# AMT Bridle Pro<sup>®</sup>

Nasal Tube Retaining Systems

Directions, Indications, & Contraindications For Use

Made in USA

[www.AppliedMedical.net](http://www.AppliedMedical.net)

# AMT Bridle™ Family of Nasal Tube Retaining Systems

Directions, Indications, & Contraindications For Use

## CONTENTS

	2 - 3
<b>ENGLISH</b>	4 - 7
<b>БЪЛГАРСКИ</b>	8 - 11
Český	12 - 15
<b>DANSK</b>	16 - 19
<b>NEDERLANDS</b>	20 - 23
Eesti keel	24 - 27
<b>SUOMI</b>	28 - 31
<b>FRANÇAIS</b>	32 - 35
<b>DEUTSCH</b>	36 - 39
<b>ΕΛΛΗΝΙΚΑ</b>	40 - 43
<b>MAGYAR</b>	44 - 47
<b>ITALIANO</b>	48 - 51
日本語	52 - 55
<b>LATVIEŠU</b>	56 - 59
<b>LIETUVIŠKAS</b>	60 - 63
<b>NORSK</b>	64 - 67
<b>POLSKIE</b>	68 - 71
<b>PORTUGUÊS</b>	72 - 75
<b>PORTUGUÊS</b>	76 - 79
<b>ROMÂNĂ</b>	80 - 83
<b>SLOVENSKÝ</b>	84 - 87
<b>SLOVENŠČINA</b>	88 - 91
<b>ESPAÑOL</b>	92 - 95
<b>SVENSKA</b>	96 - 99
<b>TÜRK</b>	100 - 103
عربي	104 - 107













108 - 109

Made in **America**  
Cleveland, Ohio

**AMT** Headquarters, Brecksville Ohio

# AMT Bridle™

## Family of Nasal Tube Retaining Systems

**MR** The Bridle clip, tape, and catheter tubing are **MR safe**

**MR** The Bridle probe, catheter, and stylet are **MR unsafe**

**DEHP** Not made with **DEHP**

**NR** Not made with **natural rubber latex**

**MD** **Medical Device**

**S** **Single Use Only**

**Caution:** Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

A thorough nasal exam is recommended before AMT Bridle™ or AMT Bridle Pro® placement in patients to ensure that adequate bone remains in the posterior septum to support the AMT Bridle System. A nasal exam should also be conducted after AMT Bridle™ or AMT Bridle Pro® placement to ensure that the umbilical tape or bridle tubing loop has passed behind the nasal septum rather than through a septal perforation.

**WARNING:** It is critical that the device is properly secured around the vomer bone and that a nasal exam is performed after placement to ensure placement is correct. If the device is not properly secured around the vomer bone, such as through a septal deviation, excessive tension applied to the device may cause extreme damage to or removal of the septum.

**NOTE:** System is supplied non-sterile, for single use only. Please inspect all contents of the kit for damage. If damaged, do not use product.

**INTENDED USE -** The Bridle Family is intended to be used by pediatric and adult patients to secure nasal feeding tubes in place in order to optimize nutrient delivery and reduce feeding tube pull out.

**INDICATIONS FOR USE -** The AMT Bridle™ or AMT Bridle Pro® - Nasal Tube Retention System is indicated to prevent inadvertent displacement or removal of Nasogastric/Nasointestinal (NG/NI) tubes.

**CONTRAINDICATIONS FOR USE -** This device is contraindicated for patients with nasal airway obstructions or abnormalities, and facial and/or cranial fractures. Do not use on patients with thrombocytopenia (<100 k/u) or immediately post septoplasty. Do not use on patients with a graft vomer bone. Extreme caution should be used with premature infants and neonatal patients. Do not use on patients that may pull on the AMT Bridle™ or AMT Bridle Pro® to such a degree as to cause serious injury.

### CLINICAL BENEFITS

Clinical benefits to be expected when using the AMT Bridle™ or AMT Bridle Pro® include but are not limited to:

- Dramatically reduces nasal tube pullout
- Minimizes interruption of tube feedings
- No tape or sutures necessary for securement
- Optimal nutrition leads to optimal recovery
- Savings on tubes, x-rays, and nursing time

### PERFORMANCE CHARACTERISTICS

Performance characteristics of the AMT Bridle™ or AMT Bridle Pro® include but are not limited to:

- MR safe after proper placement
- No patient sedation required
- May be placed in less than a minute
- Unique design with firm grip to secure nasal tubes from sizes 5F to 18F
- Secured without messy adhesive tape or sutures
- Designed for use with ANY brand of nasal tube

**WARNING:** Excessive traction on the AMT Bridle™ or AMT Bridle Pro® may cause tube displacement or nasal injury; an alternative means of securing the nasal tube should be explored.

**WARNING:** Placement of the bridle may be very difficult or impossible for pediatric patients who are nasally intubated.


**NOTE:** No patient sedation is required for placement of the AMT Bridle™ or AMT Bridle Pro®. However, sedation may be used at the healthcare provider's discretion if appropriate.

**NOTE:** These devices are intended to be placed only by qualified professionals.


**WARNING:** This device is intended for single use. Do not reuse or reprocess this medical device. Doing so may compromise biocompatibility characteristics, device performance and/or material integrity; any of which may result in potential patient injury, illness and/or death.

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

It is recommended that the bridle is placed prior to any nasal tubes for pediatric patients.

**Made in America**  
**Cleveland, Ohio**

## Kit Contents:

- Retrieval Probe
- Stylet Guide
- Flexible Bridle Catheter
- Retaining Clip
- Clip Opening Tool
- Lubricant (Not Shown)

## Clip Example:

**Pro Range Clip**

### Pro Range Clip

- 5-6F – Clip is PURPLE
- 8-10F – Clip is TEAL
- 12-14F – Clip is BLUE
- 16-18F – Clip is YELLOW

**Standard Clip/  
Pro Clip**

### Pro Clip

- 8F – Clip is WHITE
- 10F – Clip is TEAL
- 12F – Clip is BLUE

### Pro Clip

- 14F – Clip is YELLOW
- 16F – Clip is YELLOW
- 18F – Clip is Lt. BLUE

**Standard Clip/  
Pro Clip**

### Standard Clip

- 8F – Clip is WHITE
- 10F – Clip is TEAL
- 12F – Clip is BLUE

### Standard Clip

- 14F – Clip is YELLOW
- 16F – Clip is YELLOW
- 18F – Clip is Lt. BLUE

- 1** **NOTE:** It is preferred to have the patient arranged in supine position prior to placement. The AMT Bridle Pro® may be placed before or after the nasal tube. **Important: Lubricant is not necessary for use with this AMT Bridle Pro® system.** However, lubricant or water may be used on the probe or top of the AMT Bridle Pro® tubing if desired. Insert the probe into nostril opposite the nasal tube until the first rib is at the bottom of the nostril. Make adjustments for smaller patients.

**NOTE:** The stylet and catheter should be inserted along the floor of the nostril, towards the back of the throat – not towards the patient's eyes.

**NOTE FOR PEDIATRIC PATIENTS:** It is recommended that the AMT Bridle Pro® is placed prior to the nasal tube for pediatric patients. Placement of the AMT Bridle Pro® after the nasal tube may be possible, but may be more difficult due to the limited intranasal area.

**Diagrams Page 2-3:**

**1: Insert probe to first rib**

**2: Nasal Tube**

- 3** If necessary, gently twist the probes from side to side and/or up and down to encourage contact between the magnets. If no contact has occurred, then advance both the bridle tubing and the probe. **Important: Once contact has occurred, remove the orange stylet completely from the bridle tubing.**

**WARNING:** Moving or non-cooperative (especially pediatric) patients may cause additional risks during placement. Nasal anesthetic spray may be used to ease patient comfort. If appropriate, sedation may be used to assist placement. In these cases, please consult physician.

**Diagrams Page 2-3:**

**10: It may be necessary to advance or manipulate the probe to achieve magnet contact.**

• Remove stylet **AFTER** magnet connection

- 5** Cut and dispose of the bridle tube section containing the magnet, black markings and the blue probe.

If the nasal tube has not been placed, insert now per the nasal tube manufacturer's instructions.

Wipe down the nasal tube and bridle tubing if needed.

**Diagrams Page 2-3:**

**13: Bridle Tubing**

- 2** Insert the AMT Bridle Pro® tubing into the opposite nostril to approximate the magnets. Pull back on the orange stylet about 1cm (1/2") until the magnets connect. Lining up the depth indicator marks on the probe and catheter may assist in magnet connection. The "click" of the magnets may be heard or felt. Equal lengths of both probes (minus the orange stylet) should be exposed.

**NOTE FOR PEDIATRIC PATIENTS:** Magnet connection may be more difficult for pediatric patients due to the smaller area for probe manipulation or due to soft tissue constraints related to the patient's intranasal anatomy.

**Diagrams Page 2-3:**

**3: Stylet Handle**

**4: Bridle Tubing**

**5: Insert Catheter**

**6: Vomer**

**7: Probe With Magnet**

**8: Catheter With Magnet**

**9: Palate**

- 4** Slowly withdraw the probe and allow the bridle tubing to advance through the nose. Continue until the two black markings on the bridle tubing are completely pulled up and through one nostril and at least a couple inches outside the opposite nostril. This creates a loop or "bridle" around the vomer bone. If bridle tubing does not advance out of the opposite nostril, remove bridle tubing, replace the stylet and start over at Step 1.

**Diagrams Page 2-3:**

**11: Withdraw Probe**

**12: Catheter and bridle tubing will advance to form loop**

- 6** **Important: Slide the clip up the AMT Bridle Pro® tubing into position.** The clip should be located close to nostril, approximately 1cm (1/2") or one "finger width" from nostril.

**NOTE:** The clip should not touch the nostril.

**PRO RANGE CLIP:**

Place the nasal tube into the clip's circular nasal tube region.

**PRO CLIP: 8, 10, 12FR**

Place the nasal tube into the channel or groove.

**PRO CLIP: 14, 16, 18FR**

Place the nasal tube into the channel or groove.

**Diagrams Page 2-3:**

**14: Nasal Tube**

**15: Clip**

**16: Nasal Tube Location (Top Clip)**

**17: Nasal Tube Location (Middle Clip)**

**18: Nasal Tube Location (Bottom Clip)**

**7 PRO RANGE CLIP:**

Ensure nasal tube fits securely. Place loose strand of bridle tubing between the clear flats below the circular region of the clip.

**PRO CLIP: 8, 10, 12FR**

Ensure nasal tube fits securely in channel. Place loose strand of bridle tubing in the clip's hinge.

**PRO CLIP: 14, 16, 18FR**

Ensure nasal tube fits securely in channel. Place loose strand of bridle tubing between the softer inner portion of the clip and the outer, more rigid plastic section.

**Diagrams Page 2-3:**

**19: Pro Range Clip**

**20: Nasal Tube**

**21: Clip**

**22: Pro Clip 8, 10, 12FR**

**23: Nasal Tube**

**24: Clip**

**25: Pro Clip 14, 16, 18FR**

**26: Nasal Tube**

**27: Clip**

- 9** After the clip has been fully closed, tie the two strands of bridle tubing together in a simple knot below the clip.

**Diagrams Page 2-3:**

**28: 1 cm**

**29: Gently pull on bridle tubing ends to verify clip is closed**

**30: Hold tube fixed**

- 11 Clip Opening Tool:** If the clip must be opened, place the side of the opening tool into the clip. Push in and turn slightly to open.

**Diagrams Page 2-3:**

**32: Clip Opening Tool**

**8 PRO RANGE CLIP:**

Close the clip by gripping it near the base of the clip and firmly snapping it closed.

**PRO CLIP: 8, 10, 12FR**

Close the clip by folding over the plastic edge and firmly snapping it closed.

**PRO CLIP: 14, 16, 18FR**

Close the clip by folding over the plastic edge and firmly snapping it closed.

**NOTE:** Ensure proper position of the nasal tube, bridle tubing and clip prior to closure.

**ALSO NOTE:** Some nasal tubes are more rigid than others and may require a higher force to close clip.

**WARNING:** Do not use any foreign object to open or close the clip as this may damage the clip, making it less secure and/or ineffective. If damaged, use a new clip.

- 10** Cut the excess bridle tubing 1cm (1/2") below the knot.

**Note** placement of the AMT Bridle Pro® in the patient's chart.

**WARNINGS:** Throughout the useful life of the AMT Bridle Pro®, the clip must be visually inspected for signs of damage or loose attachment to the nasal tube. Monitor the clip and nasal tube position regularly for nasal tube migration.

**Diagrams Page 2-3:**

**31: Cut excess bridle tubing**

- 12 REMOVAL:** To remove the AMT Bridle Pro® and nasal tube: cut only ONE strand on the bridle tubing. Gently pull both the AMT Bridle Pro and nasal tube out of the nose.

To remove only AMT Bridle Pro®: Cut only ONE strand of the bridle tubing and open the clip. Gently pull on the opened clip to remove the Bridle tubing from nose.

**NOTE:** After removal, dispose of the device following facility protocol, local disposal guidelines, or through conventional waste.

**CAUTION:** Great care must be taken to prevent accidental swallowing as the bridle tubing can easily slide back into the nose.

**NOTE:** Recommended for up to 30 days continuous use.

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© 2021 Applied Medical Technology, Inc. **C-4614-Rev D 02/2021**
**Applied Medical Technology, Inc.**