Removal Tool

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
Removal Tool

For Use with MiniONE® Non-Balloon Button

DIRECTIONS FOR USE

INTENDED USE

The Removal Tool is intended to ease the placement and removal of Non-Balloon low profile feeding devices. The Removal Tool is intended for use by qualified clinicians for patients of all age groups (infants through elderly).

INDICATIONS FOR USE

The Low Profile Non-Balloon Feeding Device is indicated for use in patients who require long term feeding, are unable to tolerate oral feeding, who are at low risk for aspiration, require gastric decompression and/or medication delivered directly into the stomach through a secured (initial placement) or formed (replacement) stoma.

NOTE: This device is an accessory to the MiniONE® Non-Balloon Button and MiniONE® Capsule Non-Balloon Button.

CONTRAINDICATIONS

At this time there are no contraindications for use.

COMPLICATIONS

Complications may include, but are not limited to, obstruction, peritonitis, minor wound infections at the stoma site, leakage of gastric contents, gastrocolic fistula, and sepsis.

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

CLINICAL BENEFITS

Clinical Benefits to be expected when using the Removal Tool include but are not limited to:

- Eases the placement and removal of non-balloon low profile feeding devices
- Allows for the insertion of obturated Non-Balloon Buttons, which provide direct access to the stomach for feeding, decompression, and channeling medication
- Decreases the risk of patient discomfort and injury during the removal of non-balloon low-profile feeding devices

PERFORMANCE CHARACTERISTICS

Performance Characteristics of the Removal Tool include but are not limited to:

- Insertion/Removal tool distends internal bolster to a smaller profile
- Snap arm component assists in easier removal
- Various sizes to accommodate different device lengths

DEVICE PLACEMENT PROCEDURE

WARNING: PRIOR TO PLACEMENT, PLEASE INSPECT ALL CONTENTS OF THE KIT FOR DAMAGE. IF THE PACKAGE IS DAMAGED OR STERILE BARRIER IS BREACHED, DO NOT USE THE PRODUCT.

At the time of placement, test the direction of the stoma using a soft catheter of the appropriate dimension (e.g., a nasogastric tube). *Insertion at an improper angle may disrupt the stoma tract, possibly causing stomach separation from the abdominal wall, which can lead to tissue damage, infection, sepsis and/or peritonitis.*

1. Remove the existing gastrostomy tube from the stoma using the original manufacturer's instructions. (For removal see "Removal of the MiniONE® Non-Balloon Button".)
2. Insert T-Handle into the MiniONE® Non-Balloon Button by first rotating reinforcer up away from the tip of the T-Handle and then lining up the key with the key opening along the button centerline (see **Fig A1**). Insert the T-Handle fully and turn clockwise 1/2 turn (see **Fig A2**).

DEVICE PLACEMENT PROCEDURE

3. Lubricate the rod of the Snap Arm, and then slowly slide the Snap Arm through the T-Handle & into the MiniONE[®] Non-Balloon Button elongating it (see **Fig A3**). Verify rod tip is centered in the button dome.
4. Ratchet the Snap Arm down all the way to the bottom of the T-Handle until it stops.
5. Lubricate the tip of the elongated MiniONE[®] Non-Balloon Button and the stoma site with a water soluble lubricant.
6. Holding the T-Handle, gently guide the device through the stoma and into the stomach until the external flange is flush to the skin (see **Fig A4**).
7. Release and remove the Snap Arm by briefly pushing it in (while holding the T-Handle)—then squeezing the triangular ends of the Snap Arm and sliding it out (see **Fig A5**).
8. After removing the Snap Arm, disengage T-Handle by carefully rotating the metal reinforce up and turning T-Handle counterclockwise until it stops and remove it (see **Fig A6**).
9. Verify that the external flange of the MiniONE[®] Non-Balloon Button rotates freely.

NOTE: After the MiniONE[®] Non-Balloon Button has been placed, stomach contents may initially leak out of the anti-reflux valve in the feeding port area. This is a normal occurrence and occasionally happens due to the anti-reflux valve being activated by the obturator rod. To stop stomach contents from leaking, reactivate the anti-reflux valve by attaching a feeding set adapter (or the T-Handle). Once the valve has been reactivated, remove the feeding adapter (or T-Handle) and check again for stomach content leakage.

REMOVAL OF THE MINIONE[®] NON-BALLOON BUTTON

WARNING: SPONTANEOUS CLOSURE OF STOMA MAY OCCUR AS EARLY AS TWENTY-FOUR (24) HOURS AFTER REMOVAL. INSERT A NEW DEVICE IF ENTERAL FEEDING BY THIS ROUTE IS STILL INTENDED. IF CLOSURE IS DESIRED, APPLY A DRESSING OVER THE STOMA SITE. REMOVAL SHOULD NOT TAKE PLACE UNTIL THE STOMA IS FULLY ESTABLISHED, WHICH MAY TAKE TWO MONTHS OR MORE. REMOVAL IS TO BE PERFORMED BY A QUALIFIED CLINICIAN.

1. Lubricate the stoma site underneath the existing MiniONE[®] Non-Balloon Button with a water soluble lubricant. (A topical anesthetic gel may also be applied to the stoma and surrounding skin.)
 2. Using the T-Handle with metal reinforcer, first rotate reinforcer up away from the tip of the T-Handle (see **Fig B1**).
 3. Insert T-Handle into the MiniONE[®] Non-Balloon Button by lining up the key with the key opening along the button centerline (see **Fig B2**). Insert the T-Handle fully and turn clockwise 1/2 turn (see **Fig B3**).
 4. With the T-Handle in the 1/2 turn position, rotate the metal reinforcer down - sliding it underneath the bolster of the MiniONE[®] Non-Balloon Button (see **Fig B4**).
 5. Lubricate the rod of the Snap Arm, and then slowly slide the Snap Arm through the T-Handle and into the MiniONE[®] Non-Balloon Button elongating it (see **Fig B5**).
- NOTE:** Resistance should be observed when ratcheting the Snap Arm down. If no resistance is observed when ratcheting down the Snap Arm, remove the Snap Arm & T-Handle and remove the device via traction.
6. Ratchet the Snap Arm all the way down to the bottom of the T-Handle until it stops against the metal reinforcer arm (see **Fig B6**).
 7. Gently remove the device from the stoma site (see **Fig B7**).
 8. Replace as necessary following the new device's Directions for Use.

DEVICE REMOVAL PROCEDURE

Additional Methods for Removal:

- a. The MiniONE[®] Non-Balloon Button can also be traction removed without elongating the device.
 - i. Simply lubricate the stoma site underneath the existing MiniONE[®] Non-Balloon Button with a water soluble lubricant. (A topical anesthetic gel may also be applied to the stoma and surrounding skin.)
 - ii. Grasp the tabs, apply constant traction and carefully remove the device.
- b. Endoscopic removal is one other option.
 - i. Insert endoscope, cut shaft and retrieve dome.
 - ii. Do NOT let the dome pass through the intestinal tract.

WARNING: IF THE TUBE IS RESISTANT TO REMOVAL, LUBRICATE THE STOMA SITE WITH A WATER-SOLUBLE LUBRICANT. ROTATE THE TUBE GENTLY AND PUSH IT IN ABOUT 2.5 CM. CAREFULLY WORK THE TUBE FREE. NEVER USE FORCE TO REMOVE THE TUBE.

- c. If traction provides too much stress to the patient, the catheter can be cut at skin level and the internal bolster removed endoscopically.

WARNING: DO NOT ALLOW THE INTERNAL BOLSTER (CUT CATHETER SECTION) TO PASS THROUGH THE INTESTINAL TRACT.

WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, RE-STERILIZE OR REPROCESS THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.

NOTE: The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

STORAGE

Store in a controlled room temperature environment.

THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

 Sterilized using Ethylene Oxide	 Prescription Only
 Single Use Only	 Do not resterilize
 Not made with DEHP (Di(2-ethylhexyl) Phthalate)	 Not made with natural rubber latex.
 Do not use if package is damaged and consult instructions for use	 Medical Device

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© 2024 Applied Medical Technology, Inc.
Applied Medical Technology, Inc.