Capsule Monarch® & Capsule Dome

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Capsule Monarch® / Capsule Dome

Pre-loaded Replacement Gastrostomy Tube

DIRECTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician. Components are supplied sterile: For single user only. Do not re-use or re-sterilize.

INTENDED USE

The Capsule Monarch® or Capsule Dome G-Tube is intended to be used as an effective replacement gastrointestinal feeding device. The Capsule Monarch® or Capsule Dome G-Tube has a non-balloon internal bolster packaged in a dissolvable capsule to aid in the ease of insertion. The Capsule Monarch® or Capsule Dome G-Tube are intended to be used by clinicians and caregivers for infant, child, adolescent, adult, and elderly patients.

INDICATIONS FOR USE

The Monarch® / Dome G-Tube (Pre-Loaded Replacement G-Tube – with Capsule) is to be used as a percutaneous replacement gastrostomy tube for a patient with a well established gastrostomy tract. This device will assist in providing nutrition directly into the stomach through an established stoma in a patient who is unable to consume nutrition by conventional means. The Monarch® / Dome G-Tube can also deliver medication and allow for decompression of the stomach.

CONTRAINDICATIONS

Contraindications for this device are those specific to patients where there is: Lack of adherence of the stomach to the abdominal wall, lack of a well established gastrostomy site, stoma irritation, evidence of infection, presence of multiple fistulous tracts, uncertainty as gastrostomy tract direction, or evidence of tissue granulation. **The product must NEVER be used in the vasculature.**

KIT CONTENTS

Capsule Monarch® or Capsule Dome G-Tube, Clamp, Y-Port Adapter*, ENFit® Single-Port Adapter*, Dual ENFit® Y-Port Accessory*, Guidewire* (*included with some configurations)

COMPLICATIONS

Contact your physician if you experience any of the following:

Nausea, vomiting, abdominal bloating, or diarrhea • Pain, bleeding, and/or inflammation at the G-tube site • Crusting at the stoma site • Skin around the stoma site that is red, discolored, or raw • Stoma site drainage and/or puss that is white, yellow, or green and may smell bad • Repetitive leakage of food or stomach contents • Distended stomach • Fever • Tube clog • Migration of the tubing resulting in the device extending longer or shorter from the stoma • Extubation resulting in the device becoming removed from the stoma and unable to be replaced easily • Distinct indentation at the G-tube site or a distinct gap between the device and skin.

Potential complications when using the Capsule Monarch®/Dome G-Tube include but are not limited to:

Aspiration • Abscess, wound infection and skin breakdown • Hypergranulation tissue • Buried bumper syndrome • Pressure necrosis • Gastrointestinal bleeding and/or ulcerations • Ileus or gastroparesis • Intraoperative leakage • Bowel and gastric volvulus • Peritonitis • Gastrocolic fistula • Sepsis • Obstruction

For additional information, please reference our Enteral Nutrition and Troubleshooting Guide on the Resources Page of our website www.appliedmedical.net

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

CLINICAL BENEFITS, PERFORMANCE CHARACTERISTICS

Clinical Benefits to be expected when using the Capsule Monarch® or Capsule Dome G-Tube include but are not limited to:

- Provides direct access to the stomach for feeding, decompression, and channeling medication
- Ideal for patients who do not tolerate balloon buttons due to anatomy or differences in gastric environment
- Capsule Monarch® internal bolster occupies less gastric space – ideal for patients with limited intraluminal capacity
- Internal silicone retention bolster may last longer than balloon devices
- Adjustable external bolster to accommodate various patient sizing needs
- Made from medical-grade silicone to reduce irritation and improve patient comfort
- Encapsulated bolster reduces pain during insertion
- Potential to reduce hypertrophic granulation tissue

Performance Characteristics of the Capsule Monarch® or Capsule Dome G-Tube include but are not limited to:

- Internal silicone retention bolster holds the device in place and prevents the device from being pulled out of the stomach
- Self-lubricating capsule falls away, rapidly dissolves, and passes after placement
- External bolster minimizes contact and allows skin to breathe
- Capsule Monarch® internal bolster designed to conform to gastric wall
- Guidewire compatible
- Radiopaque for x-ray confirmation of device placement

DEVICE MATERIALS

The patient contacting materials in the Capsule Monarch® G-Tube and Capsule Dome G-Tube include: Medical-grade silicone (65%) • Medical-grade thermoplastic (5%) • Medical-grade silicone w/ a medical-grade thermoplastic barb (29%) • Medical-grade silicone pad print ink (1%) • Medical-grade suture (not present after device is in place) • Medical-grade adhesive (not present after device is in place) • Medical-grade cellulose (not present after device is in place)

NOTE: The Capsule Monarch® / Dome G-Tube is made with a medical-grade cellulose capsule that rapidly dissolves and passes after the device is placed.

STORAGE

The Capsule Monarch® or Capsule Dome G-Tube should be stored under the following conditions: Keep away from direct source of heat e.g. radiators, sunlight • Protect from moisture

OTHER PRECAUTIONS

The Capsule Monarch® or Capsule Dome G-Tube should only be used by or under the supervision of the personnel trained in percutaneous gastrostomy tube replacement. A thorough understanding of the technical principles, clinical applications, and risks associated with percutaneous gastrostomy tube placement is recommended prior to using this device.

NOTE: Device placement and removal is to be performed by a qualified clinician.

PLACEMENT PROCEDURE

CAUTION: Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

WARNING: KEEP THE CAPSULE MONARCH® OR CAPSULE DOME G-TUBE AWAY FROM ALL MOISTURE UNTIL IT IS READY TO BE PLACED IN PATIENT.

1. Remove the existing gastrostomy tube per the manufacturer's directions for use.
2. Select a device with equivalent French Size (or smaller) than the product just removed (Consult with the physician to determine the size to be used). Remove the Monarch® / Dome G-Tube from the package.
3. Cleanse the skin around the stoma site.
4. Remove the cover from the Capsule Monarch® or Capsule Dome G-Tube by sliding it off the end of the capsule.
5. Lubricate the tip of the Capsule Monarch® or Capsule Dome G-Tube with a water soluble lubricant, thoroughly covering the capsule with lubricant. Also be sure to place lubricant around the stoma site.

WARNING: THE CAPSULE MONARCH® OR CAPSULE DOME G-TUBE MUST BE PLACED IMMEDIATELY AFTER LUBRICATION TO PREVENT THE CAPSULE FROM DISSOLVING PREMATURELY.

WARNING: DO NOT USE OIL OR PETROLEUM BASED LUBRICANT.

6. Grasping the tube mid-way with thumb and finger, place the capsule tip of the Capsule Monarch® or Capsule Dome G-Tube into the established stoma and gently advance tube through the stoma (see figure 1).
7. Apply gentle pressure as the Capsule Monarch® or Capsule Dome G-Tube is guided through the tract until the 5 cm marking is visible just above the stoma site.

WARNING: WHEN GUIDING INTO THE TRACT, IF ANY RESISTANCE IS EXPERIENCED, DISCONTINUE INSERTION. DO NOT USE EXCESSIVE FORCE TO AVOID DAMAGING THE TRACT OR GASTRIC WALL.

8. To release the inner bolster, grasp the tube portion to secure in place, while steadily pulling on the suture removal tab (see figure 2).
- NOTE:** If steady force is applied to the suture removal tab and the bolster will not deploy, ensure that the encapsulated dome is fully inserted into the stomach. The device can be held within the stomach for a moment to allow the capsule to soften, making for easier deployment of the inner dome.
9. Pull the suture tab until the thread is completely removed. Discard the thread and the suture pull tab.
 10. While keeping the device in the desired position and holding the tube with one hand, use your other hand to gently remove the obturator and discard it.

NOTE: The capsule will remain in the stomach, where it will dissolve within the gastro-intestinal tract.

11. Gently withdraw the tube until slight tension is felt from the internal bolster contacting the internal stomach wall. Gently slide the external bolster into the desired position to ensure secure placement of the tube. Allow for slight in-out play of the tube by sliding the external bolster back to the nearest centimeter marking. Document in the patient's chart the Bolster location on the G-Tube. The Capsule Monarch® or Capsule Dome G-Tube should rotate easily.
12. Install the clamp and Y-Port as shown in Figure 3.

WARNING: CREATING TOO SNUG A FIT MAY CAUSE EMBEDDING WITH EROSION INTO THE GASTRIC WALL. RESULTING IN TISSUE NECROSIS, INFECTION, PERITONITIS, SEPSIS, AND ASSOCIATED SEQUELAE. THE EXTERNAL BOLSTER SHOULD NOT BE SUTURED INTO PLACE.

VERIFY TUBE IS WITHIN STOMACH BEFORE FEEDING

1. Aspirate for gastric contents. Spontaneous return of gastric contents should occur. If gastric contents are not present, attach the 60 ml catheter tip syringe to the Capsule Monarch® or Capsule Dome G-Tube and irrigate with approximately 10 ml of water. Aspirate again for gastric contents. Spontaneous return of gastric contents indicates that proper placement has been achieved.
 2. Radiographic examination (X-Ray) may be performed to evaluate for location of the internal bolster.
 3. A flexible guidewire may be temporarily placed in the center of the device to further facilitate radiographic examination.
 4. Fiberoptic exam may be performed through the lumen of the Capsule Monarch® or Capsule Dome G-Tube using a 3 mm Endo-scope.
 5. After air and/or gastric contents are observed, flush with water.
- When placement is confirmed, gastric decompression, medication, or feeding channeling may begin.

PATIENT CARE INSTRUCTIONS

- The Capsule Monarch® or Capsule Dome G-Tube should be replaced periodically for optimal performance. Frequent performance checks are recommended. Clogging and/or reduced flow are indicators of diminished performance.
- Be sure the internal bolster is within the stomach, and the Capsule Monarch® or Capsule Dome G-Tube rotates freely, and has approximately 1 cm of play (in-out), before feedings begin.
- The stoma area should be cleansed daily with a mild soap and water. The stoma site should be clean and dry at all times.
- To avoid clogging of the tube, flush the Capsule Monarch® or Capsule Dome G-Tube Device with at least 10 ml of water to rinse the feeding passage of any particulate after each feeding.

CAUTION: Tubing should be monitored for possible inward migration or unintentional extubation.

WARNING: NEVER INJECT AIR INTO THE G-TUBE.

WARNING: TUBE MIGRATION COULD RESULT IN THE FOLLOWING: INABILITY TO FEED, OBSTRUCTION, PERITONITIS, INFECTION, AND ASSOCIATED SEQUELAE.

WARNING: IN CASE OF FEVER, GASTRIC DISTENTION, INFECTION, BLOCKAGE, OR TISSUE NECROSIS, PATIENTS SHOULD SEE THEIR PHYSICIAN IMMEDIATELY.

MEDICATION/ NUTRITION CHANNELING

The Capsule Monarch® or Capsule Dome G-Tube Device has been designed to be a feeding/medication/decompression access into the stomach. Applications other than those indicated in these instructions, are not recommended. Not to be used intravascularly.

WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO COMPATIBLE ENTERAL DEVICES. DO NOT USE FOR NOT ENTERAL APPLICATIONS.

WARNING: WHEN USING A BOLUS STYLE CONNECTOR, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO BREATHING SYSTEM, LIMB CUFF, AND NEURAXIAL CONNECTORS.

WARNING: ONLY TIGHTEN BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.

WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.

NUTRITION CHANNELING:

1. The Capsule Monarch® or Capsule Dome G-Tube device may be used for syringe, gravity or pump feeding, or decompression.
2. Attach the opposite end of the connector being used. Catheter and luer connectors are compression fit- firmly insert while using a slight rotation into the device to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the device while avoiding excessive force to lock in place. Once connected, feeding can begin.
3. When feeding is complete, flush with 5-10ml water. Unscrew a rotating connector with a counter-clockwise turn. Catheter and luer connectors can be removed using moderate tension. Snap plug in place to keep lumen clean.

MEDICATION CHANNELING:

CAUTION: Medication should be channeled separately from feeding, one dose at a time.

- Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in water before channeling through the feeding tube. Never crush enteric coated medication or mix medication with formula.
- Using a catheter tip or ENFit® syringe flush the tube with the prescribed amount of water.

PROPER FLUSHING GUIDELINES

Proper flushing techniques and frequency can help prevent tubing blockage, clogs, and tubing failure. Follow these flushing guidelines to maintain optimal device flow conditions:

- Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain potency.
- Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube before and after medication channeling and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

UNCLOGGING A DEVICE:

First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a catheter tip or ENFit® syringe to the port. Fill the syringe with warm water and gently push and pull the syringe plunger to free the clog. It make take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced.

- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract. Use a 30 to 60 ml catheter tip syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.

CAUTION: Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.

MRI SAFETY INFORMATION

The Capsule Monarch® G-Tube and Capsule Dome G-Tube are MR Safe once placed.

REMOVAL

WARNING: REMOVAL IS RECOMMENDED TO BE PERFORMED BY A QUALIFIED CLINICIAN. SPONTANEOUS CLOSURE OF STOMA MAY OCCUR AS EARLY AS TWENTY-FOUR (24) HOURS AFTER REMOVAL. INSERT A NEW DEVICE IF ENTERAL FEEDING BY THIS ROUTE IS STILL INTENDED. IF CLOSURE IS DESIRED, APPLY A DRESSING OVER THE STOMA SITE.

1. Grasp the tube near the skin and place the other hand around the stoma site. Gently pull upward on the gastrostomy tube with a firm constant motion. See Figure 4.

WARNING: IF THE TUBE IS RESISTANT TO REMOVAL, THE CLINICIAN IS RECOMMENDED TO LUBRICATE THE STOMA SITE WITH A WATER-SOLUBLE LUBRICANT. ROTATE THE TUBE GENTLY AND PUSH IT IN ABOUT AN INCH. NEVER USE FORCE TO REMOVE THE TUBE. SEDATION MAY BE REQUIRED TO COMPLETE THE PROCEDURE.

2. If traction provides too much stress to the patient, the catheter can be cut at skin level and the internal bolster removed endoscopically.

WARNING: DO NOT ALLOW THE INTERNAL BOLSTER (CUT CATHETER SECTION) TO PASS THROUGH THE INTESTINAL TRACT.

CAUTION: AMT recommends that a balloon style device is kept on hand as a spare in case the Capsule Monarch® or Capsule Dome G-Tube device becomes pulled from the stoma or other failure occurs. Temporary placement of a spare device would prevent the stoma site from closing while scheduling a device replacement.

WARNING: AFTER REMOVAL, THE DEVICE MAY BE A POTENTIAL BIOHAZARD. THE REMOVED DEVICE SHOULD BE HANDLED AND DISPOSED OF ACCORDING TO YOUR LOCAL, STATE, AND FEDERAL LAWS AND REGULATIONS. YOUR HEALTHCARE PROFESSIONAL MAY BE ABLE TO ADVISE THE MOST ACCEPTABLE METHOD OF DISPOSAL.

TROUBLESHOOTING

Long-term device performance and functionality are dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your feeding device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

<p>A tear has formed: Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.</p>
<p>Tubing has reduced flow or has become clogged: Tubing can become blocked due to not properly flushing after each use, use of thick or improperly crushed medications, use of thick feeds/formulas, gastric reflux, and/or fungus growth. If clogged, refer to the PROPER FLUSING GUIDELINES section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced.</p>
<p>Plug will not stay closed: Ensure that the plug is being firmly and fully pressed into the interlock connector. If the plug is not staying closed, check the plug and feed-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.</p>
<p>Device has pulled from stoma: Excessive force may have been applied to the device during use. Avoid excessive force while using the device. Device will need to be placed promptly to avoid closure of the stoma site.</p>
<p>Gastric leakage: Gastric leakage may occur if the sliding external bolster is not properly adjusted or if the internal bolster is not against the stomach wall. Gently withdraw the tube until slight tension is felt from the internal bolster contacting the internal stomach wall. Gently slide the external bolster into the desired position so that it sits 1-2 mm from the skin</p>
<p>Foul smell coming from the device: Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed and stoma site should be gently cleaned with soap and warm water. If foul smell does not go away, it is recommended that you contact your healthcare professional.</p>
<p>Device has become discolored: The device can become discolored over days to months of use. This is normal depending on the types of feeds and medications being used with the device.</p>

DEVICE LONGEVITY

Feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness.

Exact device longevity cannot be predicted. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 3-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, trauma to the device, contact with sharp or abrasive objects, and overall tube care.

For optimal performance, it is recommended that the Capsule Monarch® or Capsule Dome G-Tube device be changed every 6 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure. If devices are failing or performance is degrading earlier than the typical range for device longevity, it is recommended that you speak with your healthcare professional regarding eliminating common factors that can lead to early device degradation.

WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, RE-STERILIZE OR REPROCESS THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.

THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

EN	STERILE EO	Sterilized Using Ethylene Oxide	Rx Only Prescription Only
 Do Not Resterilize	 Single Use Only	 Do not use if package is damaged and consult instructions for use	 MR Safe
 Not made with natural rubber latex	 Medical Device	 Not made with DEHP (Di(2-ethylhexyl) Phthalate)	
 Keep away from direct source of heat	 Keep Dry	ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.	

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Applied Medical Technology, Inc.