# AMT MiniONE® Enteral Extension Set

CE  
2862

- ES** Equipo de extensión enteral MiniONE®
- DE** MiniONE® Enterale Verlängerungsset
- IT** Set di estensione enterale MiniONE®
- FR** Kit de rallonge entérale MiniONE®
- NO** MiniONE® enteralt forlengelsessett
- SV** MiniONE® förlängningssett
- PT-BR** Conjunto de extensão enteral MiniONE®
- NL** MiniONE® enterale uitbreidingsset
- DA** MiniONE® enteralt forlængersæt
- FI** Enteraalinen MiniONE® -jatkosarja
- AR** مجموعة التمديد المعوية MiniONE®
- BG** Коремният удължителен комплект MiniONE®
- CS** Enterální rozšiřující sada MiniONE®
- EL** Σετ Εντερικής Επέκτασης MiniONE®
- ET** MiniONE® Pikendusmoodul
- HU** MiniONE® enterális bővítőkészletet
- JA** MiniONE® 経腸延長セットは
- LV** MiniONE® enterālais paplašinājuma komplekts
- PL** Zestaw przedłużający MiniONE®
- RO** Setul extensie enteral MiniONE®
- SL** Podaljševalni komplet MiniONE®
- TR** MiniONE® Enteral Uzatma Seti
- PT** Conjunto de Extensão Enteral MiniONE®
- SK** MiniONE® enterálna predžovacia súprava
- LT** „MiniONE®“ enterinis ilginamasis rinkinys
- MS** Set Sambungan Enteral AMT MiniONE®

Rx ONLY

C4623-D 03/2024

Made in the U.S.A.

FIG A-1



FIG A-2

MIC-KEY®  
Equivalent

FIG B-1 - G-JET® (Jejunal Port™)


FIG B-2 - MiniONE® Button


FIG C



FIG D

**MIC-KEY® Equivalent**  
 Equivalencia MIC-KEY®  
 Äquivalent zu MIC-KEY®  
 Equivalente MIC-KEY®  
 Équivalent du MIC-KEY®  
 MIC-KEY® ekvivalent  
 Kompatibel med MIC-KEY®  
 Equivalente a MIC-KEY®  
 Equivalent van MIC-KEY®  
 MIC-KEY® alternatív  
 Vastaava MIC-KEY®  
 ما يعادل MIC-KEY®  
 Эквивалент на MIC-KEY®  
 Ekvivalent MIC-KEY®

MIC-KEY® Equivalent  
 MIC-KEY® Equivalent  
 MIC-KEY® egyenértékű  
 MIC-KEY® 同等  
 MIC-KEY® ekvivalents  
 Odpowiednik MIC-KEY®  
 Echivalent MIC-KEY®  
 Ekvivalent MIC-KEY®  
 MIC-KEY® Muadili  
 MIC-KEY® Equivalent  
 Ekvivalent MIC-KEY®  
 MIC-KEY® ekvivalentas  
 Setara MIC-KEY®

**\*1 Jejunal Port**

	Jejunal port	空腸ポート
Puerto yeyunal	Hjäsuoilportti	Tukšās zarnas ports
Jejunaler Port	منفذ صائمي منفذ	Port jelitowy
Porta digiunale	порт за тънки черва	Port tub gastric
Orifice jéjunal	lačníkový otvor	Jejunalni port
Jejunal port	Jejunal Port	Jejunal Giriş
Jejunal port	Jejunaalne port	Tuščiosios žarnos jungtis
Porta jejunal	Jejunal Port	Jejunálny port
Jejunale poort	jejunal kikötő	Tapak Jejunum

## INDICATIONS, CONTRAINDICATIONS

**INTENDED USE:** The AMT MiniONE® Enteral Extension Set is an effective device used in the continuous or bolus feeding and medication channeling, which connects an enteral feeding device to an enteral feed source. It is intended to be used by clinicians and caregivers for infant, child, adolescent, adult, and elderly patients. The device aids patients with a wide array of medical conditions that make maintaining an adequate nutritional state a challenge.

**INDICATIONS FOR USE:** The MiniONE® Enteral Extension Set is intended for use as an extension set for the MiniONE® Gastrostomy Button, G-JET® Gastrojejunal Feeding Tube (Jejunal Port), Mic-Key® Secure-Lok® feeding devices, or other compatible enteral feeding devices. The extension set is for administering feeding, medications, and decompression to compatible devices.

**CONTRAINDICATIONS:** At this time there are no contraindications for use associated with use of the AMT MiniONE® Enteral Extension Set.

**COMPLICATIONS:** At this time there are no complications associated with use of the AMT MiniONE® Enteral Extension Set.

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

**CLINICAL BENEFITS:** Clinical benefits to be expected when using the AMT MiniONE® Enteral Extension Set include but are not limited to: Provides a connection method between a low profile gastrostomy device and syringes, other feed sets, or feed pump tubing • Allows for connection so that nutrition, medication, and decompression can be channeled through the device to a patient

**PERFORMANCE CHARACTERISTICS:** Performance characteristics of the AMT MiniONE® Enteral Extension Set include but are not limited to: Connect a low profile gastrostomy device to syringes, other feed sets, or tubing from a feeding pump • Plugs and clamps on devices provide additional means of leakage protection • Available in multiple tube sizes to help prevent potential blockages/clogs • Glow green adapter provides better visibility in the dark for easier nighttime feeds • Device has ENFit® compatible connection options

## EXTENSION SET INSTRUCTIONS

**WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO ENTERAL DEVICES.**

**WARNING: IF USING A LEGACY-STYLE (NON-ENFIT®) EXTENSION SET, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO THE FOLLOWING SYSTEMS: BREATHING, CONES AND SOCKETS OF ANESTHETIC AND RESPIRATORY EQUIPMENT, INTRAVENOUS, LIMB CUFF, NEURAXIAL CONNECTORS, NIPPLES OF RESPIRATORY THERAPY EQUIPMENT, URINARY, AND TEMPERATURE SENSOR CONNECTORS OF RESPIRATORY HUMIDIFYING EQUIPMENT.**

1. Prior to use, inspect contents for damage. If damaged, do not use. Obtain another package.
2. The extension set may be used for syringe, gravity or pump feeding or decompression.
3. Ensure the clamp is closed and attach the extension set (Fig A 1 / 2) to the button by lining up the dark line on the extension set connector with the dark line on the button. Fully press the extension set connector into the button. Turn 3/4 to the right (clock-wise) to lock the extension set in place (Fig B 1 / 2).

## EXTENSION SET INSTRUCTIONS

4. Attach the opposite end of the extension set to the connector being used (Fig C). If using a bolus or luer connector, firmly insert the connector while using a slight rotation into the extension set to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the extension set while avoiding excessive force to lock in place. Once connected, open the clamp to allow flow.

**WARNING: ONLY TIGHTEN BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.**

**WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.**

**WARNING: IF EXTENSION SET IS NOT PROPERLY PLACED AND LOCKED, LEAKAGE MAY OCCUR. WHEN TWISTING THE CONNECTOR, HOLD THE BUTTON IN PLACE TO AVOID ROTATION WHILE PLACING THE EXTENSION SET. NEVER INSERT EXTENSION SET CONNECTOR INTO A BALLOON INFLATION PORT (FIG D); THIS MAY CAUSE THE INTERNAL BALLOON TO DEFLATE. REFER TO BUTTON DIRECTIONS FOR USE.**

5. When feeding is complete, flush with 5-10ml water. Close the clamp to prevent leakage while removing the extension set. Remove the extension set by holding the button in place, grasp the extension set connector, twist 3/4 turn to the left (counter clock-wise) to unlock, and gently remove the extension set from the button (Fig B 1/2). Snap button plug in place to keep lumen clean.
6. The extension set should be cleansed with soap and water and thoroughly rinsed. Store in a clean, dry place. Never use a dishwasher for cleaning.

## DEVICE LONGEVITY

Extension sets are meant to be periodically replaced for optimal performance and cleanliness. AMT recommends that the extension set be changed at least every 2 weeks or as often as indicated by your healthcare professional. Device performance and functionality can degrade over time depending on usage and environmental conditions. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, frequency of use, trauma to the device, and overall extension set care.

Device should be replaced if leakage, cracks, tears, excessive residue build-up, mold, fungus, or other signs of failure are noted. Some usage types can wear down device components quicker than others. If tubing is becoming hardened, replacement is recommended in near future to avoid failure.

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

**NOTE:** Refer to button Directions for Use for the circumstances in which you should consult a healthcare professional.

## THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

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© 2024 Applied Medical Technology, Inc.
MIC-KEY® is a registered trademark of Avent, Inc.
MIC-KEY® er et registreret varemærke og tilhører Avent, Inc.
MIC-KEY® on Avent, Inc.:n rekisteröity tavaramerkki.
MIC-KEY® je registrovaná obchodní známka společnosti Avent, Inc.
MIC-KEY® on Avent, Inc. registreeritud kaubamärk.