AMT Balloon G-Tube

Balloon G-Tube Replacement Gastrostomy Feeding Device

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DEHP

AMT Balloon G-Tube

Replacement Gastrostomy Feeding Device

DIRECTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician. The device is provided STERILE for single use. This device is intended to connect to compatible enteral devices only. **Do not use for non-enteral applications.**

INTENDED USE

The AMT Balloon G-Tube is an effective replacement gastrointestinal feeding device. The AMT Balloon G-Tube is intended to be used by clinicians and caregivers for infant, child, adolescent, adult, and elderly patients.

INDICATIONS FOR USE

The AMT Balloon G-Tube is to be used as a percutaneous replacement gastrostomy tube for a patient with a well-established gastrostomy tract. This device will assist in providing nutrition directly into the stomach through an established stoma in a human patient who is unable to consume nutrition by conventional means. The AMT Balloon G-Tube can also deliver medication and allow for decompression of the stomach.

CONTRAINDICATIONS

Contraindications for this device are those specific to patients where there is: lack of adherence of the stomach to the abdominal wall, lack of a well established gastrostomy site, stoma irritation, evidence of infection, presence of multiple fistulous tracts, uncertainty as to gastrostomy tract direction, or evidence of tissue granulation. **The product must NEVER be used in the vasculature.**

NOTICE:

- Please contact a professional health care giver or physician for explanation of the warnings, care, and use of the device.
- This device has been designed to provide feeding/medication/decompression access into the stomach. Other applications not advised.

KIT CONTENTS

Balloon G-Tube with External Bolster, Clamp, ENFit® adapter (included with some configurations)

COMPLICATIONS

Contact your physician if you experience any of the following: Nausea, vomiting, abdominal bloating, or diarrhea • Pain, bleeding, and/or inflammation at the G-tube site • Crusting at the stoma site • Skin around the stoma site that is red, discolored, or raw • Stoma site drainage and/or puss that is white, yellow, or green and may smell bad • Repetitive leakage of food or stomach contents • Distended stomach • Fever • Tube clog • Migration of the tubing resulting in the device extending longer or shorter than the stoma • Extubation resulting in the device becoming removed from the stoma and unable to be replaced easily • Distinct indentation at the G-tube site or a distinct gap between the device and skin.

Potential complications when using the AMT Balloon G-Tube include but are not limited to: Aspiration • Abscess, wound infection and skin breakdown • Hypergranulation tissue • Buried bumper syndrome • Pressure necrosis • Gastrointestinal bleeding and/or ulcerations • Ileus or gastroparesis • Intraperitoneal leakage • Bowel and gastric volvulus • Peritonitis • Gastrocolic fistula • Sepsis • Obstruction

For additional information, please reference our Enteral Nutrition and Troubleshooting Guide on the Resources Page of our website www.appliedmedical.net

WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, RE-STERILIZE OR REPROCESS THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

CLINICAL BENEFITS, PERFORMANCE CHARACTERISTICS

Clinical Benefits to be expected when using the Balloon G-Tube include but are not limited to:

- Provides direct access to the stomach for feeding, decompression, and channeling medication
- Balloon design aids in easy, at home replacement when necessary
- Balloon design results in a larger gastric seal surface which may help to reduce leakage and granulation tissue
- Made from medical-grade silicone to reduce irritation and improve patient comfort
- Adjustable external bolster to accommodate various patient sizing needs

Performance Characteristics of the Balloon G-Tube include but are not limited to:

- Balloon fill volume can be easily adjusted for a custom fit
- Unique "apple" shaped balloon designed to reduce leakage and reduce the risk of accidental pull-outs from the stoma site
- Inflated balloon protects the distal tip of the tube, lowering the risk of gastric ulceration
- External bolster minimizes contact and allows skin to breathe
- Guidewire compatible
- Radiopaque for x-ray confirmation of device placement

DEVICE MATERIALS

Materials in the AMT Balloon G-Tube include: Medical-grade silicone (89%) • Medical-grade thermoplastic (10%) • Medical-grade silicone pad print ink (1%)

DIRECTIONS FOR USE

The AMT Balloon G-Tube is designed for placement in a well-established gastrostomy site with no evidence of infection, fistulous tracts, or separation of the gastric wall from the fistulous tract. A thorough understanding of the technical principles, clinical applications, and risks associated with gastrostomy tube replacement is necessary before using this product. It is recommended that the AMT Balloon G-Tube French size selected be similar in size to the device being removed (unless otherwise indicated by a physician).

PLACEMENT PROCEDURE

The device is provided as sterile, for Single Patient Use. The patient or health care giver should be made aware of notes, warnings, precautions, care instructions, and information regarding this device. Inspect the package for any breach to ensure a sterile product. If the seal appears breached or there is apparent damage, do not use, and obtain another kit.

CAUTION: Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

NOTE: Device replacement can be performed by the healthcare professional or at home by the patient/caregiver. Do not attempt to replace the device until first discussing the procedure with your healthcare professional.

1. Remove the existing gastrostomy feeding device per the manufacturer's Directions For Use.

NOTE: The AMT Balloon G-Tube French size selected should be similar in size to the device being removed (unless otherwise indicated by a physician).

2. Remove the AMT Balloon G-Tube from the package and inspect for any visible defects. If defects are observed, do not use device. Obtain another package.

3. Inflate balloon through the balloon inflation port using a Luer slip syringe with distilled or sterile water to the recommended fill volume.

NOTE: If the balloon is asymmetrical, gently roll it between your fingers or your hands until the balloon frees itself and is symmetrical.

4. Reinsert the syringe and deflate all of the water from the balloon after inspection.

5. Cleanse the skin around the stoma site. Lubricate the tip of the AMT Balloon G-Tube with a water soluble lubricant. Gently guide the tube through the stoma, and into the stomach.

WARNING: DO NOT USE OIL OR PETROLEUM BASED LUBRICANT

WARNING: WHEN GUIDING INTO THE TRACT, IF ANY RESISTANCE IS EXPERIENCED, DO NOT USE EXCESSIVE FORCE TO AVOID DAMAGING THE TRACT OR GASTRIC WALL.

6. Inflate balloon through the balloon inflation port using a Luer slip syringe with distilled or sterile water to the recommended fill volume.

WARNING: NEVER USE AIR TO FILL BALLOON.

NOTE: It may be necessary to adjust the balloon volume to prevent leakage of gastric contents based on the individual patient needs.

7. Gently snug balloon up against gastric mucosa and position the external bolster to ensure secure placement of the AMT Balloon G-Tube. Allow for slight in-and out play of tube by sliding the external bolster back to the nearest cm marking. The AMT Balloon G-Tube should also rotate freely.

WARNING: CREATING TOO SNUG A FIT MAY CAUSE EMBEDDING WITH EROSION INTO THE GASTRIC WALL, RESULTING IN TISSUE NECROSIS, INFECTION, PERITONITIS, SEPSIS, AND ASSOCIATED SEQUELAE. THE EXTERNAL BOLSTER SHOULD NOT BE SUTURED INTO PLACE.

8. Check the AMT Balloon G-Tube around the stoma for gastric leakage. If the stoma leaks, increase the balloon volume by 1-2 ml through the Balloon Inflation Valve. Do not exceed the maximum fill volume. Repeat this step until the leak stops or maximum balloon fill volume is reached.

WARNING: DO NOT EXCEED MAXIMUM BALLOON FILL VOLUME.

French Size	Fill Volumes	
	Recommended Fill	Maximum Fill
12F	2-2.5 ml	3 ml
14F	3-4 ml	5 ml
16F	4-6 ml	8 ml
18F	6-8 ml	10 ml
20-24F	7-10 ml	15 ml

TUBE CARE INSTRUCTIONS

CAUTION: Tubing should be monitored for possible inward migration or unintentional extubation.

- The AMT Balloon G-Tube should be replaced periodically for optimal performance. Frequent performance checks are recommended. Clogging and/or reduced flow are indicators of diminished performance.
- Be sure the balloon is within the stomach, and the AMT Balloon G-Tube rotates freely before feedings begin.
- The stoma area should be cleansed daily with a mild soap and water. The stoma site should be clean and dry at all times.
- It is recommended that balloon volume is checked at least every two weeks, or as recommended by your healthcare professional. Aspirate all of the sterile water or saline solution from the balloon through the Balloon Fill Valve and replace the original amount prescribed.
- To avoid clogging of the feeding tube, flush the AMT Balloon G-Tube through the Feeding Port with at least 10 ml of water to rinse the feeding passage of any particulates after feeding.
- Tubing should be monitored for possible inward migration or unintentional extubation.
- The AMT Balloon G-Tube has been designed to be a feeding/medication/decompression access into the stomach. Other applications are not advised.

WARNING: NEVER INJECT AIR INTO THE G-TUBE.

WARNING: TUBE MIGRATION COULD RESULT IN THE FOLLOWING: OBSTRUCTION, INABILITY TO FEED, PERITONITIS, INFECTION AND ASSOCIATED SEQUELAE.

WARNING: IN CASE OF FEVER, GASTRIC DISTENTION, INFECTION, BLOCKAGE, OR TISSUE NECROSIS, PATIENTS SHOULD SEE THEIR PHYSICIAN IMMEDIATELY.

VERIFY TUBE IS WITHIN STOMACH BEFORE FEEDING

1. Aspirate through the Feeding Port for gastric contents. Return of gastric contents should occur. If gastric contents are not present, attach a large catheter tip syringe (attach a transition adapter or ENFit[®] syringe for certain configurations) to the AMT G-Tube Feeding Port and irrigate with 10 ml of water. Aspirate for gastric contents.
2. Radiographic examination (X-Ray) may be performed to evaluate for location of the internal bolster (balloon).
3. A flexible guidewire may be temporarily placed in the center of the device to further facilitate radiographic examination.
4. Fiber optic exam may be performed through the Feeding Port using a 3 mm Endoscope.
5. After air and/or gastric contents are observed, flush with water.

WARNING: NEVER INJECT AIR INTO THE AMT BALLOON G-TUBE FEEDING PORT OR BALLOON FILL-VALVE.

When placement is confirmed, gastric decompression or feeding channeling can begin through the Feeding Port. Refer to the manufacturer's Directions For Use for the feeding and/or medication attachments used.

MEDICATION/ NUTRITION CHANNELING

WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO COMPATIBLE ENTERAL DEVICES. DO NOT USE FOR NON-ENTERAL APPLICATIONS.

WARNING: WHEN USING A BOLUS STYLE CONNECTOR, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO INTRAVENOUS, BREATHING SYSTEM, LIMB CUFF, AND TEMPERATURE CONNECTORS.

WARNING: ONLY TIGHTEN BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.

WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.

NUTRITION CHANNELING:

1. The AMT Balloon G-Tube device may be used for syringe, gravity or pump feeding, or decompression.
2. Attach the opposite end of the connector being used. Catheter and luer connectors are compression fit- firmly insert while using a slight rotation into the device to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the device while avoiding excessive force to lock in place. Once connected, feeding can begin.
3. When feeding is complete, flush with 5-10ml water. Unscrew a rotating connector with a counter-clockwise turn. Catheter and luer connectors can be removed using moderate tension. Snap plug in place to keep lumen clean.

MEDICATION CHANNELING:

CAUTION: Medication should be channeled separately from feeding, one dose at a time.

- Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in water before channeling through the feeding tube. Never crush enteric coated medication or mix medication with formula.
- Using a catheter tip or ENFit® syringe flush the tube with the prescribed amount of water.

PROPER FLUSHING GUIDELINES

Proper flushing techniques and frequency can help prevent tubing blockage, clogs, and tubing failure. Follow these flushing guidelines to maintain optimal device flow conditions:

- Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain potency.
- Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube before and after medication channeling and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

UNCLOGGING A DEVICE:

First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a catheter tip or ENFit® syringe to the port. Fill the syringe with warm water and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced.

- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract. **Use a 30 to 60 ml catheter tip syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.**

CAUTION: Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.

REMOVAL

NOTE: Removal should only be completed by a healthcare professional or if you have been trained by a professional to remove the device.

1. Deflate the balloon through the Balloon Inflation Valve using a Luer slip tip syringe without a needle attached.
2. Gently pull to remove the AMT Balloon G-Tube from the stoma site.
3. Replace as necessary following AMT Balloon G-Tube Directions For Use.

WARNING: AFTER REMOVAL, THE DEVICE MAY BE A POTENTIAL BIOHAZARD. THE REMOVED DEVICE SHOULD BE HANDLED AND DISPOSED OF ACCORDING TO YOUR LOCAL, STATE, AND FEDERAL LAWS AND REGULATIONS. YOUR HEALTHCARE PROFESSIONAL MAY BE ABLE TO ADVISE THE MOST ACCEPTABLE METHOD OF DISPOSAL.

NOTE: Spontaneous closure of the stoma may occur as early as twenty-four (24) hours after removal. Insert a new AMT Balloon G-Tube if feeding by this route is still intended as prescribed.

WARNING: DO NOT CUT OFF THE INTERNAL BOLSTER OR TUBING, OR ALLOW TO PASS THROUGH THE INTESTINAL TRACT.

MRI SAFETY INFORMATION

The AMT Balloon G-Tube is considered MR Safe once placed.

TROUBLESHOOTING

Long-term device performance and functionality are dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your feeding device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

Gastric leakage: Gastric leakage may occur if the sliding external bolster is not properly adjusted or if the balloon is not against the stomach wall. Adjust the external bolster so that it sits 1-2 mm from the skin. Once the bolster is adjusted, check the water volume in the balloon. Do not add more water to the balloon before first withdrawing all of the water from the balloon. Fill the balloon according to the balloon inflation range provided.
Leaking of balloon volume: If balloon is deflating, fully deflate the balloon and remove from stoma. Once removed, inflate balloon with recommended fill-volume. Check balloon for leakage by gently massaging tubing and balloon. If no leakage is noted, deflate the balloon, re-insert balloon back into stoma, and re-inflate the balloon to desired fill-volume. Do not exceed maximum fill volume. Only access the balloon inflation port for balloon inflation and deflation. Other uses can result in balloon leakage or device failure. NOTE: The balloon is made of a semi-permeable material and can lose a small amount of fill-volume over time depending on environment and usage conditions.
A tear has formed: Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.
Foul smell coming from the device: Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed and stoma site should be gently cleaned with soap and warm water. If foul smell does not go away, it is recommended that you contact your healthcare professional.
Balloon inflation valve leakage: Leakage from this valve typically occurs due to residue in the fill-valve. Always use a clean syringe when inflating the balloon and only inflate water through the valve. Leakage can also occur due to valve becoming stuck if syringe is inserted too hard into valve. Insert the syringe into the balloon inflation port if leakage is occurring to attempt to reset valve. Several attempts may be necessary before valve resets.
Device has become discolored: The device can become discolored over days to months of use. This is normal depending on the types of feeds and medications being used with the device.
Balloon is misshaped: Be sure to inflate and inspect balloon prior to placement to check balloon symmetry. Balloons becoming misshaped typically occur due to excessive force or tension on the device (pulling device out of stoma while balloon is inflated). Additionally, a device that is placed too close to the pylorus can result in a misshaped balloon and blockage of the pylorus. Balloons can be gently massaged with fingers back into symmetry if not excessively misshaped. Device should be replaced if balloon is found excessively misshaped.
Balloon will not inflate or deflate: Inflation and deflation problems typically occur due to residue blocking the fill-lumen. Always use a clean syringe when inflating the balloon. In some rare occurrences, fungus can grow inside of the tubing and block the fill-lumen. Fungus growth can occur depending on patient environment and feeds/medications being channeled through the device. If balloon will not deflate, contact your healthcare professional for removal. If deflation problem is the result of fungal growth, eliminating source of fungal growth or anti-fungal medication may be required.
Tubing has reduced flow or has become clogged: Tubing can become blocked due to not properly flushing after each use, use of thick or improperly crushed medications, use of thick feeds/formulas, gastric reflux, and/or fungus growth. If clogged, refer to the PROPER FLUSHING GUIDELINES section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced.
Balloon failure: Early balloon failure can occur due to a number of patient or environmental factors, including but not limited to: gastric pH, diet, certain medications, improper balloon fill volume, placement of device, improperly channeling food/medication into balloon port, trauma, contact with a sharp or abrasive material, and overall care of the device.
Plug will not stay closed: Ensure that the plug is being firmly and fully pressed or rotated without excessive force. If the plug is not staying closed, check the plug and feed-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.

DEVICE LONGEVITY

Balloon feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness.

Exact device longevity cannot be predicted. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, and overall tube care.

For optimal performance, it is recommended that AMT Balloon G-Tube device be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure. If devices are failing or performance is degrading earlier than the typical range for device longevity, it is recommended that you speak with your healthcare professional regarding eliminating common factors that can lead to early device degradation.

NOTE: To help prevent unneeded hospital visits, it is recommended that a spare device is kept on hand at all times for replacement in case device failure occurs prior to scheduled replacement.

THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

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EN	STERILE EO	Sterilized Using Ethylene Oxide	Rx Only Prescription Only
 Do Not Resterilize	 Single Use Only	 Do not use if package is damaged and consult instructions for use	 MR Safe
 Not made with natural rubber latex	MD Medical Device	 Not made with DEHP (Di(2-ethylhexyl) Phthalate)	
ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.			

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ENFit® to zarejestrowany znak towarowy firmy Global Enteral Device Supplier Association, Inc.
ENFit® este o marcă înregistrată a Global Enteral Device Supplier Association, Inc.
ENFit® é marca comercial registrada da Global Enteral Device Supplier Association, Inc.
ENFit® er et registreret varemærke og tilhører Global Enteral Device Supplier Association, Inc.
ENFit® on Global Enteral Device Supplier Association, Inc. registreeritud kaubamärk.
ENFit® on Global Enteral Device Supplier Association, Inc:in rekisteröity tavaramerkki.