
# MINI ONE®

CE  
2862

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






# MiniONE® Non-Balloon Button

## Low Profile Gastrostomy Feeding Tube

### DIRECTIONS FOR USE

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. **Notice:** *Peel off ID sticker from tray and store for future use in patient chart or other convenient location.* The device is provided STERILE for single use. Please inspect all contents of the kit for damage. If damaged, do not use. This device is intended to connect to compatible enteral devices only. Do not use for non-enteral applications.

The MiniONE® Non-Balloon Button should only be used by or under the supervision of personnel trained in percutaneous gastrostomy tube placement. A thorough understanding of the technical principles, clinical applications, and risks associated with percutaneous gastrostomy tube placement is recommended prior to using this device. Please contact your professional health caregiver or physician for explanation of the warnings, care, and use of the device.

### INTENDED USE

The MiniONE® Non-Balloon Button provides nutrition, medication, and decompression access into the stomach through a secured (initial placement) or formed (replacement) stoma. The MiniONE® Non-Balloon Button is intended for use by clinicians and trained patients/caregivers for patients of all age groups (infant through elderly).

### INDICATIONS FOR USE

The Low Profile Non-Balloon Feeding Device is indicated for use in patients who require long term feeding, are unable to tolerate oral feeding, who are at low risk for aspiration, require gastric decompression and/or medication delivered directly into the stomach through a secured (initial placement) or formed (replacement) stoma.

### CONTRAINDICATIONS

Contraindications for placement of the MiniONE® Non-Balloon Button include but are not limited to:

**Initial Placement:** Colonic interposition • Ascites • Portal hypertension • Peritonitis • Uncorrected coagulopathy • Infection around stoma site • Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness)

**Replacement:** Lack of adherence of the stomach to the abdominal wall • Lack of a well-established gastrostomy site • Evidence of infection • Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness) • Presence of multiple stoma fistulous tracts

### COMPLICATIONS

Potential complications when using the MiniONE® Non-Balloon Button include but are not limited to:

Nausea, vomiting, abdominal bloating or diarrhea • Aspiration • Peristomal pain • Abscess, wound infection and skin breakdown • Pressure necrosis • Hypergranulation tissue • Intraperitoneal leakage • Buried bumper syndrome • Peristomal leakage • Tube clog • Gastrointestinal bleeding and/or ulcerations • Ileus or gastroparesis • Bowel and gastric volvulus • Peritonitis • Gastrocolic fistula • sepsis

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

**WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, REPROCESS OR RE-STERILIZE THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.**

### CLINICAL BENEFITS, PERFORMANCE CHARACTERISTICS

Clinical Benefits to be expected when using the MiniONE® Non-Balloon Button include but are not limited to:

- Provides direct access to the stomach for feeding, decompression, and channeling medication
- Device may be placed either during an initial placement procedure or as a replacement device
- Ideal for patients who do not tolerate balloon buttons due to anatomy or differences in gastric environment
- Internal bolster occupies less gastric space and is less likely to block the pylorus
- Internal silicone retention bolster may last longer than balloon devices
- Device comes in a variety of sizes for patient needs
- Made from medical-grade silicone to reduce irritation and improve patient comfort

Performance Characteristics of the MiniONE® Non-Balloon Button include but are not limited to:

- Low profile design sits flush against the skin
- Internal silicone retention bolster holds the device in place and prevents the button from being pulled out of the stomach
- External bolster helps maintain tube position
- Anti-reflux valve prevents backflow of stomach contents
- Glow green feed port provides better visibility in the dark for easier nighttime feeds
- Interlock allows safe and secure attachment of the feeding set
- Insertion/Removal tool distends internal bolster to a smaller profile

### DEVICE MATERIALS

Materials in the MiniONE® Non-Balloon Button include: Medical-grade silicone (92%) • Medical-grade thermoplastic (7%) • Medical-grade silicone pad print ink (1%)

## KIT CONTENTS

| Qty | Component                          | Component Usage Details   |
|-----|------------------------------------|---|
| 1   | MiniONE® Non-Balloon Button        | Provide access into stomach for feeding, medication, decompression.     |
| 2   | 4 x 4 inch Gauze                   | Clean/dry stoma area during or after placement.                         |
| 2   | Feeding Set                        | For channeling feed / medication into button or decompression.          |
| 1   | Large Syringe (Catheter or Enfit®) | For channeling feed / medication / decompression with feeding set.      |
| 1   | Small Syringe (Luer Slip)          | For channeling feed/medication.   |
| 1   | Snap Arm Assembly Tool             | Used to elongate dome with T-Handle for placement or removal of device. |
| 1   | Placement T-Handle                 | Used with Snap Arm Assembly to place the device.                        |
| 1   | Removal T-Handle                   | Used with Snap Arm Assembly to remove the device.                       |

**NOTE:** Kit components listed above may vary slightly depending on kit configuration ordered. Additional or different feeding sets, syringes, or transitional adapters may be present for some kit types.

## TYPE OF USAGE – INITIAL PLACEMENT VS REPLACEMENT

The MiniONE® Non-Balloon Button may be placed either during an initial placement procedure or as a replacement device.

### **PERFORMING AN INITIAL PLACEMENT**

If the patient does not currently have a stoma tract for placement of the MiniONE® Non-Balloon Button device, a new stoma tract will need to be created. This process can only be completed by a healthcare professional. Follow the instructions below starting with the **SUGGESTED INITIAL PLACEMENT PROCEDURES** section in order to properly locate a proper stoma location, secure the stomach to the abdominal wall, create a stoma tract, dilate and measure the stoma, and place the device.

### **REPLACING A FEEDING DEVICE IN AN ESTABLISHED STOMA SITE**

If it is determined that a device placed in an existing stoma needs replaced (due to non-optimal functionality or as a pro-active replacement), the current device can be removed from the stoma and the new device can be inserted in the same stoma site. If performing a device replacement, skip straight to the **DEVICE PLACEMENT PROCEDURE** section for the proper method of placing the new device.

**CAUTION:** Selection of the correct size feeding device is critical for safety and comfort of the patient. An inappropriately sized device can cause necrosis, buried bumper syndrome and/or hypergranulation tissue. If patient size/weight has changed since the device was placed, fit is too tight or too loose, or it has been over six months since the device was last measured, it is recommended that you speak to your healthcare professional to ensure that device size does not need changed.

## SUGGESTED INITIAL PLACEMENT PROCEDURES

Initial placement of the MiniONE® Non-Balloon Button may be performed percutaneously under laparoscopic, fluoroscopic, or endoscopic guidance. Follow the instructions below for the operator's preferred method of placement.

**NOTE:** The insertion site for infants and children should be high on the greater curvature to prevent occlusion of the pylorus.

### **OPTION A: RADIOLOGIC PLACEMENT**

1. Place the patient in the supine position.
2. Prep and sedate the patient according to clinical protocol.
3. Ensure that the left lobe of the liver is not over the fundus or the body of the stomach.
4. Identify the medial edge of the liver by CT scan or ultrasound.
5. Glucagon 0.5 to 1.0 mg IV may be administered to diminish gastric peristalsis.

**CAUTION:** Consult Glucagon instructions for rate of IV injection and recommendations for use with insulin dependent patients.

6. Insufflate the stomach with air using a nasogastric catheter, usually 500 to 1,000 ml or until adequate distention is achieved. It is often necessary to continue air insufflation during the procedure, especially at the time of needle puncture and tract dilation, to keep the stomach distended so as to appose the gastric wall against the anterior abdominal wall.
7. Choose a catheter insertion site in the left sub-costal region, preferably over the lateral aspect or lateral to the rectus abdominis muscle (Note: The superior epigastric artery courses along the medial aspect of the rectus) and directly over the body of the stomach toward the greater curvature. Using fluoroscopy, choose a location that allows as direct a vertical needle path as possible. Obtain a cross table lateral view prior to placement when interposed colon or small bowel anterior to the stomach is suspected.

**NOTE:** PO/NG contrast may be administered the night prior or an enema administered prior to placement to pacify the transverse colon.

8. Prep and drape according to facility protocol.
9. Follow instructions in the **GASTROPEXY PLACEMENT** section.

### **OPTION B: ENDOSCOPIC PLACEMENT**

1. Perform routine Esophagogastroduodenoscopy (EGD). Once the procedure is complete and no abnormalities are identified that could pose a contraindication to placement of the tube, place the patient in the supine position and insufflate the stomach with air.
2. Transilluminate through the anterior abdominal wall to select a gastrostomy site that is free of major vessels, viscera and scar tissue. The site is usually one third the distance from the umbilicus to the left costal margin in the midclavicular line.
3. Depress the intended insertion site with a finger. The endoscopist should clearly see the resulting depression on the anterior surface of the gastric wall.
4. Prep and drape the skin at the selected insertion site.
5. Follow instructions in the **GASTROPEXY PLACEMENT** section.

## GASTROPEXY PLACEMENT

**WARNING: INITIAL PLACEMENT OF THE MINIONE® NON-BALLOON BUTTON REQUIRES THAT A GASTROPEXY IS PERFORMED TO AFFIX THE STOMACH WALL TO THE ANTERIOR ABDOMINAL WALL. DO NOT USE THE INTERNAL DOME OF THE DEVICE AS A GASTROPEXY DEVICE. AN EARLY DEVICE FAILURE MAY PREVENT THE STOMACH FROM ATTACHING TO THE ANTERIOR ABDOMINAL WALL.**

**CAUTION:** It is recommended to perform a three point gastropexy in a triangle configuration to ensure attachment of the gastric wall to the anterior abdominal wall.

1. Place a skin mark at the tube insertion site. Define the gastropexy pattern by placing three skin marks equidistant from the tube insertion site and in a triangle configuration.

**WARNING: ALLOW ADEQUATE DISTANCE BETWEEN THE INSERTION SITE AND GASTROPEXY PLACEMENT TO PREVENT INTERFERENCE OF THE T-FASTENER AND INTERNAL DOME.**

## GASTROPEXY PLACEMENT

2. Localize the puncture sites with 1% lidocaine and administer local anesthesia to the skin and peritoneum.
3. Place the first T-Fastener and confirm intragastric position. Repeat the procedure until all three T-Fasteners are inserted at the corners of the triangle.
4. Secure the stomach to the anterior abdominal wall and complete the procedure.
5. Follow instructions in the **CREATE THE STOMA TRACT** section.

**CAUTION: After the MiniONE® Non-Balloon Button is placed, do not remove the gastropexy sutures until the stoma site fully heals and the stomach is fully affixed to the anterior abdominal wall.**

## CREATE THE STOMA TRACT

1. Create the stoma tract with the stomach still insufflated and in apposition to the abdominal wall. Identify the puncture site at the center of the gastropexy pattern. With fluoroscopic guidance confirm that the site overlies the distal body of the stomach below the costal margin and above the transverse colon.

**CAUTION: Avoid the epigastric artery that courses at the junction of the medial two-thirds and lateral one-third of the rectus muscle.**

**WARNING: TAKE CARE NOT TO ADVANCE THE PUNCTURE NEEDLE TOO DEEPLY IN ORDER TO AVOID PUNCTURING THE POSTERIOR GASTRIC WALL, PANCREAS, LEFT KIDNEY, AORTA OR SPLEEN.**

2. Anesthetize the puncture site with local injection of 1% lidocaine down to the peritoneal surface.
3. Insert a .038" compatible introducer needle at the center of the gastropexy pattern into the gastric lumen directed toward the pylorus.

**NOTE:** The best angle of insertion is a 45 degree angle to the surface of the skin.

### **OPTION A. RADIOLOGIC PLACEMENT ONLY (4A & 5A):**

- 4A. Use fluoroscopic visualization to verify correct needle placement. Additionally, to aid in verification, a water filled syringe may be attached to the needle hub and air aspirated from the gastric lumen.

**NOTE:** Contrast may be injected upon return of air to visualize gastric folds and confirm position.

- 5A. Advance a guidewire, up to .038", through the needle and into the stomach. Confirm position.

### **OPTION B. ENDOSCOPIC PLACEMENT ONLY (4B & 5B):**

- 4B. Use endoscopic visualization to verify correct needle placement.
- 5B. Advance a guidewire, up to .038", through the needle into the stomach. Using endoscopic visualization, grasp the guidewire with atraumatic forceps.
6. Remove the introducer needle, leaving the guidewire in place and dispose of according to facility protocol.
7. Follow instructions in the **DILATION** section.

## DILATION

1. Use a #11 scalpel blade to create a small skin incision that extends alongside the guidewire, downward through the subcutaneous tissue and fascia of the abdominal musculature. After the incision is made, dispose of according to facility protocol.
2. Advance a dilator over the guidewire and dilate the stoma tract to the desired size.
3. Remove the dilator over the guidewire, leaving the guidewire in place.
4. Measure the Stoma Length with the AMT Stoma Measuring Device by following the **MEASURE THE STOMA LENGTH** section.

**CAUTION: Avoid excessive dilation of stoma tract as this can increase the possibility of the dome pulling through the stoma.**

## MEASURE THE STOMA LENGTH

**CAUTION: Selection of the correct size feeding device is critical for the safety and comfort of the patient. Measure the length of the patient's stoma with the AMT Balloon Stoma Measuring Device. The shaft length of the selected feeding device should be the same as the length of the stoma. An inappropriately sized feeding device can cause necrosis, buried bumper syndrome, and/or hypergranulation tissue.**

1. Please refer to Manufacturer's Directions for Use for Stoma Measuring Device being used. It is recommended that AMT's Balloon Stoma Measurement Device be used for measurement of the stoma.
2. Be sure to select the appropriate size MiniONE® Non-Balloon Button for the abdominal wall thickness measured. If the measurement appears to be between two sizes, always select the next larger size MiniONE® Non-Balloon Button. Once placed, the external flange should rotate easily.

**WARNING: UNDER-SIZING THE DEVICE MAY CAUSE EMBEDDING WITH EROSION INTO THE GASTRIC WALL, TISSUE NECROSIS, INFECTION, SEPSIS AND ASSOCIATED SEQUELAE.**

## DEVICE PLACEMENT PROCEDURE

**WARNING: PRIOR TO PLACEMENT, PLEASE INSPECT ALL CONTENTS OF THE KIT FOR DAMAGE. IF THE PACKAGE IS DAMAGED OR STERILE BARRIER IS BREACHED, DO NOT USE THE PRODUCT.**

At the time of placement, test the direction of the stoma using a soft catheter of the appropriate dimension (e.g., a nasogastric tube). *Insertion at an improper angle may disrupt the stoma tract, possibly causing stomach separation from the abdominal wall, which can lead to tissue damage, infection, sepsis and/or peritonitis.*

1. Remove the existing gastrostomy tube from the stoma using the original manufacturer's instructions. (For MiniONE® Non-Balloon removal see "**Removal of the MiniONE® Non-Balloon Button**".)
2. Insert T-Handle into the MiniONE® Non-Balloon Button by lining up the key with the key opening along the button centerline (see **Fig A1**). Insert the T-Handle fully and turn clockwise 1/2 turn (see **Fig A2**).
3. Lubricate the rod of the Snap Arm, and then slowly slide the Snap Arm through the T-Handle & into the MiniONE® Non-Balloon Button elongating it (see **Fig A3**). Verify rod tip is centered in the button dome.
4. Ratchet the Snap Arm down all the way to the bottom of the T-Handle until it stops.
5. Lubricate the tip of the elongated MiniONE® Non-Balloon Button and the stoma site with a water soluble lubricant.
6. Holding the T-Handle, gently guide the device through the stoma and into the stomach until the external flange is flush to the skin (see **Fig A4**).
7. Release and remove the Snap Arm by briefly pushing it in (while holding the T-Handle)—then squeezing the triangular ends of the Snap Arm and sliding it out (see **Fig A5**).
8. After removing the Snap Arm, disengage T-Handle by carefully rotating counterclockwise until it stops and remove it (see **Fig A6**).
9. Verify that the external flange of the MiniONE® Non-Balloon Button rotates freely.

**NOTE:** After the MiniONE® Non-Balloon Button has been placed, stomach contents may initially leak out of the anti-reflux valve in the feeding port area. This is a normal occurrence and occasionally happens due to the anti-reflux valve being activated by the obturator rod. To stop stomach contents from leaking, reactivate the anti-reflux valve by attaching a feeding set adapter (or the T-Handle). Once the valve has been reactivated, remove the feeding adapter (or T-Handle) and check again for stomach content leakage.

## PLACEMENT VERIFICATION

1. Attach the feeding set to the MiniONE<sup>®</sup> Non-Balloon Button by lining up the dark line on the feeding set connector with the dark line on the interlock of the MiniONE<sup>®</sup> Non-Balloon Button. Press feeding set connector to insert into the MiniONE<sup>®</sup> Non-Balloon Button. Turn 3/4 to the right (clock-wise) to lock the feeding set into place.
2. Attach a catheter tip syringe with 5ml of water into the feeding set. Aspirate for contents. Spontaneous return of gastric contents should occur.
3. After air and/or gastric contents are observed, flush the device with water.

**WARNING: NEVER INJECT AIR INTO THE MINIONE<sup>®</sup> NON-BALLOON BUTTON.**

4. When placement is confirmed, gastric decompression or feeding can begin.

### **ADDITIONAL OPTIONS FOR DETERMINATION OF BUTTON LOCATION:**

1. Radiographic examination (x-ray) may be performed to evaluate the location of the Button's radiopaque dome.
2. A fiberoptic exam can also be done, by using a 2.5mm or smaller endoscope inserted down through the lumen of the button.
3. A flexible guidewire can also be used in this manner and a radiographic examination can then be performed. Do NOT force guidewire.
4. Aspiration may be done to check pH for acidic fluid. Tubing should be monitored for possible migration or unintentional extubation.

## PLACEMENT CARE INSTRUCTIONS

1. The MiniONE<sup>®</sup> Non-Balloon Button should be replaced periodically for optimal performance. If in place for three months or longer, frequent performance checks are recommended. Clogging and/or reduced flow are indications of diminished performance. Refer to the **TROUBLESHOOTING** section for additional signs of reduced performance or failure.
2. Be sure the internal bolster is within the stomach and the MiniONE<sup>®</sup> Non-Balloon Button rotates freely **before the feedings begin**. Inject water or saline into the feeding port to determine proper function of the valve.
3. The stoma area should be cleansed daily with a mild soap and water. The stoma site should be clean and dry at all times. The MiniONE<sup>®</sup> Non-Balloon Button should be rotated daily for site hygiene.

## FEEDING SET INSTRUCTIONS FOR USE

**WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO ENTERAL DEVICES. DO NOT USE FOR NON-ENTERAL APPLICATIONS.**

**WARNING: IF USING A LEGACY-STYLE (NON-ENFIT<sup>®</sup>) EXTENSION SET, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO THE FOLLOWING SYSTEMS: BREATHING, CONES AND SOCKETS OF ANESTHETIC AND RESPIRATORY EQUIPMENT, INTRAVENOUS, LIMB CUFF, NEURAXIAL CONNECTORS, NIPPLES OF RESPIRATORY THERAPY EQUIPMENT, URINARY, AND TEMPERATURE SENSOR CONNECTORS OF RESPIRATORY HUMIDIFYING EQUIPMENT.**

1. Inspect contents for damage. If damaged, do not use. Obtain another package.
2. The feeding set may be used for syringe, gravity, pump feeding or decompression.
3. Ensure the feeding set clamp is closed and attach the feeding set to the button by lining up the dark line on the feeding set connector with the dark line on the button. Fully press the feeding set connector into the button. Turn 3/4 to the right (clock-wise) to lock the feeding set in place.
4. Attach the opposite end of the feeding set to the connector being used. If using a bolus or luer connector, firmly insert the connector while using a slight rotation into the feeding set to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the feeding set while avoiding excessive force to lock in place. Once connected, open the clamp to allow flow.

**WARNING: ONLY TIGHTEN BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.**

**WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.**

**WARNING: IF FEEDING SET IS NOT PROPERLY PLACED AND LOCKED, LEAKAGE MAY OCCUR. WHEN TWISTING THE CONNECTOR HOLD THE BUTTON IN PLACE TO AVOID ROTATION WHILE PLACING THE FEEDING SET.**

5. When feeding is complete, flush with 5-10ml water. Close the clamp to prevent leakage while removing the feeding set. Remove the feeding set by holding the button in place, grasp the feeding connector, twist 3/4 turn to the left (counter clock-wise) to unlock, and gently remove the feeding set from the button. Snap button plug in place to keep lumen clean.
6. The feeding set should be cleansed with soap and water and thoroughly rinsed. Store in a clean, dry place. Never use a dishwasher for cleaning.

## MEDICATION CHANNELING

Liquid medication is preferred when available. If solid medication is required, consult with your physician regarding whether it is safe to crush the medication. If safe, medication should be crushed as fine as possible (into powder form) and dissolved into water before channeling medication through device. Channeling solid medication through the tubing that is not properly crushed can lead to blockage in the tubing. Never crush enteric coated medication or mix medication with formula. Flush device with water after channeling medication.

## DECOMPRESSION

Decompression can be performed if there is a need to check for residuals, or if the patient tends to retain air during feeding causing bloating and discomfort. To decompress, connect the feeding set to the interlock connector and drain stomach contents into a container. After decompression is complete, flush the feeding set with water and remove the feeding set from the device.

## PROPER FLUSHING GUIDELINES

Proper flushing techniques and frequency can help prevent tubing blockage, clogs, and tubing failure. Follow these flushing guidelines to maintain optimal device flow conditions:

- Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain potency.
- Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube before and after medication channeling and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

## PROPER FLUSHING GUIDELINES

### UNCLOGGING A DEVICE:

First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a catheter tip syringe to a feeding set and attach into the interlock connector. Fill the syringe with warm water and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced.

**CAUTION: Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.**

## REMOVAL OF THE MINIONE® NON-BALLOON BUTTON

**WARNING: SPONTANEOUS CLOSURE OF STOMA MAY OCCUR WITHIN TWENTY-FOUR (24) HOURS AFTER REMOVAL. INSERT A NEW DEVICE IF ENTERAL FEEDING BY THIS ROUTE IS STILL INTENDED. IF CLOSURE IS DESIRED, APPLY A DRESSING OVER THE STOMA SITE. REMOVAL SHOULD NOT TAKE PLACE UNTIL THE STOMA IS FULLY ESTABLISHED, WHICH MAY TAKE TWO MONTHS OR MORE. REMOVAL IS RECOMMENDED TO BE PERFORMED BY A QUALIFIED CLINICIAN.**

1. Lubricate the stoma site underneath the existing MiniONE® Non-Balloon Button with a water soluble lubricant. (A topical anesthetic gel may also be applied to the stoma and surrounding skin.)
2. Using the T-Handle with metal reinforcer, first rotate reinforcer up away from the tip of the T-Handle (see **Fig B1**).
3. Insert T-Handle into the MiniONE® Non-Balloon Button by lining up the key with the key opening along the button centerline (see **Fig B2**). Insert the T-Handle fully and turn clockwise 1/2 turn (see **Fig B3**).
4. With the T-Handle in the 1/2 turn position, rotate the metal reinforcer down—sliding it underneath the bolster of the MiniONE® Non-Balloon Button (see **Fig B4**).
5. Lubricate the rod of the Snap Arm, and then slowly slide the Snap Arm through the T-Handle and into the MiniONE® Non-Balloon Button elongating it (see **Fig B5**).

**NOTE:** Resistance should be observed when ratcheting the Snap Arm down. If no resistance is observed when ratcheting down the Snap Arm, remove the Snap Arm & T-Handle and remove the device via traction.

6. Ratchet the Snap Arm all the way down to the bottom of the T-Handle until it stops against the metal reinforcer arm (see **Fig B6**).
7. Gently remove the device from the stoma site (see **Fig B7**).
8. Replace as necessary following the new device's Directions for Use. (For Non-Balloon placement instructions—see previous instructions.)

### Additional Methods for Removal:

- A. The MiniONE® Non-Balloon Button can also be traction removed without elongating the device.
  - i. Simply lubricate the stoma site underneath the existing MiniONE® Non-Balloon Button with a water soluble lubricant. (A topical anesthetic gel may also be applied to the stoma and surrounding skin.)
  - ii. Grasp the tabs, apply constant traction and carefully remove the device.
- B. Endoscopic removal is another option.
  - i. Insert endoscope, cut shaft and retrieve dome.
  - ii. Do NOT let the dome pass through the intestinal tract.

**WARNING: IF THE TUBE IS RESISTANT TO REMOVAL, LUBRICATE THE STOMA SITE WITH A WATER-SOLUBLE LUBRICANT. ROTATE THE TUBE GENTLY AND PUSH IT IN ABOUT 2.5 CM. CAREFULLY WORK THE TUBE FREE. NEVER USE FORCE TO REMOVE THE TUBE.**

- C. If traction provides too much stress to the patient, the catheter can be cut at skin level and the internal bolster removed endoscopically.

**WARNING: DO NOT ALLOW THE INTERNAL BOLSTER (CUT CATHETER SECTION) TO PASS THROUGH THE INTESTINAL TRACT.**

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

## DEVICE LONGEVITY

### MiniONE® Non-Balloon Button:

Low-profile feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness.

Exact device longevity cannot be predicted. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 3-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, trauma to the device, contact with sharp or abrasive objects, incorrect stoma length measurement, and overall tube care.

For optimal performance, it is recommended that MiniONE® Non-Balloon Button device be changed every 6 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure. If devices are failing or performance is degrading earlier than the typical range for device longevity, it is recommended that you speak with your healthcare professional regarding eliminating common factors that can lead to early device degradation. Also reference the **TROUBLESHOOTING** section for additional information regarding device performance issues.

**NOTE:** To help prevent unneeded hospital visits, it is recommended that a spare device is kept on hand at all times for replacement in case device failure occurs prior to scheduled replacement.

### Feeding Sets:

Feeding sets are meant to be periodically replaced for optimal performance and cleanliness. AMT recommends that the feeding set be changed at least every 2 weeks or as often as indicated by your healthcare professional. Device performance and functionality can degrade over time depending on usage and environmental conditions. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, frequency of use, trauma to the device, and overall feeding set care.

Device should be replaced if leakage, cracks, tears, excessive residue build-up, mold, fungus, or other signs of failure are noted. Some usage types can wear down device components quicker than others. If tubing is becoming hardened, replacement is recommended in near future to avoid failure.

## MRI SAFETY INFORMATION

The MiniONE® Non-Balloon Button is considered MR Safe once placed in the stoma site when none of the MR Unsafe components are attached.

### Components:

The MiniONE® Non-Balloon Button, Placement T-Handle, feeding sets, syringes, and gauze are **MR Safe**.

The Snap Arm Assembly Tool and Removal T-Handle are **MR Unsafe**.

## STORAGE

Store all device components in a dry, clean, controlled room temperature environment while not in use.

## TROUBLESHOOTING

Long-term device performance and functionality are dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your feeding device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

**Anti-reflux valve leakage or blockage:** Leakage/blockage of the anti-reflux valve typically occurs due to residue (feed, medication, gastric contents, etc) becoming stuck in the valve area, preventing the valve from fully closing. Make sure device is flushed after each use. Leakage can also be caused due to excessive pressure in stomach (see **DECOMPRESSION** section). The valve can also become inverted in rare cases. Insert feeding set into port to reset valve if this occurs. Avoid excessive pressure during channeling through the device, as this can cause a tear to occur in the valve over time.

**A tear has formed:** Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.

**Tubing has reduced flow or has become clogged:** Tubing can become blocked due to not properly flushing after each use, use of thick or improperly crushed medications, use of thick feeds/formulas, gastric reflux, and/or fungus growth. If clogged, refer to the **PROPER FLUSHING GUIDELINES** section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced.

**Plug will not stay closed:** Ensure that the plug is being firmly and fully pressed into the interlock connector. If the plug is not staying closed, check the plug and feed-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.

**Device has pulled from stoma:** The device may be undersized or excessive force was applied to the device during use. Avoid excessive force while using the device and secure feeding set if it is a source of tension. If device length has not been checked recently, resizing the stoma length may be necessary. Device will need to be placed by a medical professional promptly to avoid closure of the stoma site.

**Interlock failure or cracked:** The interlock has been designed to withstand extreme forces without detaching or cracking. However, the strength of the bond and material can reduce over prolonged use depending on medications and feeds used through the device. Excessive forces over an extended period of use can also reduce bond / material strength. The device should be replaced if interlock is found cracked, leaking, or separating from the device.

**Device fit is too tight or too loose:** A device that is not correctly fitted for the stoma site can lead to leakage, embedding with erosion into the gastric wall, tissue necrosis, infection, sepsis and associated sequelae, or device dislodgement. If the device is not fitting properly, it is recommended that the stoma site is measured to ensure the correct device length is being used.

**Foul smell coming from the device:** Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed and stoma site should be gently cleaned with soap and warm water. If foul smell does not go away, it is recommended that you contact your healthcare professional.

**Device has become discolored:** The device can become discolored over days to months of use. This is normal depending on the types of feeds and medications being used with the device.

**Device difficult to remove:** Ensure that the correct removal tool length is being used for removal. If removal of the device is not possible using the preferred removal method, traction or endoscopic removal of the device may be necessary. Follow the instructions in the **REMOVAL OF THE MINIONE® NON-BALLOON BUTTON** section. If dome is cut during the removal, **DO NOT** allow the dome to pass through the intestinal tract.

## THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

|   |   |   |   |                   |
|---|---|---|---|-------------------|
| EN  | <b>STERILE EO</b>   | Sterilized Using Ethylene Oxide   | Rx Only   | Prescription Only |
|  Do Not Resterilize                 |  Single Use Only |  Do not use if package is damaged and consult instructions for use |  MR Safe |                   |
|  Not made with natural rubber latex |  Medical Device  |  Not made with DEHP (Di(2-ethylhexyl) Phthalate)                   |   |                   |
| ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.                                  | MIC-KEY® feedset compatible   |   |   |                   |
|   | MIC-KEY® is a registered trademark of Avent, Inc.   |   |   |                   |
|   | MiniONE® is a registered trademark of Applied Medical Technology, Inc.                              |   |   |                   |

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