AMT Initial Placement Dilator Set

CE
2862

- ES AMT Colocación Inicial Conjunto Dilatador
- DE AMT Erstplatzierung Dilatorsatz
- IT AMT Posizionamento Iniziale set di dilatatori
- FR AMT Dilatateur Placement Initial Fixé
- NO AMT Innledende Plassering Utvidersettet
- SV AMT Ursprungliga Placering Dilatatorn Som
- AR مجموعة موسعات الموضعة الابتدائية
- JA AMT初期配置ダイレータセット
- PL Zestaw do demontażu początkowego AMT
- PT-BR Conjunto de dilatador para colocação inicial AMT
- RO Set de dilatare pentru plasarea inițială AMT
- DA AMT-dilatatorsæt til indledende placering
- FI AMT ensiasetus laajennin -pakkaus
- EL AMT Σετ αρχικής τοποθέτησης διαστολέα
- ET AMT dilaatorikomplekt esmapaigalduseks
- PT Conjunto Dilatador de Colocação Inicial de AMT
- NL AMT Dilator Set voor initiële plaatsing
- BG Удължителен комплект за първоначално поставяне AMT
- HU AMT Kezdeti elhelyezés tágitókészlet
- LV AMT sākotnējās ievietošanas paplašinātāja komplekts
- LT AMT pirminio įdėjimo plėtiklio rinkinys
- SK Súprava dilatátora na počiatocné zavedenie AMT
- SL Komplet dilatatorjev za prvo vstavitev AMT
- TR AMT İlk Yerleşirme Dilatör Seti
- CS Sada dilatátorů k úvodnímu umístění AMT

R_x ONLY
STERILE EO

MD

KIT CONTENTS (Compatible with 12-18F Devices)

(1) Guidewire	(4) Dilators: 7F, 12F, 16F, 20F
(1) Guidewire Compatible Introducer/Stiffener	(2) Introducer Needle, 18 gauge

KIT CONTENTS (Compatible with 10F Devices)

(1) Guidewire	(3) Dilators: 7F, 12F, 14F
(1) Guidewire Compatible Introducer/Stiffener	(2) Introducer Needle, 18 gauge

SUPPLEMENTAL COMPONENTS (Not included in the kit)

Endoscope	Sterile Water or Saline
Local Anesthetic	Water Soluble Lubricant
Scalpel	U-Stitches (or alternate means of securement)
Hemostats	

INDICATIONS FOR USE

CAUTION: Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

INTENDED USE

The AMT Initial Placement Dilator Set is an effective tool to incrementally increase the size of a gastric stoma. The AMT Initial Placement Dilator Set is intended to be used by clinicians for patients of all ages (infant to elderly).

INDICATIONS

Initial placement of a percutaneous gastrostomy device may be indicated for patients with a functioning gut who require long-term feeding. This includes patients in whom malnutrition already exists, or may result, secondary to concurrent conditions.

CONTRAINDICATIONS

Contraindications for placement of a gastrostomy tube include: Nonfunctioning gut • Infection on or around the GI site • Inability to identify transabdominal illumination • Multiple surgical procedures near the GI site

WARNING: AFTER GI DEVICE PLACEMENT, PROPER POSITIONING OF THE BALLOON OR NON-BALLOON TYPE INTERNAL BOLSTER AGAINST THE GASTRIC MUCOSA MUST BE ENDOSCOPICALLY VERIFIED. TENSION ON THE DEVICE SHOULD BE AVOIDED IN ORDER TO MINIMIZE THE RISK OF COMPLICATIONS. DEVICE SHOULD BE CAPABLE OF ROTATING WITHOUT DIFFICULTY.

COMPLICATIONS

At this time there are no complications associated with the use of the AMT Initial Placement Dilator Set.

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

CLINICAL BENEFITS

Clinical benefits to be expected when using the AMT Initial Placement Dilator Set include but are not limited to: Patients will benefit from the proper primary placement of a gastrostomy tube allowing for optimal channeling of nutritional and medication • Reduces potential trauma associated with stoma tract development • Clinicians will benefit by having most of the necessary tools to properly and easily place a gastrostomy tube during initial placement surgeries • Guidewire compatibility helps to ensure the ease of proper placement of low profile feeding devices

PERFORMANCE CHARACTERISTICS

Performance characteristics of the AMT Initial Placement Dilator Set include but are not limited to: Provides the dilators needed to surgically create a stoma site • Available in different configurations to accommodate different stoma/device sizes • Provides a means of stoma tract dilation to accommodate an interventional catheter and a means of introducing an interventional catheter to a freshly formed stoma tract • Tapered dilators to help ensure proper dilation

PROCEDURE PREPARATION


NOTE: Initial Placement is to be performed by a qualified clinician.

- Inspect the kit for damage. If damaged or sterile barrier is breached, do not use product.
- Use a clinically approved method to prep and sedate the patient for endoscopic procedure.

SELECTING AND PREPARING THE GASTROSTOMY SITE:

- Once the patient is in the supine position and the endoscope is in place within the stomach, insufflate the stomach with air.
- Dim the room lights and deflect the endoscope to the anterior surface. The endoscope's light should clearly transilluminate through the anterior abdominal wall.
- Choose an insertion site that is free of major vessels, viscera, and scar tissue.

CAUTION: Attaining definitive transillumination at the insertion site is one of the most critical steps in the procedure. Do not continue if there is any question concerning an insertion site free of major vessels, viscera, and scar tissue.

- Depress the illuminated, intended insertion site. The depression on the anterior surface of the gastric wall should be clearly visible endoscopically.

CAUTION: Selection of the proper insertion site is a critical step in the procedure. Do not continue if there is inadequate visualization of the depression on the anterior surface of the gastric wall.

- Prep and drape the skin at the selected insertion site.

PLACEMENT PROCEDURE

1. SECUREMENT OF STOMACH TO ANTERIOR ABDOMINAL WALL

- At physician's discretion, an alternate method of securement may be necessary with larger adult patients due to greater thickness of the abdominal wall.
- Using two U-stitches inserted on opposite sides of the insertion site, secure the stomach against the anterior abdominal wall.
- Once the stitches are complete, take care to secure the two ends of each stitch in order to maintain the stomach against the abdominal wall.

PLACEMENT PROCEDURE

2. GUIDEWIRE PLACEMENT

While maintaining insufflation of the stomach and endoscopically visualizing the anterior mucosal surface, insert an 18-gauge needle through the peritoneum and stomach wall.

NOTE: The most desirable insertion angle is perpendicular to the surface of the abdomen and stomach. Other restrictions may apply at the discretion of the physician.

CAUTION: It is imperative that one takes extreme caution when inserting the introducer needle to ensure the posterior wall and other vessels or viscera of the stomach are not perforated.

- When the 18-gauge needle is observed within the stomach, to minimize air loss through the needle, quickly cover the needle hub.
- Pass the enclosed guidewire through the 18-gauge needle.
- Once the guidewire is observed in the stomach via the endoscope, remove the 18-gauge needle over the guidewire, taking caution not to dislodge the guidewire from the stomach.

NOTE: In the subsequent steps requiring use of the guidewire, do not pull up on the guidewire, dislodging it from the stoma.

3. STOMA TRACT CREATION AND STOMA DILATION

- Using a scalpel, create a small skin incision that extends directly alongside the guidewire. If desired, utilize the endoscope to verify the incision extends downwards through the anterior gastric wall into the stomach.
- Using the enclosed dilator set, expand the stoma tract to the necessary diameter for the chosen device size.
- Begin with a smaller size dilator. Lubricate surface of the dilator with a water-soluble lubricant.

CAUTION: It is imperative to take extreme caution when inserting dilators to ensure the posterior wall of the stomach is not perforated.

- Pass the dilator over the guidewire. Firmly but gently, use a rotating motion to twist the dilator into the stoma tract and slightly into the stomach. Utilize endoscope to visually verify dilator passage.
- Remove the first dilator using caution not to dislodge the guidewire. To minimize air loss through the stoma tract, hold a finger over the incision site.
- Repeat using larger dilators as necessary.
- Refer to **Chart 1** for recommended necessary dilation of French size for selected low profile device.

Chart 1				
Product French Size	10F	12F	14F	16-18F
Dilator	14F	16F	18F	22F

4. DEVICE PREPARATION

CAUTION: When using a low profile feed device, do not use an undersized device. Once placed, the device should rotate easily.

WARNING: UNDERSIZING THE DEVICE MAY CAUSE EMBEDDING WITH EROSION IN THE GASTRIC WALL, TISSUE NECROSIS, INFECTION, SEPSIS, AND ASSOCIATED SEQUELAE.

- If utilizing a balloon device, inflate the balloon prior to insertion to verify that it is fully functional. If the balloon is not functional, discard the product and obtain a new product.
- Lubricate the tip portion of the enclosed introducer stiffener using a water-soluble lubricant.
- Insert the introducer stiffener into the device opening until the tip is flush with the tip of the device and necessary stiffening is achieved.
- Lubricate the outer surface of the device using a water-soluble lubricant.

5. DEVICE PLACEMENT

- Pass the introducer stiffener and device over the guidewire.
- Insert device through the stoma and into the stomach. If necessary, use a firm but gentle rotating motion to insert.
- If utilizing a balloon device, inflate the balloon through the balloon fill valve using the recommend amount of sterile water or saline and a slip tip syringe. Verify endoscopically that the entire internal bolster portion of the device is within the stomach.
- For non-guidewire compatible gastrostomy devices, the guidewire may need to be removed prior to placement. Follow gastrostomy device directions for use for method of insertion and internal bolster deployment.

WARNING: IF USING A BALLOON PRODUCT, ENDOSCOPICALLY GRIPPING OR GRASPING THE BALLOON OR TIP OF THE DEVICE MAY CAUSE DEVICE FAILURE.

- If utilizing a non-balloon device, release the bolster per manufacturing instruction.
- Verify endoscopically that internal bolster is retaining device against anterior surface of the mucosal wall.
- Rotate the device to ensure that it has a snug but not overly tight fit.

CAUTION: The inflated balloon or internal bolster should be capable of rotating freely. There should be no compression of either the gastric mucosa or the skin.


- Slowly withdraw the introducer stiffener and guidewire from the device.
- CAUTION:** Care must be taken when withdrawing the introducer stiffener to ensure that the entire device is not also removed.

Tie the sutures over the device.
 Deflate the stomach and remove the endoscope.
WARNING: KIT CONTENTS ARE INTENDED FOR SINGLE USE. DO NOT REUSE, RE-STERILIZE OR REPROCESS THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.

NOTE: The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

POSTOPERATIVE PROCEDURE AND PRECAUTIONS:

- It is advisable to wait until the following morning before initiating feedings. First feeding time, prescribed rate and volume are at physician's discretion.

STERILE EO Sterilized Using Ethylene Oxide		Rx Only Prescription Only	
Do Not Resterilize	Single Use Only	Do not use if package is damaged and consult instructions for use	
Not made with natural rubber latex	Medical Device	Not made with DEHP (Di(2-ethylhexyl) Phthalate)	

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Applied Medical Technology, Inc.