# MINI ONE®

## BALLOON BUTTON

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**Rajah 4:** Port Pengembungan Belon

# MiniONE® Balloon Button

## Low Profile Gastrostomy Feeding Tube

### INSTRUCTIONS FOR USE

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**Notice:** Peel off ID sticker from tray and store for future use in patient chart or other convenient location.

The device is provided STERILE for single use. Please inspect all contents of the kit for damage. If damaged, do not use. This device is intended to connect to compatible enteral devices only. Do not use for non-enteral applications.

### INTENDED USE

The MiniONE® Balloon Button provides nutrition, medication, and decompression access into the stomach through a secured (initial placement) or formed (replacement) stoma. The MiniONE® Balloon Button is intended to be used by clinicians and trained patients/caregivers.

### INDICATIONS FOR USE

The Low Profile Balloon Feeding Device is indicated for use in patients who require long term feeding, are unable to tolerate oral feeding, who are at low risk for aspiration, require gastric decompression and/or medication delivered directly into the stomach through a secured (initial placement) or formed (replacement) stoma. The Low Profile Balloon Feeding Device is intended for all age groups.

### CONTRAINDICATIONS

Contraindications for placement of the MiniONE® Balloon Button include but are not limited to:

#### Initial Placement:

Colonic interposition • Ascites • Portal hypertension • Peritonitis • Uncorrected coagulopathy • Infection around stoma site • Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness)

#### Replacement:

Lack of adherence of the stomach to the abdominal wall • Lack of a well established gastrostomy site • Evidence of infection • Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness) • Presence of multiple stoma fistulous tracts

#### Notice:

- Please contact a professional health care giver or physician for explanation of the warnings, care, and use of the device.
- This device has been designed to provide feeding/medication/decompression access to the stomach. Other applications not advised.

### COMPLICATIONS

Potential complications when using the MiniONE® Balloon Button include but are not limited to:

Nausea, vomiting, abdominal bloating or diarrhea • Aspiration • Peristomal pain • Abscess, wound infection and skin breakdown • Pressure necrosis • Hypergranulation tissue • Intraoperative leakage • Buried bumper syndrome • Peristomal leakage • Balloon failure or dislodgement • Tube clog • Gastrointestinal bleeding and/or ulcerations • Ileus or gastroparesis • Bowel and gastric volvulus

**WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, REPROCESS OR RE-STERILIZE THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.**

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

### CLINICAL BENEFITS, PERFORMANCE CHARACTERISTICS

Clinical Benefits to be expected when using the MiniONE® Balloon Button include but are not limited to:

- Provides direct access to the stomach for feeding, decompression, and channeling medication
- Device may be placed either during an initial placement procedure or as a replacement device
- Balloon design aids in easy, at home replacement when necessary
- Soft, flexible bolster is easy to clean and less likely to create trauma to the stoma site
- Balloon design results in a larger gastric seal surface which may help to reduce leakage and granulation tissue
- Made from medical-grade silicone to reduce irritation and improve patient comfort

Performance Characteristics of the MiniONE® Balloon Button include but are not limited to:

- Low profile design sits flush against the skin
- Anti-leakage valve prevents backflow of stomach contents
- Balloon fill volume can be easily adjusted for a custom-fit
- Interlock allows safe and secure attachment of the extension set
- Unique "apple" shape balloon designed to reduce leakage and reduce the risk of accidental pull-outs from the stoma site
- Glow green feed port provides better visibility in the dark for easier nighttime feeds

### DEVICE MATERIALS

Materials in the MiniONE® Balloon Button include: Medical-grade silicone (81%) • Medical-grade thermoplastic (16%) • Stainless steel spring (2%) • Medical-grade silicone pad print ink (1%)

## KIT CONTENTS

Qty	Component	Component Usage Details
1	MiniONE <sup>®</sup> Balloon Button (Fig 1)	Provide access into stomach for feeding, medication, decompression
2	4 x 4 inch Gauze	Clean/dry stoma area during or after placement
1	Water Soluble Lubricant	Place on balloon / stoma area for easier insertion
1	12" Bolus Feeding set with Straight Adapter*	For channeling feed / medication into button or decompression
1	12" Y-Port Feeding set with Right Angle Adapter*	For channeling feed / medication into button or decompression
1	Large Syringe (Catheter or ENFit <sup>®</sup> )*	For channeling feed / medication / decompression with feeding set
1	Luer-Slip Syringe (5/6ml 12-14F), (10/12ml 16-24F)	For balloon inflation or channeling feed/medication
1	Optional Introducer / Stiffener (12-18F only)	Provides stiffness during placement of device

\*The feeding sets and large syringe are not included with the SM1-5 kit configurations of the MiniONE<sup>®</sup> Balloon Button device.

## TYPE OF USAGE – INITIAL PLACEMENT VS REPLACEMENT

The MiniONE<sup>®</sup> Balloon Button may be placed either during an initial placement procedure or as a replacement device.

### **PERFORMING AN INITIAL PLACEMENT**

If the patient does not currently have a stoma tract for placement of the MiniONE<sup>®</sup> Balloon Button device, a new stoma tract will need to be created. This process can only be completed by a healthcare professional. Follow the instructions below starting with the **SUGGESTED INITIAL PLACEMENT PROCEDURES** section in order to properly locate a proper stoma location, secure the stomach to the abdominal wall, create a stoma tract, dilate and measure the stoma, and place the device.

### **REPLACING A FEEDING DEVICE IN AN ESTABLISHED STOMA SITE**

If it is determined that a device placed in an existing stoma needs to be replaced (due to non-optimal functionality or as a pro-active replacement), the current device can be removed from the stoma and the new device can be inserted in the same stoma site. If performing a device replacement, skip straight to the **DEVICE PLACEMENT PROCEDURE** section for the proper method of placing the new device.

**NOTE:** Device replacement can be performed by the healthcare professional or at home by the patient/caregiver. Do not attempt to replace the device until first discussing the procedure with your healthcare professional.

**CAUTION:** Selection of the correct size feeding device is critical for safety and comfort of the patient. An inappropriately sized selection of the correct sized device can cause necrosis, buried bumper syndrome, and/or hypergranulation tissue. If patient size/weight has changed since the device was placed, fit is too tight or too loose, or it has been over six months since the device was last measured, it is recommended that you speak to your healthcare professional to ensure that device size does not need changed.

## SUGGESTED INITIAL PLACEMENT PROCEDURES

**WARNING: THE MINIONE<sup>®</sup> BALLOON BUTTON SHOULD ONLY BE USED BY OR UNDER THE SUPERVISION OF PERSONNEL TRAINED IN PERCUTANEOUS GASTROSTOMY TUBE PLACEMENT. A THOROUGH UNDERSTANDING OF THE TECHNICAL PRINCIPLES, CLINICAL APPLICATIONS, AND RISKS ASSOCIATED WITH PERCUTANEOUS GASTROSTOMY TUBE PLACEMENT IS RECOMMENDED PRIOR TO USING THIS DEVICE.**

Initial placement of the MiniONE<sup>®</sup> Balloon Button may be performed percutaneously under laparoscopic, fluoroscopic, or endoscopic guidance. Follow the instructions below for the operator's preferred method of placement.

**NOTE:** The insertion site for infants and children should be high on the greater curvature to prevent occlusion of the pylorus when the balloon is inflated.

### **OPTION A: RADIOLOGIC PLACEMENT**

1. Place the patient in the supine position.
2. Prep and sedate the patient according to clinical protocol.
3. Ensure that the left lobe of the liver is not over the fundus or the body of the stomach.
4. Identify the medial edge of the liver by CT scan or ultrasound.
5. Glucagon 0.5 to 1.0 mg IV may be administered to diminish gastric peristalsis.

**CAUTION:** Consult glucagon instructions for rate of iv injection and recommendations for use with insulin dependent patients.

6. Insufflate the stomach with air using a nasogastric catheter, usually 500 to 1,000 ml or until adequate distention is achieved. It is often necessary to continue air insufflation during the procedure, especially at the time of needle puncture and tract dilation, to keep the stomach distended so as to appose the gastric wall against the anterior abdominal wall.
7. Choose a catheter insertion site in the left sub-costal region, preferably over the lateral aspect or lateral to the rectus abdominal muscle (N.B. the superior epigastric artery courses along the medial aspect of the rectus) and directly over the body of the stomach toward the greater curvature. Using fluoroscopy, choose a location that allows as direct a vertical needle path as possible. Obtain a cross table lateral view prior to placement of gastrostomy when interposed colon or small bowel anterior to the stomach is suspected.

**NOTE:** PO/NG contrast may be administered the night prior or an enema administered prior to placement to pacify the transverse colon.

8. Prep and drape according to facility protocol.
9. Follow instructions in the **GASTROPEXY PLACEMENT** section.

### **OPTION B: ENDOSCOPIC PLACEMENT**

1. Perform routine Esophagogastroduodenoscopy (EGD). Once the procedure is complete and no abnormalities are identified that could pose a contraindication to placement of the tube, place the patient in the supine position and insufflate the stomach with air.
2. Transilluminate through the anterior abdominal wall to select a gastrostomy site that is free of major vessels, viscera and scar tissue. The site is usually one third the distance from the umbilicus to the left costal margin at the midclavicular line.
3. Depress the intended insertion site with a finger. The endoscopist should clearly see the resulting depression on the anterior surface of the gastric wall.
4. Prep and drape the skin at the selected insertion site.
5. Follow instructions in the **GASTROPEXY PLACEMENT** section.

## GASTROPEXY PLACEMENT

**WARNING: INITIAL PLACEMENT OF THE MINIONE® BALLOON BUTTON REQUIRES THAT A GASTROPEXY IS PERFORMED TO AFFIX THE STOMACH WALL TO THE ANTERIOR ABDOMINAL WALL. DO NOT USE THE RETAINING BALLOON OF THE DEVICE AS A GASTROPEXY DEVICE. AN EARLY BALLOON FAILURE MAY PREVENT THE STOMACH FROM ATTACHING TO THE ANTERIOR ABDOMINAL WALL.**

**CAUTION: It is recommended to perform a three point gastropexy in a triangle configuration to ensure attachment of the gastric wall to the anterior abdominal wall.**

1. Place a skin mark at the tube insertion site. Define the gastropexy pattern by placing three skin marks equidistant from the tube insertion site and in a triangle configuration.

**WARNING: ALLOW ADEQUATE DISTANCE BETWEEN THE INSERTION SITE AND GASTROPEXY PLACEMENT TO PREVENT INTERFERENCE OF THE T-FASTENER AND INFLATED BALLOON.**

2. Localize the puncture sites with 1% lidocaine and administer local anesthesia to the skin and peritoneum.
3. Place the first T-Fastener and confirm intragastric position. Repeat the procedure until all three T-Fasteners are inserted at the corners of the triangle.
4. Secure the stomach to the anterior abdominal wall and complete the procedure.
5. Follow instructions in the **CREATE THE STOMA TRACT** section.

**CAUTION: After the MiniONE® Balloon Button is placed, do not remove the gastropexy sutures until the stoma site fully heals and the stomach is fully affixed to the anterior abdominal wall.**

## CREATE THE STOMA TRACT

1. Create the stoma tract with the stomach still insufflated and in apposition to the abdominal wall. Identify the puncture site at the center of the gastropexy pattern. With fluoroscopic guidance confirm that the site overlies the distal body of the stomach below the costal margin and above the transverse colon.

**CAUTION: Avoid the epigastric artery that coursed at the junction of the medial two-thirds and lateral one-third of the rectus muscle.**

**WARNING: TAKE CARE NOT TO ADVANCE THE PUNCTURE NEEDLE TOO DEEPLY IN ORDER TO AVOID PUNCTURING THE POSTERIOR GASTRIC WALL, PANCREAS, LEFT KIDNEY, AORTA OR SPLEEN.**

2. Anesthetize the puncture site with local injection of 1% lidocaine down to the peritoneal surface.
3. Insert a .038" compatible introducer needle at the center of the gastropexy pattern into the gastric lumen directed toward the pylorus.

**NOTE:** The best angle of insertion is a 45 degree angle to the surface of the skin.

### **OPTION A. RADIOLOGIC PLACEMENT ONLY (4A & 5A):**

- 4A. Use fluoroscopic visualization to verify correct needle placement. Additionally, to aid in verification, a water filled syringe may be attached to the needle hub and air aspirated from the gastric lumen.

**NOTE:** Contrast may be injected upon return of air to visualize gastric folds and confirm position.

- 5A. Advance a guidewire, up to .038", through the needle and into the stomach. Confirm position.

### **OPTION B. ENDOSCOPIC PLACEMENT ONLY (4B & 5B):**

- 4B. Use endoscopic visualization to verify correct needle placement.
- 5B. Advance a guidewire, up to .038", through the needle into the stomach. Using endoscopic visualization, grasp the guidewire with atraumatic forceps.
6. Remove the introducer needle, leaving the guidewire in place and dispose of according to facility protocol.
7. Follow instructions in the **DILATION** section.

## DILATION

1. Use a #11 scalpel blade to create a small skin incision that extends alongside the guidewire, downward through the subcutaneous tissue and fascia of the abdominal musculature. After the incision is made, dispose of according to facility protocol.
2. Advance a dilator over the guidewire and dilate the stoma tract to the desired size.
3. Remove the dilator over the guidewire, leaving the guidewire in place.
4. Measure the Stoma Length with the AMT Stoma Measuring Device by following the **MEASURE THE STOMA LENGTH** section.

**CAUTION: Avoid excessive dilation of stoma tract as this can increase the possibility of the balloon pulling through the stoma.**

## MEASURE THE STOMA LENGTH

**CAUTION: Selection of the correct size feeding device is critical for the safety and comfort of the patient. Measure the length of the patient's stoma with the AMT balloon stoma measuring device (fig 2). The shaft length of the selected feeding device should be the same as the length of the stoma. An inappropriately sized feeding device can cause necrosis, buried bumper syndrome, and/or hypergranulation tissue.**

1. Please refer to Manufacturer's Directions for Use for Stoma Measuring Device being used. It is recommended that AMT's Balloon Stoma Measurement Device be used for measurement of the stoma.
2. Be sure to select the appropriate size MiniONE® Balloon Button for the abdominal wall thickness measured. If the measurement appears to be between two sizes, always select the next larger size MiniONE® Balloon Button. Once placed, the external flange should rotate easily.

**WARNING: UNDER-SIZING THE DEVICE MAY CAUSE EMBEDDING WITH EROSION INTO THE GASTRIC WALL, TISSUE NECROSIS, INFECTION, SEPSIS, ASSOCIATED SEQUELAE, AND/OR DEVICE FAILURE.**

## DEVICE PLACEMENT PROCEDURE

**CAUTION: Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.**

1. Select the correct MiniONE® Balloon Button size for placement.

**NOTE:** When replacing a device, stoma length should be periodically measured to ensure the correct MiniONE® Balloon Button size is being used. If patient size/weight has changed since device was placed, fit is too tight or loose, or it has been over six months since the device was last measured, it is recommended that you speak to your healthcare professional to ensure device length does not need changed.

2. Prior to placing the MiniONE® Balloon Button, inflate balloon through the balloon inflation port (see Fig 4) using a Luer slip syringe with distilled or sterile water to the recommended fill volume. The recommended fill-volume can be found on **Table 1** or printed above the balloon inflation port of the device. Remove the syringe and verify balloon integrity by gently squeezing the balloon to check for leaks. Visually inspect the balloon to verify symmetry. Verify that sizing information is appropriate for measured length. Reinsert the syringe and deflate all of the water from the balloon after inspection.

3. Lubricate the tip of the tube with water soluble lubricant. Do not use mineral oil or petroleum jelly. Insert optional introducer into the feed port if increased stiffness is desired during placement.


4. Gently guide the tube through the stoma and into the stomach until the external flange is flush to the skin.

5. Remove Introducer (if used in step 3).

6. Inflate balloon with distilled or sterile water according to the fill volume in the chart in **Table 1**.

7. Gently lift the tabs and check for signs of gastric leakage.

**NOTE:** If leakage is observed, increase balloon volume in increments of 0.5-1 ml. **Do not exceed maximum fill volume.**

Fr Size	Minimum Fill Volume	Recommended Fill Volume	Maximum Fill Volume
12 F	2 ml	2.5 ml	3 ml
14 F	3 ml	4 ml	5 ml
16 F	4 ml	6 ml	8 ml
18 F	6 ml	8 ml	10 ml
20 F	7 ml	10 ml	15 ml
24 F	7 ml	10 ml	15 ml

## PLACEMENT VERIFICATION

1. Ensure the feeding set clamp is closed and attach the feeding set to the button by lining up the dark line on the feeding set connector with the dark line on the button. Fully press the feeding set connector into the button. Turn 3/4 to the right (clock-wise) to lock the feeding set in place.

2. Attach a large diameter catheter / screw tip syringe with 5ml of water into the feeding set. Open the feeding set clamp. Aspirate for contents. Spontaneous return of gastric contents should occur.

3. After air and/or gastric contents are observed, flush the device with water.

**WARNING: NEVER INJECT AIR INTO THE MINIONE® BALLOON BUTTON.**

**WARNING: NEVER CONNECT THE EXTENSION SET TO THE BALLOON INFLATION PORT.**

4. When placement is confirmed, gastric decompression or feeding channeling can begin.

## FEEDING SET INSTRUCTIONS FOR USE

**WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO COMPATIBLE ENTERAL DEVICES. DO NOT USE FOR NON-ENTERAL APPLICATIONS.**

**WARNING: IF USING A LEGACY-STYLE (NON-ENFIT®) EXTENSION SET, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO THE FOLLOWING SYSTEMS: BREATHING, CONES AND SOCKETS OF ANESTHETIC AND RESPIRATORY EQUIPMENT, INTRAVENOUS, LIMB CUFF, NEURAXIAL CONNECTORS, NIPPLES OF RESPIRATORY THERAPY EQUIPMENT, URINARY, AND TEMPERATURE SENSOR CONNECTORS OF RESPIRATORY HUMIDIFYING EQUIPMENT.**

1. Inspect contents for damage. If damaged, do not use. Obtain another package.

2. The feeding set may be used for syringe, gravity or pump feeding or decompression.

3. Ensure the feeding set clamp is closed and attach the feeding set to the button by lining up the dark line on the feeding set connector with the dark line on the button. Fully press the feeding set connector into the button. Turn 3/4 to the right (clock-wise) to lock the feeding set in place.

4. Attach the opposite end of the feeding set to the connector being used. If using a bolus or luer connector, firmly insert the connector while using a slight rotation into the feeding set to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the feeding set while avoiding excessive force to lock in place. Once connected, open the clamp to allow flow.

**WARNING: ONLY TIGHTEN BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.**

**WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.**

**WARNING: IF FEEDING SET IS NOT PROPERLY PLACED AND LOCKED, LEAKAGE MAY OCCUR. WHEN TWISTING THE CONNECTOR HOLD THE BUTTON IN PLACE TO AVOID ROTATION WHILE PLACING THE FEEDING SET. NEVER INSERT FEEDING SET CONNECTOR INTO A BALLOON INFLATION PORT; THIS MAY CAUSE THE INTERNAL BALLOON TO DEFLATE.**

5. When feeding is complete, flush with 5-10ml water. Close the clamp to prevent leakage while removing the feeding set. Remove the feeding set by holding the button in place, grasp the feeding connector, twist 3/4 turn to the left (counter clock-wise) to unlock, and gently remove the feeding set from the button. Snap button plug in place to keep lumen clean.

6. The feeding set should be cleansed with soap and water and thoroughly rinsed. Store in a clean, dry place. Never use a dishwasher for cleaning.

## PLACEMENT CARE INSTRUCTIONS

1. The MiniONE<sup>®</sup> Balloon Button should be replaced periodically for optimal performance. If in place for three months or longer, frequent performance checks are recommended. Clogging and/or reduced flow are indications of diminished performance. Refer to the **TROUBLESHOOTING** section for additional signs of reduced performance or failure.
2. Be sure the balloon is within the stomach and the MiniONE<sup>®</sup> Balloon Button rotates freely **before the feedings begin**. Inject distilled or sterile water into the feeding port to determine proper function of the valve.
3. The stoma area should be cleansed daily with a mild soap and water. The stoma site should be clean and dry at all times. The MiniONE<sup>®</sup> Balloon Button should be rotated daily for site hygiene.

## MEDICATION CHANNELING

Liquid medication is preferred when available. If solid medication is required, consult with your physician regarding whether it is safe to crush the medication. If safe, medication should be crushed as fine as possible (into powder form) and dissolved into water before channeling medication through device. Channeling solid medication through the tubing that is not properly crushed can lead to blockage in the tubing. Never crush enteric coated medication or mix medication with formula. Flush device with water after channeling medication.

## DECOMPRESSION

Decompression can be performed if there is a need to check for residuals, or if the patient tends to retain air during feeding causing bloating and discomfort. To decompress, connect the feeding set to the interlock connector and drain stomach contents into a container. After decompression is complete, flush the feeding set with water and remove the feeding set from the device.

## PROPER FLUSHING GUIDELINES

Proper flushing techniques and frequency can help prevent tubing blockage, clogs, and tubing failure. Follow these flushing guidelines to maintain optimal device flow conditions:

- Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain patency.
- Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube before and after medication channeling and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

**UNCLOGGING A DEVICE:** First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a syringe to an extension set and attach into the interlock connector. Fill the syringe with warm water and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced.

**CAUTION:** Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.

## MRI SAFETY INFORMATION

Non-clinical testing demonstrated that the MiniONE<sup>®</sup> Balloon Button is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 1,000-gauss/cm (10-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the MiniONE<sup>®</sup> Balloon Button is expected to produce a maximum temperature rise of 1.6°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the MiniONE<sup>®</sup> Balloon Button extends approximately 20-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

## DEVICE LONGEVITY

**MiniONE<sup>®</sup> Balloon Button:** Low-profile balloon feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness.

Exact device longevity cannot be predicted. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, incorrect stoma length measurement, and overall tube care.

For optimal performance, it is recommended that MiniONE<sup>®</sup> Balloon Button device be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure. If devices are failing or performance is degrading earlier than the typical range for device longevity, it is recommended that you speak with your healthcare professional regarding eliminating common factors that can lead to early device degradation. Also reference the **TROUBLESHOOTING** section for additional information regarding device performance issues.

**NOTE:** To help prevent unneeded hospital visits, it is recommended that a spare device is kept on hand at all times for replacement in case device failure occurs prior to scheduled replacement.

**Feeding Sets:** Feeding sets are meant to be periodically replaced for optimal performance and cleanliness. AMT recommends that the feeding set be changed at least every 2 weeks or as often as indicated by your healthcare professional. Device performance and functionality can degrade over time depending on usage and environmental conditions. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, frequency of use, trauma to the device, and overall feeding set care. Device should be replaced if leakage, cracks, tears, excessive residue build-up, mold, fungus, or other signs of failure are noted. Some usage types can wear down device components quicker than others. If tubing is becoming hardened, replacement is recommended in near future to avoid failure.

## REMOVAL OF THE MINIONE® BALLOON BUTTON

1. Deflate the balloon by attaching the Luer slip tip syringe into the balloon inflation port (Fig. 4) to withdraw the water.
2. Once the balloon is fully deflated, gently remove the device from the stoma site.
3. Replace as necessary following the MiniONE® Balloon Button directions for use.

**NOTE:** Spontaneous closure of the stoma may occur as early as 24 hours after removal. Insert new device if enteral feeding by this route is still intended.

**WARNING: DO NOT CUT OFF THE INTERNAL BOLSTER OR TUBING, OR ALLOW TO PASS THROUGH THE INTESTINAL TRACT.**

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

## TROUBLESHOOTING

Long-term device performance and functionality are dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your feeding device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

<p><b>Anti-reflux valve leakage or blockage:</b> Leakage/blockage of the anti-reflux valve typically occurs due to residue (feed, medication, gastric contents, etc) becoming stuck in the valve area, preventing the valve from fully closing. Make sure device is flushed after each use. Leakage can also be caused due to excessive pressure in stomach (see <b>DECOMPRESSION</b> section). The valve can also become inverted in rare cases. Insert extension set into port to reset valve if this occurs. Avoid excessive pressure during channeling through the device, as this can cause a tear to occur in the valve over time.</p>
<p><b>Leaking of balloon volume:</b> If balloon is deflating, fully deflate the balloon and remove from stoma. Once removed, inflate balloon with recommended fill-volume. Check balloon for leakage by gently massaging tubing and balloon. If no leakage is noted, deflate the balloon, re-insert balloon back into stoma, and re-inflate the balloon to desired fill-volume. Do not exceed maximum fill volume. Only access the balloon inflation port for balloon inflation and deflation. Other uses can result in balloon leakage or device failure. <b>NOTE:</b> The balloon is made of a semi-permeable material and can lose a small amount of fill-volume over time depending on environment and usage conditions.</p>
<p><b>A tear has formed:</b> Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.</p>
<p><b>Tubing has reduced flow or has become clogged:</b> Tubing can become blocked due to not properly flushing after each use, use of thick or improperly crushed medications, use of thick feeds/formulas, gastric reflux, and/or fungus growth. If clogged, refer to the <b>PROPER FLUSHING GUIDELINES</b> section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced.</p>
<p><b>Foul smell coming from the device:</b> Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed and stoma site should be gently cleaned with soap and warm water. If foul smell does not go away, it is recommended that you contact your healthcare professional.</p>
<p><b>Balloon inflation valve leakage:</b> Leakage from this valve typically occurs due to residue in the fill-valve. Always use a clean syringe when inflating the balloon and only inflate water through the valve. Leakage can also occur due to valve becoming stuck if syringe is inserted too hard into valve. Insert the syringe into the balloon inflation port if leakage is occurring to attempt to reset valve. Several attempts may be necessary before valve resets.</p>
<p><b>Device has become discolored:</b> The device can become discolored over days to months of use. This is normal depending on the types of feeds and medications being used with the device.</p>
<p><b>Balloon is misshaped:</b> Be sure to inflate and inspect balloon prior to placement to check balloon symmetry. Balloons becoming misshaped typically occur due to excessive force or tension on the device (pulling device out of stoma while balloon is inflated). This can occur if the device is too short for the patient's stoma site. Additionally, a device that is placed too close to the pylorus can result in a misshaped balloon and blockage of the pylorus. Balloons can be gently massaged with fingers back into symmetry if not excessively misshaped. Device should be replaced if balloon is found excessively misshaped.</p>
<p><b>Balloon will not inflate or deflate:</b> Inflation and deflation problems typically occur due to residue blocking the fill-lumen. Always use a clean syringe when inflating the balloon. In some rare occurrences, fungus can grow inside of the tubing and block the fill-lumen. Fungus growth can occur depending on patient environment and feeds/medications being channeled through the device. If balloon will not deflate, contact your healthcare professional for removal. If deflation problem is the result of fungal growth, eliminating source of fungal growth or anti-fungal medication may be required.</p>
<p><b>Interlock failure or cracked:</b> The interlock has been designed to withstand extreme forces without detaching or cracking. However, the strength of the bond and material can reduce over prolonged use depending on medications and feeds used through the device. Excessive forces over an extended period of use can also reduce bond / material strength. The device should be replaced if interlock is found cracked, leaking, or separating from the device.</p>
<p><b>Device fit is too tight or too loose:</b> Balloon fit can be adjusted by changing the balloon inflation volume within the balloon inflation range in <b>Table 1</b>. If balloon fit is too loose, increase balloon fill-volume above recommended but not above max fill-volume. If balloon fit is too tight, decrease fill-volume below recommended but not below min fill-volume. If device does not properly fit with the fill-volume range, a new device length will be needed.</p>
<p><b>Balloon failure:</b> Early balloon failure can occur due to a number of patient or environmental factors, including but not limited to: gastric pH, diet, certain medications, improper balloon fill volume, placement of device, improperly channeling food/medication into balloon port, trauma, contact with a sharp or abrasive material, incorrect stoma length measurement, and overall care of the device.</p>
<p><b>Plug will not stay closed:</b> Ensure that the plug is being firmly and fully pressed into the interlock connector. If the plug is not staying closed, check the plug and feed-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.</p>

## THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

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<b>EN</b>	<b>STERILE EO</b>	Sterilized Using Ethylene Oxide	Rx Only	Prescription Only
 Do Not Resterilize	 Single Use Only	 Do not use if package is damaged and consult instructions for use	 MR Conditional	
 Not made with natural rubber latex	<b>MD</b> Medical Device	 Not made with DEHP (Di(2-ethylhexyl) Phthalate)		
ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.	MIC-KEY® feedset compatible			
	MIC-KEY® is a registered trademark of Avent, Inc.			
	MiniONE® is a registered trademark of Applied Medical Technology, Inc.			

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MiniONE® on Applied Medical Technology, Inc. registreeritud kaubamärk.
MIC-KEY® az Avent, Inc. bejegyzett márkaneve.
A MiniONE® az Applied Medical Technology, Inc. bejegyzett védjegye.
MIC-KEY® ir reģistrēta Avent, Inc. preču zīme.
MiniONE® ir reģistrēta Applied Medical Technology, Inc. preču zīme.
„MIC-KEY®“ yra registruotasis „Avent, Inc.“ prekių ženklas.