

TLC® Self-Retaining Retractor System

- ES** Sistema Retractor Autorretenedor TLC®
- DE** TLC® - Selbsthaltendes Retraktionssystem
- IT** Auto-Retrazione Divaricatore TLC® Sistema
- FR** Autostatique Système D'écarteur TLC®
- NO** TLC® Selvhøldende Sårhakesystemet
- SV** TLC® - Självhållande Retraktorsystemet
- PT** Autoestáticos Sistema de Afastadores TLC®
- NL** TLC® Zelfdragend Retractorsysteem
- DA** TLC® Selv-Retinerende Retractor Systemet
- FI** TLC® -Itsekiinnittyvä Järjestelmä
- EL** TLC® Αυτοσυγκρατούμενο Σύστημα Συγκράτησης
- AR** منظومة مبيعاد ذاتية التثبيت
- BG** Самозадържаща Система за екартиране TLC®
- CS** Samorozvírací Retraktorový Systém TLC®
- ET** TLC® Isefikseeruv Retraktorisüsteem
- HU** TLC® seberterpesztő rendszer
- JA** TLC® 自留 Retractor System
- LV** TLC® Pašstipriņoša Spriegotāja Sistēma
- LT** TLC® Savaimė Besilaikanti Retraktorius Sistema
- PL** Samoretrakcyjny System Retrakcyjny TLC®
- PT-BR** Autorretido Sistema Retrator TLC®
- RO** Autoportant Sistem Retractor TLC®
- SK** TLC® Samonosný Retrakčný Systém
- SL** Samodržalni Retrakcijski Sistem TLC®
- TR** TLC® Kendinden Tespitlii Retraktör Sistemi

Rx ONLY
STERILE | EO

MD

TLC® Self-Retaining Retractor System

INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

The device is provided STERILE for single use. Prior to use, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

DESCRIPTION:

The TLC® Retractor System is a sterile, single use, self-retaining retraction system.

INTENDED USE:

The TLC® Retractor System is intended to aid in tissue retraction during a surgical procedure. The TLC® Retractor System is intended for use during surgical operations for patients ranging in age from infant to elderly.

INDICATIONS FOR USE:

The TLC® Retractor System is intended to aid in tissue retraction during a surgical procedure.

CONTRAINDICATIONS:

At this time, there are no contraindications for use of the TLC® Retractor System.

COMPLICATIONS:

At this time, there are no complications associated with the use of the TLC® Retractor System.

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

CLINICAL BENEFITS:

Clinical Benefits to be expected when using the TLC® Retractor System include but are not limited to: Allows surgeons to retract and secure delicate soft tissue for a wide variety of procedures • Self-retaining design creates an organized surgical site which allows the surgeons to focus on the procedure while providing O.R. staff the freedom to perform other tasks • Molded disposable plastic ring is cost-effective, reducing labor, cleaning and re-sterilization costs • Ability to perform procedures without stay-sutures or manual retraction • Increases the surgeon's visualization, organization, and control of surgical sites • Increases control with uniform tension that allows surgeons better tissue control as the retractor is precisely positioned around the surgical site

PERFORMANCE CHARACTERISTICS:

Performance Characteristics of the TLC® Retractor System include but are not limited to: Designed with hinges that adjust to individual patient anatomy which creates faster deployment • Designed without thumbscrews that would need to be adjusted or would obstruct the procedure • Low profile design allows for greater visualization and unimpeded access of the surgical site • Ribbed stays provide a secure fit, virtually eliminating slippage and risk of popping out • Available in a variety of shapes and sizes to cover a broad spectrum of surgical procedures • One-piece construction eliminates the risk of detachment • Concave structure makes stays easy to access and easy to place

WARNINGS:

1. Use of excessive force may cause tissue damage or create other hazards.
2. Applying unbalanced forces to the retractor system may cause the retractor system to fail and cause impairment to surgery or other hazards.
3. Metal components of the retractor system including retractor hook stays may be sharp and may damage or perforate skin, gloves, implants, vessels, tissues, organs, or other surgical tools. Use caution when palpating tissue or using surgical implants while using this product.
4. Improper use may cause stays to suddenly break under tension and may lead to patient or operator injury. Do not use stays in retractor frames which have sharp, serrated, or other rigid uneven surfaces or stay notches.
5. This device is intended for single use. Do not reuse, re-sterilize or reprocess this medical device. Doing so may compromise biocompatibility characteristics, device performance and/or material integrity; any of which may result in potential patient injury, illness and/or death.

INSTRUCTIONS FOR USE


NOTE: The TLC® Retractor System is for use by qualified healthcare professionals during a surgical procedure.

1. When using the Retractor Frame and Stays, bend the hinges of the Retractor Frame as needed and position the Retractor Frame around the surgical site. Holding the Retractor Frame with one hand, fix one Hook or Rake Stay into the appropriate tissue, apply appropriate retraction, and then place the tail into the Retractor Frame notch. Repeat on the opposite side using balanced tension. Place additional Hook or Rake Stays as needed while maintaining balanced forces.
2. When using Deavers, place the end of the large curl into the surgical site to the desired tissue depth. Gently pull on the small curl to retract the tissue to the desired exposure.
3. When using Straps, place one end of the Strap in the Retractor Frame near the hinges. Applying firm tension, place the other end of the Strap into the opposite side of the Retractor Frame.

NOTE: The device may be a potential biohazard. The device can be disposed of by following facility protocol and in accordance with local, state, and federal laws and regulations.

THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

STERILE EO Sterilized Using Ethylene Oxide		Rx ONLY Prescription Only	
Do Not Resterilize	Single Use Only	Do not use if package is damaged and consult instructions for use	
Not made with natural rubber latex	Medical Device	Not made with DEHP (Di(2-ethylhexyl) Phthalate)	

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Applied Medical Technology, Inc.