CE  
2862

# NutraGlide®

## Nasal Feeding Tube

ES	NutraGlide® Sonda de alimentación nasal
DE	NutraGlide® Nasale Ernährungssonde
IT	NutraGlide® Sondino di alimentazione nasale
FR	NutraGlide® Sonde d'alimentation nasale
NO	NutraGlide® Nasal matesonde
SV	NutraGlide® Nässond
PT-BR	NutraGlide® Tubo de alimentação nasal
NL	NutraGlide® Nasale voedingssonde
DA	NutraGlide® Nasal ernæringslange
FI	NutraGlide® Nenäruokintaletku
AR	أنبوب التغذية الأنفي NutraGlide®
BG	NutraGlide® Назогастрална сонда
CS	NutraGlide® Nazální výživová sonda
ET	NutraGlide® Nasogastraalsond
EL	NutraGlide® Ρινικός σωλήνας σίτισης
HU	NutraGlide® Orron keresztüli táplálócső
JA	NutraGlide® 鼻腔栄養チューブ
LV	NutraGlide® Nazāla barošanas zonde
LT	NutraGlide® Nosies kaniulė
PT	NutraGlide® Tubo de alimentação nasal
PL	NutraGlide® Zgłębnik donosowy do żywienia
RO	NutraGlide® Tub de alimentare pe cale nazală
SK	NutraGlide® Hadička na podávanie výživy nosom
SL	NutraGlide® Nazalna cevka za hranjenje
TR	NutraGlide® Nazal Besleme Tüpü


Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Contact: 1-800-869-7382 [cs@appliedmedical.net](mailto:cs@appliedmedical.net) [www.appliedmedical.net](http://www.appliedmedical.net)  
Applied Medical Technology Inc. 8006 Katherine Boulevard Brecksville OH, 44141 USA


Rx Only

**EN:** Legacy Connector  
**ES:** Conector heredado  
**DE:** Legacy-Anschluss  
**IT:** Connettore Legacy  
**FR:** Connecteur pré-existant  
**NO:** Eldre type kontakt  
**SV:** Äldre koppling  
**PT-BR:** Conector legado  
**NL:** Legacy aansluiting  
**DA:** Ældre kontakt  
**FI:** Perinteinen liitin  
**AR:** مقبس التوصيل التقليدي  
**BG:** Традиционен конектор

**CS:** Starší konektor  
**ET:** Universaalne ühendus  
**EL:** Σύνδεσμος Legacy  
**HU:** Régi típusú csatlakozó  
**JA:** レガシーコネクタ  
**LV:** Vēsturiskais savienotājs  
**LT:** „Legacy“ jungtis  
**PT:** Conector Legacy  
**PL:** Starsze złącze  
**RO:** Conector vechi  
**SK:** Pôvodný konektor  
**SL:** Starejši konektor  
**TR:** Eski Tip Konektör

**EN:** ENFit® Connector  
**ES:** Conector ENFit®  
**DE:** ENFit®-Anschluss  
**IT:** Connettore ENFit®  
**FR:** Connecteur ENFit®  
**NO:** ENFit®-kontakt  
**SV:** ENFit®-koppling  
**PT-BR:** Conector ENFit®  
**NL:** ENFit® aansluiting  
**DA:** ENFit®-kontakt  
**FI:** ENFit®-liitin  
**AR:** مقبس توصيل ENFit®  
**BG:** Конектор ENFit®

**CS:** Konektor ENFit®  
**ET:** ENFit® ühendus  
**EL:** Σύνδεσμος ENFit®  
**HU:** ENFit® csatlakozó  
**JA:** ENFit®コネクタ  
**LV:** ENFit® savienotājs  
**LT:** „ENFit®“ jungtis  
**PT:** Conector ENFit®  
**PL:** Złącze ENFit®  
**RO:** Conector ENFit®  
**SK:** Konektor ENFit®  
**SL:** Konektor ENFit®  
**TR:** ENFit® Konektör

## INSTRUCTIONS FOR USE

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. The NutraGlide® is intended for single use. This device is intended to connect to compatible enteral devices only. Do not use for non-enteral applications.

### DESCRIPTION

The NutraGlide® is a nasal feeding tube with a Y-Port, widening lumen diameter, lubricious coating and optional stylet.

### INTENDED USE

The NutraGlide® is intended to deliver enteral nutrition, medication, and fluid to neonatal, pediatric and adult patients via the nasogastric or naso-intestinal route. The NutraGlide® Nasal Feeding Tube is intended to be placed by trained healthcare professionals and used by clinicians and trained caregivers/users.

**WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE OR REPROCESS THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.**

### INDICATIONS FOR USE

The NutraGlide® Nasal Feeding Tube is indicated for the administration of nutrition, fluids, and medications in neonatal, pediatric, and adult patients.

### CONTRAINDICATIONS

The NutraGlide® is contraindicated for use in patients with anatomical abnormalities or diseases of the nose, throat, or esophagus that would prohibit safe placement of the device.

### COMPLICATIONS

Potential complications when using the NutraGlide® include but are not limited to: Patient discomfort • Nausea • Vomiting • Diarrhea • Pneumothorax • Gastrointestinal bleeding or ulceration • Gastrointestinal or esophageal perforation • Aspiration and/or aspiration pneumonia • Airway Obstruction • Tissue irritation or necrosis • Contamination

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

### CLINICAL BENEFITS, PERFORMANCE CHARACTERISTICS

Clinical Benefits to be expected when using the NutraGlide® include but are not limited to:

Allows patients who cannot eat to receive nutrition • Variety of tube lengths allows the channeling of enteral feed or medications to the stomach, duodenum or jejunum • Tube design reduces occurrence of clogs which minimizes the number of replacements • Tube design may minimize patient trauma and discomfort during placement • Variety of tube sizes (FR size) available to meet a wide range of patient needs

Performance Characteristics of the NutraGlide® include but are not limited to:

Guidewire compatible throughout the entire length of the device • Tube lumen that increases in diameter as it advances through the gastrointestinal tract to resist clog formation and ease clog removal • Centimeter markings to assist in tube placement and check for tube migration • Radiopaque polyurethane tube for x-ray placement verification • MRI Safe once stylet is removed • Y-Port to allow a second delivery option while the main port is connected • Lubricious coating for easy stylet removal (if used) • ENFi® compatible option

### DEVICE MATERIALS

Materials in the NutraGlide® include: Medical-grade elastomeric thermoplastic, non-PVC (78% - 82%) • Medical-grade thermoplastic (16% - 20%) • Medical-grade adhesive (1%) • Medical-grade surface treatment (1%) • Laser-engraved onto the medical-grade elastomeric thermoplastic tubing (<1%)

### DEVICE PLACEMENT PROCEDURE

**CAUTION: This feeding tube is to be placed by trained healthcare professionals only.**

**CAUTION: Prior to placement, please inspect all contents of the kit for damage. If the package or device is damaged, do not use the product.**

1. Read all directions and warnings prior to tube placement.
2. Place patient in sitting or semi-Fowler's position to aid in device placement. Do not lean patient forward. The patient's head and neck should not be extended.

**CAUTION: Do not use a stylet during placement if the device was not packaged with one.**

3. Estimate insertion depth for gastric placement. Using the printed marks on the tube itself, measure from the tip of the patient's nose to the earlobe, and from the earlobe to the xiphoid process.

**WARNING: ESTIMATION OF TUBING LENGTH IS CRITICAL. DO NOT INSERT EXCESS TUBING AS THIS MAY RESULT IN KINKING OF TUBE.**

4. Determine the preferred nostril for insertion.
5. Provide patient with water and straw to help with tube placement.
6. For guidewire assisted placement only:

**NOTE:** Stylet cannot be used with guidewire assisted placement. Carefully remove the stylet by applying gentle traction. If resistance occurs, inject up to 10 mL of water into the tube to aid in stylet removal. Refer to Table 1 for recommended water volume to be used. Wait 60 seconds after water injection before reattempting stylet removal.

- A. Place tube onto guidewire by inserting guidewire through distal opening in tube.
- B. Continue with guidewire assisted placement following the remaining instructions below.

Table 1	
Lubricious Coating Activation	
Fr Size	Recommended Volume
5F	≤ 1 mL
6F	≤ 3 mL
8F	≤ 5 mL
10F	≤ 8 mL

## DEVICE PLACEMENT PROCEDURE

7. Apply water based lubricant to the outside of the tube to ease placement and limit patient discomfort.
8. Insert the feeding tube, being careful not to twist or pull the stylet (if used) during placement to avoid damaging the tube. Direct tube posteriorly, aiming tip parallel to nasal septum and superior surface of hard palate. Advance tube to nasopharynx. Encourage patient to sip water to help guide tube into esophagus and towards stomach.

**WARNING: COUGHING OR RESPIRATORY DISTRESS MAY INDICATE PASSAGE OF THE TUBE INTO THE AIRWAY. IF THIS IS SUSPECTED, REMOVE DEVICE IMMEDIATELY AND REINSERT.**

**WARNING: CAUTION SHOULD BE TAKEN IF ANY TYPE OF ENDOTRACHEAL DEVICE IS IN PLACE, AS THESE MAY FACILITATE THE PASSAGE OF THE FEEDING TUBE INTO THE AIRWAY. PLACEMENT OF THE FEEDING TUBE INTO THE TRACHEA OR LUNGS MAY RESULT IN SERIOUS INJURY.**

**WARNING: WHEN GUIDING THE FEEDING TUBE INTO PLACE, DISCONTINUE INSERTION IF ANY RESISTANCE IS EXPERIENCED AND NOTIFY CLINICIAN IMMEDIATELY. DO NOT USE EXCESSIVE FORCE UPON INSERTION.**

9. Utilize the printed centimeter marks on the tube to aid insertion and check tube migration.
10. Proper placement of the feeding tube must be confirmed prior to delivery of substances through the tube. Follow facility protocol to ensure proper tube position (i.e. aspiration of gastric contents, X-Ray, etc.).
11. The stylet (if used) can be removed after tube position has been confirmed. Carefully remove the stylet by applying gentle traction. If resistance occurs, inject up to 10 mL of water into the tube to aid in stylet removal. Refer to Table 1 for recommended water volume to be used. Wait 60 seconds after water injection before reattempting stylet removal.

**WARNING: NEVER REINSERT STYLET WHEN TUBE IS IN PATIENT.**

12. If gastric placement is desired, secure the tube. AMT recommends using the AMT Bridle™ or AMT Bridle Pro® for tube securement. Document the tube length reading on the tube at the opening of the nare (or nostril).
13. If small bowel placement is desired, secure tube to cheek loosely (using tape or AMT CINCH®), leaving a loop of excess tube between nose and securement point. The tube tip will move into the small bowel as a result of gastric peristalsis. Placing the patient into Fowler's position on the right side will aid with tube migration. After migration of tube into the small bowel is complete, secure the tube. AMT recommends using the AMT Bridle™ or AMT Bridle Pro® for tube securement. Document the tube length reading on the tube at the opening of the nare (or nostril).
14. Feeding can begin per physician's order and following facility protocol.

**CAUTION: Tube position must be confirmed prior to use.**

**WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO COMPATIBLE ENTERAL DEVICES. DO NOT USE FOR NON-ENTERAL APPLICATIONS.**

**NOTE:** If using a legacy-style (non-ENFit®) device, the device has the potential to misconnect to the following systems: breathing, intravenous, limb cuff, neuralxial connectors, nipples of respiratory therapy equipment, urinary, and temperature sensor connectors of respiratory humidifying equipment.

**WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.**

## MEDICATION CHANNELING

Liquid medication is preferred when available. If solid medication is required, consult with your physician regarding whether it is safe to crush the medication. If safe, medication should be crushed as finely as possible (into powder form) and dissolved into water before channeling medication through device. Channeling solid medication through the tubing that is not properly crushed can lead to blockage in the tubing. Never crush enteric coated medication or mix medication with formula. Flush device with water after channeling medication.

## PROPER FLUSHING GUIDELINES

Proper flushing techniques and frequency can help prevent tubing blockage, clogs, and tubing failure. Follow these flushing guidelines to maintain optimal device flow conditions:

- Use room temperature water for tube flushing.
- The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants.

**NOTE:** Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain patency.

- Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, and before and after every intermittent feeding.
- Flush the feeding tube before and after medication channeling and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.

**UNCLOGGING A DEVICE:**

Connect a syringe filled with warm water to the Feed Port and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact the healthcare professional, as the tube may need to be replaced.

**CAUTION: Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.**

## DEVICE REMOVAL

**NOTE:** The NutraGlide® has satisfied biocompatibility testing as a device for long-term use per ISO 10993-1. The feeding tube should be monitored, regularly assessed, and replaced when clinically indicated based on functionality and patient condition. The feeding tube should be monitored and regularly assessed for any resistance of flow and cleanliness of the port as well as a regular patient assessment for any nasopharynx trauma such as the occurrence of pain or the development of bleeding.

1. Flush the feeding tube with water.
2. Seal the Y-Port using the attached plugs or clamp the feeding tube. This will reduce the risk of fluid aspiration into the airway during removal.
3. Gently withdraw the device from the patient and discard the feeding tube following facility protocol.

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

## MRI SAFETY INFORMATION

The NutraGlide® is MRI Safe once stylet (if used) is removed. A patient with an indwelling device can be scanned safely.

## STORAGE

Store device in a dry, room temperature environment until ready to use.

## DEVICE LONGEVITY

Nasal feeding tubes are meant to be periodically replaced for optimal performance, functionality, and cleanliness. The NutraGlide® has satisfied biocompatibility testing as a device for long-term use per ISO 10993-1. The feeding tube should be monitored, regularly assessed, and replaced when clinically indicated based on functionality and patient condition. The feeding tube should be monitored and regularly assessed for any resistance of flow and cleanliness of the port as well as a regular patient assessment for any nasopharynx trauma such as the occurrence of pain or the development of bleeding. Device performance and functionality can degrade over time depending on many factors, including: gastric pH, medications, trauma to the device, and overall tube care. Device should be replaced if signs of failure are noted.

## TROUBLE SHOOTING

Long-term device performance and functionality are dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your feeding device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

**Tubing has reduced flow or has become clogged:** Tubing can become blocked due to not properly flushing after each use, use of thick or improperly crushed medications, use of thick feeds/formulas, gastric reflux, and/or fungus growth. If clogged, refer to the **PROPER FLUSHING GUIDELINES** section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced.

**Device has become discolored:** The device can become discolored over days to months of use. This is normal depending on the types of feeds and medications being used with the device.

**A tear has formed:** Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.

**Foul smell coming from the device:** Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed. If foul smell does not go away, it is recommended that you contact your healthcare professional.

**Plug will not stay closed:** Ensure that the plug is being firmly and fully pressed or rotated without excessive force. If the plug is not staying closed, check the plug and feed-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.

## THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

 <p>Not made with DEHP (Di(2-ethylhexyl) Phthalate)</p>	 <p>Single Use Only</p>	 <p>MR Safe</p>	 <p>Medical Device</p>
 <p>Not made with natural rubber latex</p>	 <p>Non-Sterile</p>	<p>Rx Only Prescription Only</p>	
<p>ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.</p>			
<p>The ENFit® connection is intended to reduce misconnections with other non-enteral healthcare applications.</p>			
<p>NutraGlide® is a registered trademark of Applied Medical Technology, Inc.</p>			

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ENFit®, Global Enteral Device Supplier Association, Inc. şirketinin tescilli ticari markasıdır.
NutraGlide®, Applied Medical Technology, Inc. şirketinin tescilli ticari markasıdır.