Remora® Smoke and Fluid Evacuation System

ES	Sistema de evacuación de fluidos y humos Remora®
DE	Remora® Absaugsystem für Rauch und Flüssigkeiten
IT	Sistema di evacuazione di fumo e fluidi Remora®
FR	Système d'évacuation de la fumée et des fluides Remora®
NO	Remora® røyk- og væskeevakueringsystem
SV	Remora® evakueringsystem för rök och vätska
PT-BR	Sistema de Evacuação de Fluidos e Vapores Remora®
NL	Remora®-systeem van rook- en warmteafvoer
DA	Remora® røg- og væskeevakueringsystem
FI	Remora® Savun- ja nesteenpoistojärjestelmä
AR	نظام تفريغ الدخان والسوائل Remora®
BG	Remora® Евакуационна система за дим и течности
CS	Remora® Systém odsávání kouře a tekutin
ET	Remora® suitsu ja vedelike eemaldussüsteem
EL	Σύστημα εκκένωσης καπνού και υγρών Remora®
HU	Remora® Füst- és Folyadékkelvezető Rendszer
JA	Remora® 煙および液体排出システム
LV	Remora® dūmu un šķidrumu sūkņēšanas sistēma
LT	„Remora®“ dūmų ir skysčių pašalinimo sistema
PT	Remora® Sistema de Evacuação de Fumo e Fluidos
PL	System odprowadzania dymu i płynów Remora®
RO	Sistem de evacuaire a fumului și fluidelor Remora®
SK	Systém Remora® na vypúšťanie dymu a tekutín
SL	Remora® Sistem za odvajanje dima in tekočin
TR	Remora® Duman ve Sıvı Tahliye Sistemi
KO	Remora® 매연 및 체액 배출 시스템
MS	Sistem Pemindahan Asap dan Cecair Remora®

C4406-I 02/2024

Made in USA

Remora[®] Smoke and Fluid Evacuation System

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. The device is provided STERILE for single use. Please inspect all contents of the kit for damage. If damaged or sterile barrier is breached, do not use.

INTENDED USE:

The Remora[®] is intended to remove surgical smoke and fluid caused by use of an electrosurgical tool. The Remora[®] is intended to be used by surgeons for patients of all ages. The device aids surgeons with a wide array of electrosurgical procedures that present risk to the patient and clinicians from surgical smoke created during the procedures.

INDICATIONS FOR USE:

The Remora[®] is an accessory device used to provide an air suction path for the reduction of gases, body fluids, and smoke associated with the use of an electric cauterization system.

CONTRAINDICATIONS:

At this time, there are no contraindications for use of the Remora[®].

COMPLICATIONS:

At this time there are no complications associated with use of the Remora[®].

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

CLINICAL BENEFITS:

Clinical benefits to be expected when using the Remora[®] include but are not limited to:

- Protects patient, surgeon, and surgical staff from electrosurgical smoke
- Simultaneously removes potentially harmful smoke and liquids directly from the surgical site at point of origin
- Improves conductivity of electrosurgical instruments to potentially result in a more efficient procedure by removing fluid from the surgical site
- Minimalist design allows for improved visibility of the surgical field (especially in deep surgical pockets or cavities)
- Increases efficiency, assists in keeping the operating room free of toxic smoke, and improves safety as carcinogens and other toxins are removed from the OR

PERFORMANCE CHARACTERISTICS:

Performance Characteristics of the Remora[®] include but are not limited to:

- Designed for universal compatibility - slips over most electrosurgical pencils and has a universal suction adapter which plugs into most O.R. suction sources
- Designed with unimpeded maneuverability
- Insulates the pencil tip
- Alternate ports prevent blockage and reduces trauma to tissue
- Device is OSHA and AORN smoke evacuation compliant
- Available in a standard and a micro configuration to accommodate a variety of surgical procedures
- Available in a two sizes to accommodate different tips (needle, blade, or extended) and electrosurgical pencil handles

CONTENTS:

(1) Remora[®], (1) Tube Retainer, (1) Needle Tip Adapter (short tube), (1) Clamp (included with some configurations)

INSTRUCTIONS:

NOTE: Remora[®] is for use by surgeons to aid in electrosurgical procedures.

NOTE: The Remora[®] is a smoke and fluid evacuation system that was designed to attach to (or work with) almost any electrosurgical pencil and has a universal suction adapter which plugs into most O.R. suction sources.

1. Securely fit the Remora[®] 2mm-3mm past the end of the electrode.

2. If necessary slide the needle tip adapter over the electrode, before placing the Remora[®] to achieve a secure fit.

3. Place the tube retainer around the Remora[®] tube and electrocautery pencil opposite the electrode. Insert the Remora[®] connection adapter into an operating room suction system.

WARNING: The Remora[®] must be used in accordance with the operating parameters and limitations of the applicable electrocauterization system and suction system.

WARNING: This device is intended for single use. Do not reuse, re-sterilize or reprocess this medical device. Doing so may compromise biocompatibility characteristics, device performance and/or material integrity; any of which may result in potential patient injury, illness and/or death.

NOTE: The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

R_x ONLY Prescription Only	STERILE EO Sterilized Using Ethylene Oxide
 Not made with natural rubber latex	 Not made with DEHP (Di(2-ethylhexyl) Phthalate)
 Do not re-sterilize	 Do not use if package is damaged and consult instructions for use
 Single Use Only	MD Medical Device

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© 2024 Applied Medical Technology, Inc.
Applied Medical Technology, Inc.