# AMT micro G-JET®

Micro Transgastric-Jejunal Feeding Device





Directions for Use - Tube Use and Care Instructions

For placement instructions for hospital professionals please see [Device Placement Directions for Use](#) document.

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**EN:** The AMT micro G-JET<sup>®</sup> Transgastric-Jejunal Feeding Device is considered compatible with MRI equipment according to the American Society for Testing and Materials (ASTM International), under the name F2503-08. Non-clinical trials have shown that patients with the micro G-JET<sup>®</sup> can safely undergo magnetic resonance imaging (MRI).

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Read the following warnings and instructions for use prior to using the AMT micro G-JET® device. Not following these guidelines can result in device failure, serious injury, and/or death. The AMT micro G-JET® is considered safe for the intended use when used within the guidelines listed throughout these instructions.

**PRIOR TO USE:**

- Components are supplied sterile: for single user only. Do not re-use or re-sterilize.
- Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

**DEVICE PLACEMENT AND REPLACEMENT:**

- Due to the presence of the tubing support throughout the entire length of the tubing, do not cut the device to custom length. Do not cut the distal end of the tubing to create a customized jejunal length. Doing so will eliminate the soft, tapered tip of the device and expose the tubing support.
- A gastropexy must be performed to affix the stomach to the anterior abdominal wall, the feeding tube insertion site identified, stoma tract dilated and measured prior to initial tube insertion to ensure patient safety and comfort. The length of the tube should be sufficient to be placed beyond the ligament of treitz. Do not use the retention balloon of the feeding tube as a gastropexy device. The balloon may burst and fail to attach the stomach to the anterior abdominal wall.
- Selection of the correct size feeding device is critical for the safety and comfort of the patient. Measure the length of the patient's stoma with the stoma measuring device. The shaft length of the selected feeding device should be the same as the length of the stoma. An inappropriately sized feeding device can cause necrosis, buried bumper syndrome and/or hypergranulation tissue.
- Removal is recommended to be performed by a qualified clinician. Spontaneous closure of stoma may occur as early as twenty-four (24) hours after removal. Insert a new device if enteral feeding by this route is still intended. If closure is desired, apply a dressing over the stoma site.

**GENERAL USE:**

- Do not exceed the maximum rated balloon volume. Do not use air. Do not inject contrast into the balloon.
- Never connect the jejunal port to suction. Do not measure residuals from the jejunal port. Do not use continuous or high intermittent suctions. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.
- Do not channel medications through the jejunal feeding port. This will clog the tube. To help avoid clogging, use a feeding pump to deliver formula through the jejunal lumen. Never allow formula to stand in the tube.
- This device is intended to connect to compatible enteral devices only. Do not use for non-enteral applications.
- Do not use excessive force or pressure when channeling feed or medication, or to attempt to clear a clog in the tubing. This can lead to damage of the tubing and may result in feeding the stomach. If clog cannot be cleared or tubing is clogging on a frequent basis, this may indicate that the device needs replaced. Failure to follow this warning may result in device failure and/or serious patient injury.
- Do not use excessive force to flush the tube. Excessive force can damage the tube and can cause injury to the gastrointestinal tract. **Use a 30 to 60 ml syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.**
- This device is intended for single use. Do not reuse, re-sterilize or reprocess this medical device. Doing so may compromise biocompatibility characteristics, device performance and/or material integrity; any of which may result in potential patient injury, illness and/or death.

**WHEN USING AN EXTENSION SET OR TRANSITION ADAPTER:**

- Never connect the extension set to the balloon inflation port.
- This device has the potential to misconnect with small bore connectors of other healthcare applications. Only use this device to connect to compatible enteral devices. Do not use for non-enteral applications.
- If using a legacy-style (non-ENFit®) extension set, this device has the potential to misconnect to the following systems: breathing, cones and sockets of anesthetic and respiratory equipment, intravenous, limb cuff, neuraxial connectors, nipples of respiratory therapy equipment, urinary, and temperature sensor connectors of respiratory humidifying equipment .
- Only tighten by hand. Never use excessive force or a tool to tighten a rotating connector. Improper use can lead to cracking, leakage, or other failure.
- Ensure device is connected to an enteral port only and NOT to an IV set.

**BALLOON INFLATION VOLUME**

Inflate the balloon with sterile or distilled water within the fill-volume range below.

Table 1			
FR Size	Min Volume	Recommended Volume	Max Volume
14 F	3 ML	4 ML	5 ML

- 14F sized tubes: Inflate the balloon with 3-5 ml of sterile or distilled water.
- CAUTION: Do not exceed 5 ml total balloon volume. Do not use air. Do not inject contrast into the balloon.**

**KIT CONTENTS**

**CAUTION: Components are supplied Sterile: For single user only. Do not re-use or re-sterilize. For enteral nutrition and /or medication only.**

**KIT CONTENTS**

(1) micro Transgastric-jejunal Feeding Device	(1) Small Syringe (luer slip)
(1) Jejunal Feeding Set with Glow Green™ Right Angle Adapter	(1) Small Syringe with adapter (For direct port access) <b>(AMT P/N: 4-7000)</b>
(1) Gastric Feeding Set with White Right Angle Adapter	(1) Guidewire Introducer
(1) Large Syringe (Catheter or ENFit®)	(2) 4" x 4" (10 cm x 10 cm) Gauze

## DESCRIPTION, INDICATIONS, AND DURATION OF USE

**DESCRIPTION:** The AMT micro Transgastric-Jejunal Feeding Device (**Fig. A**) provides for simultaneous gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. It enters the stomach through a gastric stoma. The tube is held in place (within the stoma tract) by means of an inflatable balloon (**Fig. A-1**) and a low-profile external bolster (**Fig. A-2**). The external bolster contains two ports, one labeled "JEJUNAL" (Glow Green™) and one labeled "GASTRIC" (white). Incorporated within each port is a one-way valve which is opened by attaching the appropriate AMT Extension Set. AMT Extension Sets (**Figs. B & C**) are color coded and connect to the feeding and drainage ports. The JEJUNAL (Glow Green™) port is used for feeding into the small intestine. The GASTRIC (white) port is used to drain the stomach with the use of low intermittent suction or gravity drainage. A third port labeled "BAL" is used to inflate and deflate the balloon.

**INTENDED USE:** The AMT micro G-JET® is intended for gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. The AMT micro G-JET® is intended to be placed by qualified healthcare professionals. The AMT micro G-JET® is intended to be used by clinicians and caregivers/users.

**INDICATIONS FOR USE:** The AMT micro Transgastric-Jejunal Feeding Device is indicated for use in adult, adolescent, child, and infant patients over 6 kg who cannot absorb adequate nutrition through the stomach, who have intestinal motility problems, gastric outlet obstruction, severe gastroesophageal reflux, are at risk of aspiration, or in those who have had previous esophagectomy or gastrectomy. The use of this tube is also clinically indicated when simultaneous gastric decompression and jejunal feeding are needed. This includes patients in whom malnutrition already exists, or may result, secondary to concurrent conditions.

**CONTRAINDICATIONS:** Contraindications for placement of a transgastric-jejunal feeding device include, but are not limited to ascites, colonic interposition, portal hypertension, peritonitis and morbid obesity.

**COMPLICATIONS:** The following complications may be associated with any transgastric-jejunal feeding device; any of which may result in potential patient injury, illness, and/or death. Skin Breakdown • Infection • Hypergranulation Tissue • Stomach or Duodenal Ulcers • Intraabdominal Leakage • Pressure Necrosis • Intussusception • Gastric or Intestinal Perforation • Peritonitis • Device Migration

**NOTE:** Verify package integrity. Do not use if package is damaged or sterile barrier compromised.

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

**CLINICAL BENEFITS:** Clinical benefits to be expected when using the AMT micro G-JET® include but are not limited to: Provides access to both the stomach and jejunum for feeding, decompression, or channeling medication • Device may be placed either during an initial placement procedure or as a replacement device • Four point stabilization of the external bolster reduces motion that can cause granulation tissue • Provides leakage protection via one-way valves to ensure patients retain nutrition and/or medication • Mutually exclusive Gastric and Jejunal ports help minimize the possibility of accidental feeding of the stomach, which could lead to choking and aspiration • Channeling of medication and feeding to the patient is immediate and sustains life • Made from medical-grade silicone to reduce skin irritation and improve patient comfort • Designed to reduce digestive problems associated with sustained pyloric opening, such as reflux and potential aspiration of bile or chyme

**PERFORMANCE CHARACTERISTICS:** Performance Characteristics of the AMT micro G-JET® include but are not limited to: Low profile design sits flush against the skin • The only low profile gastric-jejunal feeding device with an 8F jejunal segment • 8F jejunal segment is well-suited for the narrow lumen of the small intestine in pediatric patients • AMT's unique "apple" shape balloon, adapted for smaller stomachs, designed to reduce leakage and reduce the risk of accidental pull-outs from the stoma site • Designed with a gentle, atraumatic tip • Anti-kink coil spans the entire 8F jejunal segment to help prevent kinking and maintain flow • Glow Green technology provides better visibility in the dark for easier for nighttime feeds • Tapered transition for a more secure fit • Created with dual safety straps that allow independent access to the Gastric and/or Jejunal ports

**DEVICE MATERIALS:** The AMT micro G-JET® is composed of the following materials: Medical-grade silicone (64%) • Medical-grade thermoplastic (31%) • Stainless Steel Spring (4%) • Medical-grade silicone pad print ink (1%) •

**DEVICE DURATION OF USE:** Device meets biocompatibility standard ISO 10993 and is approved for permanent contact (greater than 30 days).

### AMT MICRO G-JET®

Low-profile balloon feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness. Exact device longevity cannot be predicted. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, incorrect stoma length measurement, and overall tube care.

For optimal performance, it is recommended that the AMT micro G-JET® device be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure. If devices are failing or performance is degrading earlier than the typical range for device longevity, it is recommended that you speak with your healthcare professional regarding eliminating common factors that can lead to early device degradation. Also reference the **TROUBLESHOOTING** section for additional information regarding device performance issues.

### EXTENSION SETS

Extension sets are meant to be periodically replaced for optimal performance and cleanliness. AMT recommends that the extension set be changed at least every 2 weeks or as often as indicated by your healthcare professional. Device performance and functionality can degrade over time depending on usage and environmental conditions. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, frequency of use, trauma to the device, and overall extension set care.

Device should be replaced if leakage, cracks, tears, excessive residue build-up, mold, fungus, or other signs of failure are noted. Some usage types can wear down device components quicker than others. If tubing is becoming hardened, replacement is recommended in near future to avoid failure.

## EXTENSION SET USAGE INSTRUCTIONS

### FEEDING AND DECOMPRESSION THROUGH THE AMT MICRO TRANSGASTRIC-JEJUNAL FEEDING DEVICE

The following information explains how to properly connect the AMT Extension Sets and to feed into the small intestine and decompress the stomach.

**WARNING: NEVER CONNECT THE EXTENSION SET TO THE BALLOON INFLATION PORT.**

**WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO COMPATIBLE ENTERAL DEVICES. DO NOT USE FOR NON-ENTERAL APPLICATIONS.**

**WARNING: ONLY TIGHTEN BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.**

**WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.**

**WARNING: NEVER CONNECT THE JEJUNAL PORT TO SUCTION. DO NOT MEASURE RESIDUALS FROM THE JEJUNAL PORT.**

**WARNING: IF USING A LEGACY-STYLE (NON-ENFIT®) EXTENSION SET, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO THE FOLLOWING SYSTEMS: BREATHING, CONES AND SOCKETS OF ANESTHETIC AND RESPIRATORY EQUIPMENT, INTRAVENOUS, LIMB CUFF, NEURAXIAL CONNECTORS, NIPPLES OF RESPIRATORY THERAPY EQUIPMENT, URINARY AND TEMPERATURE SENSOR CONNECTORS OF RESPIRATORY HUMIDIFYING EQUIPMENT.**

### EXTENSION SET ASSEMBLY FOR JEJUNAL FEEDING

1. Open the jejunal feeding port cover (**Fig A-3**) located at the top of the AMT micro Transgastric-Jejunal Feeding Device.
2. Insert the AMT Glow Green™ adapter extension set (**Fig B**) into the Glow Green™ port labeled "JEJUNAL" by aligning the lock and key connector. Align the black orientation marking on the set with the corresponding black orientation line on the jejunal feeding port.
3. Lock the set into the jejunal feeding port by pushing in fully and rotating the connector **CLOCKWISE** until a slight resistance is felt (approximately 3/4 turn).

**NOTE:** DO NOT rotate the connector past the stop point.

4. Attach the opposite end of the extension set to the connector being used. If using a bolus or luer connector, firmly insert the connector while using a slight rotation into the extension set to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the extension set while avoiding excessive force to lock in place. Once connected, open the clamp to allow flow.

**JEJUNAL FEEDING**

**CAUTION: Do not channel medications through the jejunal feeding port. This will clog the tube. To help avoid clogging, use a feeding pump to deliver formula through the lumen. Never allow formula to stand in the tube.**

1. Assemble the equipment: formula, feeding set, irrigating syringe, enteral feeding pump, water for flushing the tube, and extension set.
2. Wash your hands with soap and water. Rinse and dry them thoroughly.
3. Shake the formula and wipe the container tops before opening. If you do not use all the formula, cover the open container, record the date and time, and store it in the refrigerator. Throw away formula more than 24 hours old. Do not mix new formula with old formula. There is always a chance it will spoil.
4. Pour the formula into the feeding bag.
5. Prime the 12 inch extension set by injecting a prescribed amount of water or injecting no less than the recommended priming volume indicated in the table in the "General Flushing Guidelines" section. Close the clamp.
6. Connect the AMT Extension Set (Fig. B) to the Glow Green™ port labeled "JEJUNAL". Open the clamp on the extension set.
7. Inject an additional 1-2 ml or prescribed amount of water into the AMT Extension Set by using the 35 ml irrigating syringe from your kit.
8. Attach the extension set to the feeding set connector. Connect the feeding set tubing to the enteral feeding pump. Follow the manufacturer's directions to set the pump.
9. If the feed set tubing has a clamp, open it.
10. Start the feeding pump.
11. Flush the jejunal port EVERY FOUR HOURS by injecting a prescribed amount of water or by injecting no less than the recommended flushing volume indicated in the table in the "General Flushing Guidelines" section. DO NOT USE FORCE.

**NOTE:** If you are simultaneously decompressing the stomach while feeding and you see formula in the gastric drainage, stop the feeding and notify your specialist.

12. When the feeding is complete, clamp the feeding set tubing and disconnect the feeding set tubing from the extension set.

13. Flush the extension set and feeding tube by injecting a prescribed amount of water or by injecting no less than the recommended flushing volume indicated in the table in the "General Flushing Guidelines" section. Disconnect the extension set from the feeding tube port.

**EXTENSION SET ASSEMBLY FOR GASTRIC DECOMPRESSION OR MEDICATION CHANNELING**

1. Open the gastric feeding port cover (Fig A-4) located at the top of the AMT micro Transgastric-Jejunal Feeding Device.
2. Insert the AMT white adapter extension set (Fig C) into the white port labeled "GASTRIC" by aligning the lock and key connector. Align the black orientation marking on the set with the corresponding black orientation line on the gastric port.
3. Lock the set into the gastric decompression port by pushing in fully and rotating the connector CLOCKWISE until a slight resistance is felt (approximately 3/4 turn).

**NOTE:** Do not rotate the connector past the stop point.

4. Attach the opposite end of the extension set to the connector being used. If using a bolus or luer connector, firmly insert the connector while using a slight rotation into the extension set to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the extension set while avoiding excessive force to lock in place. Once connected, open the clamp to allow flow.

**CAUTION: Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.**

**GASTRIC DECOMPRESSION**

Some specialists recommend decompressing the stomach (letting the air or fluid out) before or after feeding. FOLLOW YOUR SPECIALIST'S INSTRUCTIONS. Connect the AMT Gastric Extension Set to the white port labeled "GASTRIC". The end of the Extension Set may be left open for gravity drainage or connected to low intermittent suction. This allows stomach contents and/or gas to escape. Flush the gastric port EVERY FOUR HOURS by injecting a prescribed amount of water or by injecting no less than the recommended flushing volume indicated in the table in the "General Flushing Guidelines" section. Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.

**EXTENSION SET REMOVAL**

1. When use is complete, flush by injecting a prescribed amount of water or by injecting no less than the recommended flushing volume indicated in the table in the "General Flushing Guidelines" section. Close the clamp to prevent leakage while removing the extension set. Remove the extension set by holding the button in place, grasp the extension connector, twist 3/4 turn to the left (COUNTER CLOCKWISE) to unlock, and gently remove the extension set from the button. Snap button plug in place to keep lumen clean.
2. The extension set should be cleaned with soap and water and thoroughly rinsed. Store in a clean, dry place. Never use a dishwasher for cleaning.

**MICRO G-JET® FEEDING DEVICE – GENERAL USE AND CARE**

**GENERAL FLUSHING GUIDELINES:** Proper tube flushing is the best way to avoid clogging and maintain tube patency. The following are guidelines to avoid clogging and maintain tube patency.

- The table below shows the minimum priming/displacement volume and flushing volume that should be used to maintain tube patency based on the device size and attached extension set. Please consult your physician for recommended priming and flushing volumes. Priming should occur prior to the start of feeding or medication channeling. Flushing should occur every 4 hours during continuous feeding, anytime the feeding is interrupted, after every intermittent feeding, after medication channeling and between medications, or at least every 8 hours if the tube is not being used.

**NOTE:** Using less than the minimum recommended priming / flushing volume as indicated in the table may increase the likelihood of a tube clog during use, especially in cases of administration of multiple substances back to back. Use of very thick or sticky feeds or medications may require a greater flush volume to fully clear residue from the tubing. It is recommended that you speak with your healthcare professional regarding the best priming and flushing practice for your specific dietary and medication requirements.

- A smaller syringe may be used for flushing and priming the device. A 5 ml syringe is recommended for priming/flushing volumes less than 5 ml. Caution should be taken if using a syringe smaller than 30 ml. If using a smaller syringe and resistance is felt, stop flushing/priming immediately and switch to a larger syringe (30 to 60 ml). Continuing to flush/prime with a smaller syringe when resistance is felt can increase pressure on the tube and potentially damage smaller tubes.
- Use room temperature tap water for tube flushing. Sterile water may be appropriate where the quality of municipal water supplies is of concern. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain patency.
- Flush the feeding tube before and after medication channeling and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog. Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and to mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in warm water before channeling through the feeding tube. Never crush enteric-coated medication or mix medication with formula. Avoid using acidic irrigants such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging.
- Do not use excessive force to flush the tube. Excessive force can damage the tube and can cause injury to the gastrointestinal tract.
- Document the time and amount of water used in the patient's record. This will enable all caregivers to monitor the patient's needs more accurately.

	Minimum Recommended Priming/Displacement Volumes*		Minimum Recommended Flushing Volumes**	
	Gastric Port	Jejunal Port	Gastric Port	Jejunal Port
Direct Syringe	0.5 ml	0.5 ml	1.0 ml	1.0 ml
2" Extension + Device	1.0 ml	1.0 ml	2.0 ml	2.0 ml
12" Extension + Device	4.5 ml	2.5 ml	6.5 ml	4.5 ml

\*Priming/Displacement Volume is the amount of fluid needed to fill the entire length of the tubing and eliminate the air in the tubing.

\*\*Flushing Volume is the amount of fluid needed to clear the substances from the tube in order to avoid clogging and maintain tube patency.

**MEDICATION CHANNELING:** Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in water before channeling through the feeding tube. Never crush enteric-coated medication or mix medication with formula. Using a syringe, flush the tube with the prescribed amount of water.

**TUBING BLOCKAGE:** Obstructions are usually caused by formula accumulation inside the tube. Body fluids and medications can also clog the tube. Prevent this by thoroughly flushing the tube EVERY FOUR HOURS by injecting a prescribed amount of water or by injecting no less than the recommended flushing volume indicated in the table in the "General Flushing Guidelines" section. When you feel too much resistance and cannot inject water at all, the tube may need to be repositioned or replaced because it may be kinked. High pressure can damage the tube. Do not insert foreign objects into the tube. This may puncture the tube. If none of these measures work, call your specialist.

Tube occlusion is generally caused by: Poor flushing techniques • Failure to flush after measurement of gastric residuals • Inappropriate channeling of medication • Pill fragments • Viscous medications • Thick formulas, such as concentrated or enriched formulas that are generally thicker and more likely to obstruct tubes • Formula contamination that leads to coagulation • Reflux of gastric or intestinal contents up the tube

#### TO UNCLOG A TUBE

**WARNING: DO NOT USE EXCESSIVE FORCE OR PRESSURE WHEN CHANNELING FEED OR MEDICATION, OR TO ATTEMPT TO CLEAR A CLOG IN THE TUBING. THIS CAN CAUSE THE TUBING TO RUPTURE OR LEAD TO DAMAGE OF THE TUBING SUPPORT STRUCTURE. IF CLOG CANNOT BE CLEARED OR TUBING IS CLOGGING ON A FREQUENT BASIS, THIS MAY INDICATE THAT THE DEVICE NEEDS TO BE REPLACED. FAILURE TO FOLLOW THIS WARNING MAY RESULT IN DEVICE FAILURE AND/OR SERIOUS PATIENT INJURY (E.G. BOWEL IRRITATION OR PERFORATION).**

1. Make sure that the feeding tube is not kinked or clamped off.
2. If the clog is visible above the skin surface, gently massage or milk the tube between fingers to break up the clog.
3. Inspect the anti-reflux valve for signs of obstruction.
4. Next, place a 30-60ml syringe with warm water (use the recommended flushing volume indicated in the table in the "General Flushing Guidelines" section) into the appropriate adaptor or lumen of the tube and gently pull back on, and then depress the plunger to dislodge the clog.
5. If the clog remains, repeat step #4. Gentle suction alternating with syringe pressure will relieve most obstructions.
6. If this fails, consult with the physician. Do not use cranberry juice, cola drinks, meat tenderizer or chymotrypsin, as they can actually cause clogs or create adverse reactions in some patients. If the clog is stubborn and cannot be removed, the tube will have to be replaced.

**GRANULATION TISSUE:** A small amount of extra tissue around the stoma site is simply the body's attempt to heal the wound. This is normal. It is also normal for the extra tissue to secrete a small amount of mucus. If the tissue bleeds excessively, drains, or interferes with the external bolster, call your specialist.

**UNUSUAL SYMPTOMS:** Consult your health care practitioner if any of the following symptoms occur: Abdominal pain • Abdominal discomfort • Abdominal tenderness • Abdominal distension • Dizziness or fainting • Unexplained fever • Unusual amount of bleeding through or around the tube.

**DIARRHEA:** The flow rate may be too fast. Decrease the rate, then increase it slowly until you reach the prescribed flow rate. Diluting the formula with water may help slow the diarrhea, but check with your specialist first • Prepare fresh formula for each feeding. Discard unused formula after 24 hours • Keep your equipment clean. Prepare only enough formula for eight hours. Wash, rinse and dry the equipment after each use • If the diarrhea continues for more than three days despite your efforts, call your specialist • Never attempt to bolus feed through the jejunal lumen.

**STOMA CARE:** Cleanse the site with normal saline three times daily to remove the small amount of mucus that normally accumulates around the stoma. If the mucus dries, it may irritate the skin. Applying a small amount of water may loosen the dried material. After the stoma heals, a thorough cleansing with soap and water is best. A dressing is not necessary, and contributes to skin maceration and infection. Unless recommended by your specialist, avoid medicated ointments or powders. DO NOT use oil or petroleum based products. DO NOT ROTATE THE TUBE. Jejunal tubes must NOT be rotated because they will kink. Should the tube become kinked, call your health care professional.

**BALLOON MAINTENANCE:** It is recommended that balloon volume is checked at least every two weeks, or as recommended by your healthcare professional.

**NOTE:** Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

**Healthcare Professionals:** Insert a Luer slip syringe into the balloon inflation port (**Fig. A-5**) and withdraw the fluid while holding the tube in place. Compare the amount of water in the syringe to the amount recommended or the amount initially prescribed and documented in the patient record. If the amount is less than recommended or prescribed, refill the balloon with the water initially removed, then draw up and add the amount needed to bring the balloon volume up to the recommended and prescribed amount of water. Be aware as you deflate the balloon there may be some gastric contents that can leak from around the tube. Document the fluid volume, the amount of volume to be replaced (if any), the date and time. Wait 10-20 minutes and repeat the procedure. The balloon is leaking if it has lost fluid, and the tube should be replaced. A deflated or ruptured balloon could cause the tube to dislodge or be displaced. If the balloon is ruptured, it will need to be replaced. Secure the tube into position using tape, then follow facility protocol and/or call the physician for instructions.

**End Users:** It is recommended that balloon volume is checked at least every two weeks, or as recommended by your healthcare professional. While holding the external bolster in place, insert a luer slip syringe into the balloon inflation port (**Fig. A-5**). Withdraw the water from the balloon. The amount of water you withdraw should equal the prescribed amount. If you withdraw less water than the prescribed amount, replace the amount recommended by your specialist. Never add water before you withdraw all of the water from the balloon. **NOTE:** The use of a luer-lock syringe may make this more difficult. A luer-slip syringe is recommended.

#### DAILY CARE & MAINTENANCE CHECKLIST

##### • Assess the patient

- Assess the patient for any signs of pain, pressure or discomfort, warmth, rashes, purulent or gastrointestinal drainage.
- Assess the patient for any signs of pressure necrosis, skin breakdown or hypergranulation tissue.

##### • Clean the stoma site

- Use warm water and mild soap.
- Use a circular motion moving from the tube outwards. Clean sutures, external bolsters and any stabilizing devices using a cotton-tipped applicator.
- Rinse thoroughly and dry well.

##### • Assess the tube

- Assess the tube for any abnormalities such as damage, clogging or abnormal discoloration.

##### • Clean the feeding tube

- Use warm water and mild soap being careful not to pull or manipulate the tube excessively.
- Rinse thoroughly, dry well.

##### • Clean the jejunal, gastric and balloon ports

- Use a cotton tip applicator or soft cloth to remove all residual formula and medication.

##### • Do not rotate the external bolster

- This will cause the tube to kink and possibly lose position.

##### • Verify placement of the external bolster

- Verify that the external bolster rests 2-3 mm above the skin.

##### • Flush the feeding tube

- Flush the feeding tube with water every 4 hours during continuous feeding, anytime the feeding is interrupted, after every intermittent feeding, after medication channeling and between medications, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube after checking gastric residuals. Avoid using acidic irrigants such as cranberry juice and cola beverages to flush feeding tubes.

## DEVICE REPLACEMENT

**WARNING: REMOVAL IS RECOMMENDED TO BE PERFORMED BY A QUALIFIED CLINICIAN. SPONTANEOUS CLOSURE OF STOMA MAY OCCUR AS EARLY AS TWENTY-FOUR (24) HOURS AFTER REMOVAL. INSERT A NEW DEVICE IF ENTERAL FEEDING BY THIS ROUTE IS STILL INTENDED. IF CLOSURE IS DESIRED, APPLY A DRESSING OVER THE STOMA SITE.**

1. The AMT micro Transgastric-Jejunal Feeding Device should be replaced periodically for optimal performance. Frequent performance checks are recommended. Clogging and/or reduced flow are indicators of diminished performance.
2. The AMT micro Transgastric-Jejunal Feeding Device should only be used by or under the supervision of personnel trained in percutaneous gastrostomy tube placement. A thorough understanding of the technical principles, clinical applications, and risks associated with percutaneous gastrostomy tube placement is recommended prior to using this device.
3. Please contact a professional health care giver or physician for replacement of the AMT micro Transgastric-Jejunal Feeding Device.

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

## MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the AMT micro Transgastric-Jejunal Feeding Device is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only
- Maximum spatial gradient magnetic field of 1,000-gauss/cm (10-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg in the Normal Operating Mode

Under the scan conditions defined, the AMT micro Transgastric-Jejunal Feeding Device is expected to produce a maximum temperature rise of 2°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the AMT micro Transgastric-Jejunal Feeding Device extends approximately 20-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

## TROUBLESHOOTING

Long-term device performance and functionality are dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your feeding device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

<p><b>Anti-reflux valve leakage or blockage:</b> Leakage/blockage of the anti-reflux valve typically occurs due to residue (feed, medication, gastric contents, etc) becoming stuck in the valve area, preventing the valve from fully closing. Make sure device is flushed after each use. Leakage can also be caused due to excessive pressure in stomach (see <b>GASTRIC DECOMPRESSION</b> section). The valve can also become inverted in rare cases. Insert extension set into port to reset valve if this occurs. Avoid excessive pressure during channeling through the device, as this can cause a tear to occur in the valve over time.</p>
<p><b>Leaking of balloon volume:</b> If balloon is deflating, fully deflate the balloon. Compare the amount of water in the syringe to the amount recommended or the amount initially prescribed and documented in the patient record. Inflate balloon with recommended fill-volume. Do not exceed maximum fill volume. Wait 10-20 minutes and repeat the procedure. The balloon is leaking if it has lost fluid, and the tube should be replaced. Only access the balloon inflation port for balloon inflation and deflation. Other uses can result in balloon leakage or device failure. <b>NOTE:</b> The balloon is made of a semi-permeable material and can lose a small amount of fill-volume over time depending on environment and usage conditions.</p>
<p><b>A tear has formed:</b> Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.</p>
<p><b>Balloon inflation valve leakage:</b> Leakage from this valve typically occurs due to residue in the fill-valve. Always use a clean syringe when inflating the balloon and only inflate water through the valve. Leakage can also occur due to valve becoming stuck if syringe is inserted too hard into valve. Insert the syringe into the balloon inflation port if leakage is occurring to attempt to reset valve. Several attempts may be necessary before valve resets.</p>
<p><b>Foul smell coming from the device:</b> Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed and stoma site should be gently cleaned with soap and warm water. If foul smell does not go away, it is recommended that you contact your healthcare professional.</p>
<p><b>Plug will not stay closed:</b> Ensure that the plug is being firmly and fully pressed into the interlock connector. If the plug is not staying closed, check the plug and feed-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.</p>
<p><b>Device has become discolored:</b> The device can become discolored over days to months of use. This is normal depending on the types of feeds and medications being used with the device.</p>
<p><b>Difficulty decompressing the stomach:</b> Obstructions within the gastric lumen may reduce or prevent gastric decompression. Flush the gastric port EVERY FOUR HOURS by injecting a prescribed amount of water or by injecting no less than the recommended flushing volume indicated in the table in the "General Flushing Guidelines" section. If resistance is felt or water does not flow freely, it is recommended that you contact your healthcare professional (see <b>GASTRIC DECOMPRESSION</b> section).</p>
<p><b>Balloon will not inflate or deflate:</b> Inflation and deflation problems typically occur due to residue blocking the fill-lumen. Always use a clean syringe when inflating the balloon. In some rare occurrences, fungus can grow inside of the tubing and block the fill-lumen. Fungus growth can occur depending on patient environment and feeds/medications being channeled through the device. If balloon will not deflate, contact your healthcare professional for removal. If deflation problem is the result of fungal growth, eliminating source of fungal growth or anti-fungal medication may be required.</p>
<p><b>Balloon is misshaped:</b> Be sure to inflate and inspect balloon prior to placement to check balloon symmetry. Balloons becoming misshaped typically occur due to excessive force or tension on the device (pulling device out of stoma while balloon is inflated). This can occur if the device is too short for the patient's stoma site. Additionally, a device that is placed too close to the pylorus can result in a misshaped balloon and blockage of the pylorus. Balloons can be gently massaged with fingers back into symmetry if not excessively misshaped. Device should be replaced if balloon is found excessively misshaped.</p>
<p><b>Device fit is too tight or too loose:</b> Balloon fit can be adjusted by changing the balloon inflation volume within the balloon inflation range in <b>Table 1</b>. If balloon fit is too loose, increase balloon fill-volume above recommended but not above max fill-volume. If balloon fit is too tight, decrease fill-volume below recommended but not below min fill-volume. If device does not properly fit with the fill-volume range, a new device length will be needed.</p>
<p><b>Interlock failure or cracked:</b> The interlock has been designed to withstand extreme forces without detaching or cracking. However, the strength of the bond and material can reduce over prolonged use depending on medications and feeds used through the device. Excessive forces over an extended period of use can also reduce bond / material strength. The device should be replaced if interlock is found cracked, leaking, or separating from the device.</p>
<p><b>Balloon failure:</b> Early balloon failure can occur due to a number of patient or environmental factors, including but not limited to: gastric pH, diet, certain medications, improper balloon fill volume, placement of device, improperly channeling food/medication into balloon port, trauma, contact with a sharp or abrasive material, incorrect stoma length measurement, and overall care of the device.</p>
<p><b>Tubing has reduced flow or has become clogged:</b> Tubing can become blocked due to not properly flushing after each use, use of thick or improperly crushed medications, use of thick feeds/formulas, gastric reflux, and/or fungus growth. If clogged, refer to the <b>GENERAL USE AND CARE</b> section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced. <b>WARNING: DO NOT USE EXCESSIVE FORCE OR PRESSURE WHEN CHANNELING FEED OR MEDICATION, OR TO ATTEMPT TO CLEAR A CLOG IN THE TUBING. THIS CAN LEAD TO DAMAGE OF THE TUBING AND MAY RESULT IN FEEDING THE STOMACH. IF CLOG CANNOT BE CLEARED OR TUBING IS CLOGGING ON A FREQUENT BASIS, THIS MAY INDICATE THAT THE DEVICE NEEDS TO BE REPLACED. FAILURE TO FOLLOW THIS WARNING MAY RESULT IN DEVICE FAILURE AND/OR SERIOUS PATIENT INJURY (E.G. BOWEL IRRITATION OR PERFORATION).</b></p>

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<b>EN</b>	<b>STERILE EO</b>	Sterilized Using Ethylene Oxide	Rx Only	Prescription Only
 Do Not Resterilize	 Single Use Only	 Do not use if package is damaged and consult instructions for use	 MR Conditional	
 Not made with natural rubber latex	<b>MD</b> Medical Device	 Not made with DEHP (Di(2-ethylhexyl) Phthalate)		
ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.	Jejunal Port is MIC-KEY® feedset compatible			
	MIC-KEY® is a registered trademark of Avent, Inc.			
	The ENFit® connection is intended to reduce misconnections with other non-enteral applications.			

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MIC-KEY® on Avent Inc:in rekisteröity tavaramerkki
MIC-KEY® to zarejestrowany znak towarowy firmy Avent, Inc.
MIC-KEY® este o marcă înregistrată a Avent, Inc.
MIC-KEY® on Avent, Inc. registreeritud kaubamärk.
MIC-KEY® az Avent, Inc. bejegyzett márkaneve.
MIC-KEY® ir reģistrēta Avent, Inc. preču zīme.