AMT PEG-Assist® Initial Placement GI Kit

CE
2862

- ES** Estuche AMT PEG-Assist® de colocación inicial
- DE** AMT PEG-Assist®-Initialplatzierungskit
- IT** Kit di posizionamento iniziale AMT PEG-Assist®
- FR** Kit AMT PEG-Assist® de placement initial
- NO** AMT PEG-Assist® Innledende Plasseringssett
- SV** AMT PEG-Assist® initialinsättningsatts
- AR** طقم الجهاز الهضمي للموضعة الابتدائية AMT PEG-Assist®
- JA** AMT PEG-Assist® 初期配置GIキット
- PL** AMT PEG-Assist® Initial Placement GI Kit
- PT-BR** Conjunto de colocação inicial GI AMT PEG-Assist®
- RO** Trusă de plasare inițială GI AMT PEG-Assist®
- DA** AMT PEG-Assist® GI-sæt til førstegangsindsetting
- FI** AMT PEG-Assist® GI ensiasetus pakkaus
- EL** AMT PEG-Assist® GI Kit Αρχικής Τοποθέτησης
- ET** AMT PEG-Assist® GI-seadme esmapaigalduskomplekt
- PT** Kit GI de Colocação Inicial PEG-Assist® da AMT
- NL** AMT PEG-Assist® Initiële Plaatsing GI Kit
- BG** GI комплект за първоначално поставяне AMT PEG-Assist®
- HU** AMT PEG-Assist® kezdeti behelyezési GI készlet
- LV** AMT PEG-Assist® Sākotnējās ievietošanas gremošanas trakta komplekts
- LT** „AMT PEG-Assist®“ pirminio implantavimo GI rinkinys
- SK** Súprava AMT PEG-Assist® na úvodné GI zavedenie
- SL** GI komplet za začetno vstavitvev AMT PEG-Assist®
- TR** AMT PEG-Assist® İlk Yerleşirme GI Seti
- CS** GI souprava pro úvodní umístění AMT PEG-Assist®

R_x ONLY
STERILE EO

MD

KIT CONTENTS (Compatible with 12-18F Devices)

(1) Guidewire	(1) Guidewire Compatible Balloon Stoma Measuring Device
(2) Introducer Needle, 18 gauge	(1) Slip Tip Syringe
(1) #11 Scalpel	(2) Guidewire Compatible Introducer/Stiffener
(5) Dilators: 7, 12, 16, 18, 22F	(4) 4" x 4" (10 cm x 10 cm) Gauze

KIT CONTENTS (Compatible with 10F Devices)

(1) Guidewire	(1) Guidewire Compatible Balloon Stoma Measuring Device
(2) Introducer Needle, 18 gauge	(1) Slip Tip Syringe
(1) #11 Scalpel	(1) Guidewire Compatible Introducer/Stiffener
(3) Dilators: 7, 12, 14F	(4) 4" x 4" (10 cm x 10 cm) Gauze

SUPPLEMENTAL COMPONENTS (NOT INCLUDED IN THE KIT)

Endoscope	Sterile Water or Saline
Local Anesthetic	Needle Holders
Hemostats	U-Stitches (or alternate means of securement)

DESCRIPTION AND INDICATIONS FOR USE

CAUTION: Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

DESCRIPTION

The AMT Initial Placement Kit is specifically designed to improve upon the efficacy of the initial device placement procedure. Used with the procedure that follows, this initial placement kit has many unique features that can possibly:

- Eliminate the need for a follow-up in-office procedure for replacement of G-tube device with a low profile feeding device (LPFD).
- Allow the need for only one endoscope pass, allowing the endoscope to remain in the stomach until the entire procedure is completed.

INTENDED USE

The AMT Initial Placement Kit is an effective kit for use during the initial creation of a stoma for enteral feeding. The AMT Initial Placement Kits are intended to be used by clinicians for patients of all ages (infant to elderly).

INDICATIONS FOR USE

Initial placement of a percutaneous gastrostomy device may be indicated for patients with a functioning gut who require long-term feeding. This includes patients in whom malnutrition already exists, or may result, secondary to concurrent conditions.

CONTRAINDICATIONS

Contraindications for placement of a gastrostomy tube include but are not limited to: Nonfunctioning gut • Multiple surgical procedures near the GI site • Infection on or around the GI site • Inability to identify transabdominal illumination

WARNING: AFTER GI DEVICE PLACEMENT, PROPER POSITIONING OF THE BALLOON OR NON-BALLOON TYPE INTERNAL BOLSTER AGAINST THE GASTRIC MUCOSA MUST BE ENDOSCOPICALLY VERIFIED. TENSION ON THE DEVICE SHOULD BE AVOIDED IN ORDER TO MINIMIZE THE RISK OF COMPLICATIONS. DEVICE SHOULD BE CAPABLE OF ROTATING WITHOUT DIFFICULTY.

COMPLICATIONS

There are no complications associated with the AMT Initial Placement Kit at this time.

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device

CLINICAL BENEFITS

Clinical Benefits to be expected when using the AMT Initial Placement Kit include but are not limited to:

- Patients will benefit from the proper primary placement of a gastrostomy tube allowing for optimal channeling of nutrition and medication
- Reduces potential trauma associated with stoma tract development
- Clinicians will benefit by having most of the necessary tools to properly and easily place gastrostomy tube during initial placement surgeries
- Endoscope is able to remain in the stomach until the entire procedure is completed
- Guidewire compatibility helps to ensure the ease of proper placement of low profile feeding devices

PERFORMANCE CHARACTERISTICS

Performance Characteristics of the AMT Initial Placement Kit include but are not limited to:

- Provides a means of initiating a stoma tract in in-tact tissue
- Provides a means of stoma tract dilation to accommodate an interventional catheter and a means of introducing an interventional catheter to a freshly formed stoma tract
- Helps to ensure proper fit, comfort, and safety of gastrostomy feeding device
- Tapered dilators to help ensure proper dilation
- Available in different configurations to accommodate different stoma/device sizes

PROCEDURE PREPARATION

NOTE: Initial Placement is to be performed by a qualified clinician.


- Inspect the kit for damage. If damaged or sterile barrier is breached, do not use product.
- Use a clinically approved method to prep and sedate the patient for endoscopic procedure.

PROCEDURE PREPARATION

SELECTING AND PREPARING THE GASTROSTOMY SITE

- Once the patient is in the supine position and the endoscope is in place within the stomach, insufflate the stomach with air.
- Dim the room lights and deflect the endoscope to the anterior surface. The endoscope's light should clearly transilluminate through the anterior abdominal wall.
- Choose an insertion site that is free of major vessels, viscera, and scar tissue.

CAUTION: Attaining definitive transillumination at the insertion site is one of the most critical steps in the procedure. Do not continue if there is any question concerning an insertion site free of major vessels, viscera, and scar tissue.

- Depress the illuminated, intended insertion site. The depression on the anterior surface of the gastric wall should be clearly visible endoscopically.

CAUTION: Selection of the proper insertion site is a critical step in the procedure. Do not continue if there is inadequate visualization of the depression on the anterior surface of the gastric wall.

- Prep and drape the skin at the selected insertion site.

PLACEMENT PROCEDURE

1. SECUREMENT OF STOMACH TO ANTERIOR ABDOMINAL WALL

- At physician's discretion, an alternate method of securement may be necessary with larger adult patients due to greater thickness of the abdominal wall.
- Using two U-stitches inserted on opposite sides of the insertion site, secure the stomach against the anterior abdominal wall.
- Once the stitches are complete, take care to secure the two ends of each stitch in order to maintain the stomach against the abdominal wall.

2. GUIDEWIRE PLACEMENT

- While maintaining insufflation of the stomach and endoscopically visualizing the anterior mucosal surface, insert the 18-gauge introducer needle through the peritoneum and stomach wall.

NOTE: The most desirable insertion angle is perpendicular to the surface of the abdomen and stomach. Other restrictions may apply at the discretion of the physician.

CAUTION: It is imperative that one takes extreme caution when inserting the introducer needle to ensure the posterior wall and other vessels or viscera of the stomach are not perforated.

- When the 18-gauge needle is observed within the stomach, to minimize air loss through the needle, quickly cover the needle hub.
- Pass the enclosed guidewire through the 18-gauge needle sheath.

- Once the guidewire is observed in the stomach via the endoscope, remove the 18-gauge needle over the guidewire, taking caution not to dislodge the guidewire from the stomach.

NOTE: In the subsequent steps requiring use of the guidewire, do not pull up on the guidewire, dislodging it from the stoma.

3. STOMA TRACT CREATION AND STOMA MEASUREMENT

- Using the enclosed #11 scalpel, create a small skin incision that extends directly alongside the guidewire. If necessary, utilize the endoscope to verify that the incision extends downwards through the anterior gastric wall into the stomach.

- Utilizing the enclosed balloon stoma measuring device, obtain a stoma depth measurement. Slide the measuring disc up the tube - away from the balloon.

CAUTION: Selection of the correct size feeding device is critical for the safety and comfort of the patient. Measure the length of the patient's stoma with the stoma measuring device. The shaft length of the selected feeding device should be the same as the length of the stoma. An inappropriately sized feeding device can cause necrosis, buried bumper syndrome and/or hypergranulation tissue.

NOTE: It may be necessary to utilize a smaller size dilator prior to the measuring device in order to increase ease of insertion.

- Lubricate the measuring device and gently pass it over the guidewire and insert it through the incision until it is visible within the stomach.
- Insert a slip tip syringe containing 4ml AIR into the handle slip tip orifice and inflate balloon. Remove the syringe. Be careful not to put pressure on the valve area during the removal of the syringe, as this can cause immediate deflation.
- Pull gently on the Measuring Device until you can feel the balloon against the inside stomach wall. Slide the measuring disc down on the tube until it rests on the outside of the abdomen.
- Read the measurement value on the tube above the disc (side farther from balloon) to determine the stoma tract length. The measurement accuracy of the device is +/- 0.2cm.
- To deflate the balloon, place fingers on sides of the inflation valve. Gently squeeze to remove air from balloon. A noticeable sound may be heard when air escapes. It may be necessary to gently roll your fingers while squeezing to open the valve and deflate the balloon.

NOTE: Squeezing too hard will block the valve port, preventing balloon deflation. Only a partial squeeze - enough to open the valve - is all that is needed.






NOTE: Typical deflation using a slip tip syringe can also be used to deflate the balloon.

- Gently remove the Measuring Device.

NOTE: If measured stoma depth is between available feeding device sizes, utilize the larger size.

4. STOMA DILATION

- Using the enclosed dilator set, expand the stoma tract to the necessary diameter for chosen device size.
- Begin with a smaller size dilator. Lubricate surface of the dilator with a water-soluble lubricant.

CAUTION: It is imperative that one takes extreme caution when inserting dilators to ensure the posterior wall of the stomach is not perforated.

- Pass the dilator over the guidewire. Firmly but gently, use a rotating motion to twist the dilator into the stoma tract and slightly into the stomach. Utilize endoscope to visually verify dilator passage.
- Remove the first dilator using caution not to dislodge the guidewire. To minimize air loss through the stoma tract, hold a finger over the incision site.
- Repeat using larger dilators as necessary.
- Refer to **Chart 1** for recommended necessary dilation of French size for selected low profile device.

Chart 1				
Product French Size	10F	12F	14F	16-18F
Dilator	14F	16F	18F	22F

5. DEVICE PREPARATION

WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, RE-STERILIZE OR REPROCESS THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.

CAUTION: When using a low profile feeding device, do not use an undersized low profile device. Once placed, the device should rotate easily.

WARNING: UNDERSIZING THE DEVICE MAY CAUSE EMBEDDING WITH EROSION IN THE GASTRIC WALL, TISSUE NECROSIS, INFECTION, SEPSIS, AND ASSOCIATED SEQUELAE.

- If utilizing a balloon device, inflate the balloon prior to insertion to verify that it is fully functional. If the balloon is not functional, discard the product and obtain a new product.
- Lubricate the tip portion of the enclosed introducer stiffener using a water-soluble lubricant.

NOTE: One of the smaller dilators may also be used as a device stiffener.

- Insert the introducer stiffener into the device opening until the tip is flush with the tip of the device and necessary stiffening is achieved.
- Lubricate the outer surface of the low profile feeding device using water-soluble lubricant.

6. DEVICE PLACEMENT

- Pass the introducer stiffener over the guidewire.
- Insert the device into the stoma and then into the stomach. If necessary, use a firm but gentle rotating motion for insertion.
- If utilizing a balloon device, inflate the balloon through the balloon fill valve using the recommended amount of sterile water or saline and a slip tip syringe. Verify endoscopically that the entire internal bolster portion of the device is within the stomach.
- For non-guidewire compatible gastrostomy devices, the guidewire may need to be removed prior to placement. Follow gastrostomy device directions for use for method of insertion and internal bolster deployment.

WARNING: IF USING A BALLOON PRODUCT, ENDOSCOPICALLY GRIPPING OR GRASPING THE BALLOON OR TIP OF THE DEVICE MAY CAUSE DEVICE FAILURE.

- If utilizing a non-balloon device, release the bolster per manufacturing instruction.
- Verify endoscopically that internal bolster is retaining device against anterior surface of the mucosal wall.
- Rotate the device to ensure that it has a snug but not overly tight fit.

CAUTION: The inflated balloon or internal bolster should be capable of rotating freely. There should be no compression of either the gastric mucosa or the skin.

- Slowly withdraw the introducer stiffener and guidewire from the device.

CAUTION: Care must be taken when withdrawing the introducer stiffener to ensure that the entire device is not also removed.

- Tie the sutures over the device
- Deflate the stomach and remove the endoscope.

POSTOPERATIVE PROCEDURE AND PRECAUTIONS:

- It is advisable to wait until the following morning before initiating feedings. First feeding time, prescribed rate and volume are at physician's discretion.

NOTE: Kit contents can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

STERILE EO Sterilized Using Ethylene Oxide		Rx Only Prescription Only	
 Do Not Resterilize	 Single Use Only	 Do not use if package is damaged and consult instructions for use	
 Not made with natural rubber latex	 Medical Device	 Not made with DEHP (Di(2-ethylhexyl) Phthalate)	
PEG-Assist® is a registered trademark of Applied Medical Technology, Inc.			

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© 2024 Applied Medical Technology, Inc.
PEG-Assist® je registrovaná obchodní známka společnosti Applied Medical Technology, Inc.