# AMT G-JET®

## Low-Profile Transgastric-Jejunal Feeding Device Device Placement Directions for Use

Instructions for placement by hospital professionals only. For end user instructions please see [Tube Use and Care Instructions](#) document.

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
**Rx ONLY**  
**STERILE EO**

DEHP

**WARNING: DUE TO THE PRESENCE OF THE TUBING SUPPORT (SPRING) IN THE TUBING, DO NOT CUT THE DEVICE TO CUSTOM LENGTH. DO NOT stretch the tube as damage and/or delamination may occur between the tubing and the internal support structure.**

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Fig. D

\*NOTE: Device not included in Kit. Device sold separately. Refer to *Measuring the Stoma Length* section for instructions.


Fig. E

Fig. E-1



Fig. E-2

**CAUTION:** Components are supplied Sterile. For single user only. Do not re-use, re-sterilize, or reprocess this medical device. Doing so may compromise biocompatibility characteristics, device performance and/or material integrity; any of which may result in potential patient injury, illness and/or death.

For enteral nutrition and/or medication only.

#### KIT CONTENTS

(1) Low-Profile Transgastric-Jejunal Feeding Device	(1) Small Syringe (Luer Slip)
(1) Jejunal Feeding Set with Glow Green Right Angle Adapter	(1) Small Syringe with adapter (For direct port access) (AMT P/N: 4-7000)
(1) Gastric Feeding Set with White Right Angle Adapter	(1) Guidewire Introducer
(1) Large Syringe (Catheter or ENFit®)	(2) 4 "x 4" (10 cm x 10 cm) Gauze

#### SUPPLEMENTAL COMPONENTS (NOT INCLUDED IN THE KIT)

T-Fastener Set	Dilator Set (AMT P/N: IP-DIL)
.038" Compatible Introducer Needle	Guidewire Compatible Balloon Stoma Measurement Device (AMT P/N: 1804)
#11 Scalpel Blade	.038" Guidewire with Dispenser (AMT P/N: 4-GW-5)

**DESCRIPTION:** The AMT Transgastric-Jejunal feeding device (Fig A) provides for simultaneous gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum.

**INTENDED USE:** The AMT G-JET® is intended for gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. The AMT G-JET® is intended to be placed by qualified healthcare professionals. The AMT G-JET® is intended to be used by clinicians and caregivers/users.

**INDICATIONS FOR USE:** The AMT Transgastric-Jejunal feeding tube is indicated for use in patients who cannot absorb adequate nutrition through the stomach, who have intestinal motility problems, gastric outlet obstruction, severe gastroesophageal reflux, are at risk of aspiration, or in those who have had previous esophagectomy or gastrectomy. The use of this tube is also clinically indicated when simultaneous gastric decompression and jejunal feeding are needed. This includes patients in whom malnutrition already exists, or may result, secondary to concurrent conditions. The AMT Transgastric-Jejunal feeding tube is intended for use in adult, adolescent, child, and infant patients over 10kg.

**CONTRAINDICATIONS:** Contraindications for placement of a transgastric-jejunal feeding device include, but are not limited to ascites, colonic interposition, portal hypertension, peritonitis and morbid obesity.

**COMPLICATIONS:** The following complications may be associated with any transgastric-jejunal feeding device: Skin Break-down • Infection • Hypergranulation Tissue • Stomach or Duodenal Ulcers • Intraperitoneal Leakage • Pressure Necrosis • Jejunal Limb Migration • Bowel perforation • Tubing Support (Spring) Dislodgement • Intussusception

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

**CLINICAL BENEFITS:** Clinical Benefits to be expected when using the AMT G-JET® include but are not limited to: Provides direct access to the stomach and jejunum for feeding, decompression, or channeling medication • Device may be placed either during an initial placement procedure or as a replacement device • Four point stabilization of the external bolster reduces motion that can cause granulation tissue • Provides leakage protection via one-way valves to ensure patients retain nutrition and/or medication • Mutually exclusive Gastric and Jejunal ports help minimize the possibility of accidental feeding of the stomach, which could lead to choking and aspiration • Channeling of medication and feeding to the patient is immediate and sustains life • Made from medical-grade silicone to reduce skin irritation and improve patient comfort • Device comes in a variety of sizes for patient needs

**PERFORMANCE CHARACTERISTICS:** Performance Characteristics of the AMT G-JET® include but are not limited to: Low profile design sits flush against the skin • Unique internal tube structure – transitions from a tri-lumen design in its proximal (gastric) segment to one functional lumen in its distal (jejunal) segment • Jejunal lumen offers a larger cross-sectional area in the jejunal portion to help prevent clogging • Anti-kink technology (16F and 18F only) to help maintain flow • Tapered transition for a more secure fit • Unique "apple" shape balloon designed to reduce leakage and reduce the risk of accidental pull-outs from the stoma site • Guidewire adapter included in kit for ease of placement • Created with dual safety straps that allow independent access to the Gastric and/or Jejunal ports • Glow Green Jejunal port provides better visibility in the dark for easier for nighttime feeds

**DEVICE MATERIALS:** The AMT G-JET® is composed of the following materials: Medical-grade silicone (80%) • Medical-grade thermoplastic (13%) • Stainless steel spring (4%) • Medical-grade pad print ink (3%) • Stainless steel in tubing for 16 Fr and 18 Fr devices

## PLACEMENT PROCEDURE (FOR HEALTHCARE PROFESSIONALS)

**CAUTION:** Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

The AMT Transgastric-Jejunal feeding device may be placed percutaneously under fluoroscopic or endoscopic guidance or as a replacement to an existing device using an established stoma tract.

**CAUTION:** A gastrostomy must be performed to affix the stomach to the anterior abdominal wall, the feeding tube insertion site identified, stoma tract dilated and measured prior to initial tube insertion to ensure patient safety and comfort. The length of the tube should be sufficient to be placed beyond the ligament of treitz.

**CAUTION:** Do not use the retention balloon of the feeding tube as a gastrostomy device. The balloon may burst and fail to attach the stomach to the anterior abdominal wall.

## TUBE PREPARATION

1. Select the appropriate size Transgastric-Jejunal feeding device, remove from package and inspect for damage.

**NOTE:** The suture loop at the distal end of device is most often utilized during endoscopic placements. If the suture loop of the device is not necessary for the anticipated placement method, it is recommended to cut and gently remove the suture loop prior to placement. **DO NOT** use excessive force to remove the suture loop as this can damage tapered end of the tubing.

2. Inflate balloon through the balloon inflation port using a Luer slip syringe with distilled or sterile water to the recommended fill volume. (Fig A-5).

## TUBE PREPARATION

3. Remove the syringe and verify balloon integrity by gently squeezing the balloon to check for leaks. Visually inspect the balloon to verify symmetry. Symmetry may be achieved by gently rolling the balloon between the fingers. Reinsert the syringe and remove all the water from the balloon.
4. Using a Luer slip syringe, flush water through both the gastric and jejunal ports (**Fig A-1 & A-2**) to verify patency.

**WARNING: DO NOT CUT THE DISTAL END OF THE TUBING TO CREATE A CUSTOMIZED JEJUNAL LENGTH. DOING SO WILL ELIMINATE THE SOFT, TAPERED TIP OF THE DEVICE AND EXPOSE THE TUBING SUPPORT.**

5. Lubricate distal end of the tube with a water-soluble lubricant. Do Not Use Mineral Oil Or Petroleum Jelly.
6. Generously lubricate the jejunal lumen with water-soluble lubricant. Do Not Use Mineral Oil Or Petroleum Jelly.
7. Insert the introducer (**Fig E**) into the Jejunal port until the hub is in contact with the Jejunal feeding port (**Fig E-1**). The introducer opens the one-way valve and protects it from damage by the guidewire. When maneuvering the guidewire through the jejunal lumen and into the introducer, ensure that the tubing of the device is in-line with the introducer (**Fig E-2**). This will assist the guidewire as it passes through the introducer.

## SUGGESTED RADIOLOGIC PLACEMENT PROCEDURE

1. Place the patient in the supine position.
2. Prep and sedate the patient according to clinical protocol.
3. Insure that the left lobe of the liver is not over the fundus or the body of the stomach.
4. Identify the medial edge of the liver by CT scan or ultrasound.
5. Glucagon 0.5 to 1.0 mg IV may be administered to diminish gastric peristalsis.

**CAUTION: Consult glucagon instructions for use for rate of iv injection and recommendations for use with insulin dependent patients.**

6. Insufflate the stomach with air using a nasogastric catheter, usually 500 to 1,000 ml or until adequate distention is achieved. It is often necessary to continue air insufflation during the procedure, especially at the time of needle puncture and tract dilation, to keep the stomach distended so as to appose the gastric wall against the anterior abdominal wall.
7. Choose a catheter insertion site in the left sub-costal region, preferably over the lateral aspect or lateral to the rectus abdominis muscle (N.B. the superior epigastric artery courses along the medial aspect of the rectus) and directly over the body of the stomach toward the greater curvature. Using fluoroscopy, choose a location that allows as direct a vertical needle path as possible. Obtain a cross table lateral view prior to placement of gastrostomy when interposed colon or small bowel anterior to the stomach is suspected.

**NOTE:** PO/NG contrast may be administered the night prior or an enema administered prior to placement to pacify the transverse colon.

8. Prep and drape according to facility protocol.

## GASTROPEXY PLACEMENT

**CAUTION: It is recommended to perform a three point gastropexy in a triangle configuration to ensure attachment of the gastric wall to the anterior abdominal wall.**

1. Place a skin mark at the tube insertion site. Define the gastropexy pattern by placing three skin marks equidistant from the tube insertion site and in a triangle configuration.

**WARNING: ALLOW ADEQUATE DISTANCE BETWEEN THE INSERTION SITE AND GASTROPEXY PLACEMENT TO PREVENT INTERFERENCE OF THE T-FASTENER AND INFLATED BALLOON.**

2. Localize the puncture sites with 1% lidocaine and administer local anesthesia to the skin and peritoneum.
3. Place the first T-Fastener and confirm Intra-gastric Position. Repeat the procedure until all three T-Fasteners are inserted at the corners of the triangle.
4. Secure the stomach to the anterior abdominal wall and complete the procedure.

## CREATE THE STOMA TRACT

1. Create the stoma tract with the stomach still insufflated and in apposition to the abdominal wall. Identify the puncture site at the center of the gastropexy pattern. With fluoroscopic guidance confirm that the site overlies the distal body of the stomach below the costal margin and above the transverse colon.

**CAUTION: Avoid the epigastric artery that courses at the junction of the medial two-thirds and lateral one-third of the rectus muscle.**

**WARNING: TAKE CARE NOT TO ADVANCE THE PUNCTURE NEEDLE TOO DEEPLY IN ORDER TO AVOID PUNCTURING THE POSTERIOR GASTRIC WALL, PANCREAS, LEFT KIDNEY, AORTA OR SPLEEN.**

2. Anesthetize the puncture site with local injection of 1% lidocaine down to the peritoneal surface.
3. Insert a .038" compatible introducer needle at the center of the gastropexy pattern into the gastric lumen directed toward the pylorus.

**NOTE:** The best angle of insertion is a 45 degree angle to the surface of the skin.

4. Use fluoroscopic visualization to verify correct needle placement. Additionally, to aid in verification, a water filled syringe may be attached to the needle hub and air aspirated from the gastric lumen.

**NOTE:** Contrast may be injected upon return of air to visualize gastric folds and confirm position.

5. Advance a guidewire, up to .038", through the needle and coil in the fundus of the stomach. Confirm position.
6. Remove the introducer needle, leaving the guidewire in place and dispose of according to facility protocol.
7. Advance a .038" compatible flexible catheter over the guidewire and using fluoroscopic guidance, manipulate the guidewire into the antrum of the stomach.
8. Advance the guidewire and flexible catheter until the catheter tip is at the pylorus.
9. Negotiate through the pylorus and advance the guidewire and catheter into the duodenum and beyond the Ligament of Treitz.
10. Remove the catheter and leave the guidewire in place.

## DILATION

1. Use a #11 scalpel blade to create a small skin incision that extends alongside the guidewire, downward through the subcutaneous tissue and fascia of the abdominal musculature. After the incision is made, dispose of according to facility protocol.
2. Advance a dilator over the guidewire and dilate the stoma tract to the desired size.
3. Remove the dilator over the guidewire, leaving the guidewire in place.
4. Measure the Stoma Length with the AMT Stoma Measuring Device.

## MEASURING THE STOMA LENGTH

**CAUTION:** Selection of the correct size feeding device is critical for the safety and comfort of the patient. Measure the length of the patient's stoma with the stoma measuring device (fig D). The shaft length of the selected feeding device should be the same as the length of the stoma. An inappropriately sized feeding device can cause necrosis, buried bumper syndrome and/or hypergranulation tissue.

1. Inspect Measuring Device (Fig D) and verify balloon inflates prior to using. If any damage is noticed, discard and obtain another device for use.
2. Slide the measuring disc up the tube - away from the balloon.
3. Lubricate the balloon tip of the Measuring Device using a water soluble lubricant.
4. Gently slide the tube portion of the Measuring Device through the stoma and into the stomach.

**NOTE:** A guidewire can optionally be used to aid in guidance through the stoma tract. Measuring Device will slide over guidewires up to 0.038" in diameter.

5. Insert a slip tip syringe containing 4ml AIR into the handle side orifice and inflate balloon. Remove the syringe.
6. Pull gently on the Measuring Device until you can feel the balloon against the inside stomach wall.
7. Slide the measuring disc down on the tube until it rests on the outside of the abdomen.
8. Read the measurement value on the tube above the disc (side farther from balloon) to determine the stoma tract length. The measurement accuracy of the device is +/- 0.2cm.
9. To deflate the balloon, the thumb should be placed directly below inflation port, while index finger is placed directly above inflation port. Gently squeeze to remove air from balloon. A noticeable sound may be heard when air escapes.

**NOTE:** Squeezing too hard will block the valve port, preventing balloon deflation. Only a partial squeeze - enough to open the valve - is all that is needed.

**NOTE:** Typical deflation using a slip tipped syringe can also be used to deflate the balloon.

10. Gently remove the Measuring Device.

## TUBE PLACEMENT FOR RADIOLOGIC PROCEDURE

**NOTE:** A peel-away sheath may be used to facilitate advancement of the tube through the stoma tract.

1. Select the appropriate AMT Transgastric-Jejunal feeding device and prepare according to the directions in the "Tube Preparation" section above.

2. Advance the distal end of the tube over the guidewire until the proximal end of the guidewire exits the introducer.

**NOTE:** Direct visualization and manipulation of the introducer and guidewire may be required to pass the guidewire through the end of the introducer.

3. Hold the introducer hub and jejunal feeding port while advancing the tube over the guidewire and into the stomach.
4. Rotate the AMT Transgastric-Jejunal feeding device while advancing to facilitate passage of the tube through the pylorus and into the jejunum.
5. Advance the tube until the tip of the tube is beyond the Ligament of Treitz and the balloon is in the stomach.
6. Ensure the external bolster is flush with the skin.
7. Using a Luer slip syringe, inflate the balloon according to the **BALLOON INFLATION VOLUME** section of the *Tube Use and Care Instructions*.
8. Remove the guidewire through the introducer while holding the introducer in position. Remove the introducer.

## VERIFY TUBE POSITION

1. Verify proper tube placement radiographically to avoid potential complication (e.g. bowel irritation or perforation) and ensure the tube is not looped within the stomach or small bowel.

**NOTE:** The jejunal portion of the tube contains a radiopaque material and can be used to radiographically confirm position. Do not inject contrast into the balloon.

2. Using the small syringe with winged adapter, inject contrast dye directly into the gastric and jejunal ports to verify patency and proper tube position (Fig A-6).
3. Check for moisture around the stoma. If there are signs of gastric leakage, check the tube position and external bolster placement. Add fluid as needed in 1-2 ml increments.

**CAUTION:** Do not exceed total balloon volume indicated.

4. Assure that the external bolster is not placed too tightly against the skin and rests 2-3 mm above the abdomen.
5. Document the date, the type, the size and lot number of the tube, the fill volume of the balloon, skin condition and patient tolerance to the procedure. Start feeding and medication administration per physician orders and after confirmation of proper tube placement and patency.

## RADIOLOGIC PLACEMENT THROUGH AN ESTABLISHED GASTROSTOMY TRACT

1. Insert a floppy-tipped guidewire (up to .038" diameter) through the existing gastrostomy tube under fluoroscopic guidance. The distal end of the guidewire should be visible within the stomach.
2. Remove the existing gastrostomy tube over the guidewire, leaving the distal end of the guidewire partially coiled within the stomach.
3. Using the AMT Guidewire Compatible Balloon Stoma Measuring Device, properly measure the length of the stoma according to the instructions outlined in the "Measuring the Stoma Length" section above.
4. After the measuring device is removed, navigate the guidewire through the pylorus and advance the guidewire into the duodenum.
5. Refer to the instructions outlined in the sections entitled "Tube Placement for Radiologic Procedure" and "Verify Tube Position" to complete device placement.

## SUGGESTED ENDOSCOPIC PLACEMENT PROCEDURE

1. Perform routine Esophagogastroduodenoscopy (EGD). Once the procedure is complete and no abnormalities are identified that could pose a contraindication to placement of the tube, place the patient in the supine position and insufflate the stomach with air.
2. Transilluminate through the anterior abdominal wall to select a gastrostomy site that is free of major vessels, viscera and scar tissue. The site is usually one third the distance from the umbilicus to the left costal margin at the midclavicular line.
3. Depress the intended insertion site with a finger. The endoscopist should clearly see the resulting depression on the anterior or surface of the gastric wall.
4. Prep and drape the skin at the selected insertion site.
5. Place the device according to the directions in the "Gastropexy Placement" section.

## CREATE THE STOMA TRACT AND DILATION

1. Create the stoma tract with the stomach still insufflated and in apposition to the abdominal wall. Identify the puncture site at the center of the gastropexy pattern. With endoscopic guidance confirm that the site overlies the distal body of the stomach below the costal margin and above the transverse colon.

**CAUTION: Avoid the epigastric artery that courses at the junction of the medial two-thirds and lateral one-third of the rectus muscle.**

**WARNING: TAKE CARE NOT TO ADVANCE THE PUNCTURE NEEDLE TOO DEEPLY IN ORDER TO AVOID PUNCTURING THE POSTERIOR GASTRIC WALL, PANCREAS, LEFT KIDNEY, AORTA OR SPLEEN.**

2. Anesthetize the puncture site with local injection of 1% lidocaine down to the peritoneal surface.
3. Insert a .038" compatible introducer needle at the center of the gastropexy pattern into the gastric lumen directed toward the pylorus.

**NOTE:** The best angle of insertion is a 45 degree angle to the surface of the skin.

4. Use endoscopic visualization to verify correct needle placement.
5. Advance a guidewire, up to .038", through the needle into the stomach. Using endoscopic visualization, grasp the guidewire with atraumatic forceps.
6. Remove the introducer needle, leaving the guidewire in place and dispose of according to facility protocol.
7. Dilate the stoma tract according to the directions in the "Dilation" section listed above.

## TUBE PLACEMENT FOR ENDOSCOPIC PROCEDURE

1. Measure the stoma length according to the directions in the "Measuring the Stoma Length" section listed above.
2. Select the appropriate sized AMT Transgastric-Jejunal feeding device and prepare according to the directions in the "Tube Preparation" section listed above.
3. Advance the distal end of the tube over the guidewire until the proximal end of the guidewire exits the introducer.

**NOTE:** Direct visualization and manipulation of the introducer and guidewire may be required to pass the guidewire through the end of the introducer.

4. Hold the introducer hub and jejunal port while advancing the tube over the guidewire and into the stomach.
5. Using endoscopic guidance, grasp the suture loop or the tip of the tube with atraumatic forceps.
6. Advance the AMT Transgastric-Jejunal feeding device through the pylorus and upper duodenum. Continue to advance the tube using the forceps until the tip is positioned beyond the Ligament of Treitz and the balloon is in the stomach.
7. Release the tube and withdraw the endoscope and forceps in tandem, leaving the tube in place.
8. Ensure that the external bolster is flush with the skin.
9. Using a Luer slip syringe, inflate the balloon according to the **BALLOON INFLATION VOLUME** section of the *Tube Use and Care Instructions*.
10. Remove the guidewire through the introducer while holding the introducer in place.
11. Remove the introducer.
12. Verify the tube position according to the directions in the "Verify Tube Position" section listed above.

## ENDOSCOPIC PLACEMENT THROUGH AN EXISTING GASTROSTOMY TRACT

1. Following established protocol, perform routine Esophagogastroduodenoscopy (EGD). Once the procedure is complete and no abnormalities are identified that could pose a contraindication to placement of the tube, place the patient in the supine position and insufflate the stomach with air.
2. Manipulate the endoscope until the indwelling gastrostomy tube is in the visual field.
3. Insert a floppy-tip guidewire through the indwelling gastrostomy tube and remove the tube.
4. Place tubing according to the "Tube Placement for Endoscopic Procedure" section above.
5. Verify proper tube placement according to prior instructions in "Verify Tube Position" section.

**WARNING: DO NOT STRETCH THE TUBE AS DAMAGE AND/OR DELAMINATION MAY OCCUR BETWEEN THE TUBING AND THE INTERNAL SUPPORT STRUCTURE.**

**NOTE:** The device and the other kit components can be disposed of by following local disposal guidelines or through facility protocol.

**Jejunal Port is MIC-KEY® feedset compatible**

MIC-KEY® is a registered trademark of Avent, Inc.

ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.

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**Applied Medical Technology, Inc.**