# MiniACE® Button

## Low Profile Antegrade Continence Enema Button

CE  
2862

- ES** MiniACE® Botón. Botón de Enema de Continencia Anterógrado de Bajo Perfil.
- DE** MiniACE® Button. Flacher, antegrader Kontinenzeinlauf-Button.
- IT** MiniACE® Bottone. Pulsante per clistere di continenza anterograda a basso profilo.
- FR** Bouton MiniACE®. Bouton de lavement de continence antérograde à profil bas.
- NO** MiniACE® Knapp. Lavprofil-antegrads inkontinensklyster-knapp.
- SV** MiniACE® Knapp. Lågprofilsknapp för antegrad kontinensirrigation.
- PT-BR** Botão MiniACE®. Botão de Enema de Continência Anterógrada de Baixo Perfil.
- NL** MiniACE®-Knop. Vooruitbewegende continentiëklysmaknop met laag profiel.
- DA** MiniACE® Knap. Fremadrettet kontinenslavematknop med lav profil.
- FI** MiniACE®-Nappi. Matalan profiilin antegradinen huhtelunappi.
- AR** MiniACE® زر. حصر الحقنة الشرجية التقدمي غير المرئي.
- BG** Бутон MiniACE®. Бутон за нископрофилна антеградна задържаща клизма.
- CS** MiniACE® Knoflík. Knoflík s nízkým profilem pro kontinentní nálevovou apendikocékostomii.
- ET** MiniACE® Nupp. Low Profile Antegrade Continence Enema Button (Madalaprofiililine ACE nupp).
- EL** MiniACE® Κομπτί. Κομπτί κλύσματος συγκράτησης ορθοδομής χαμηλού προφίλ.
- HU** MiniACE® Gomb. Alacsony profilú Antegrade continencia beöntési gomb.
- JA** MiniACE® ボタン。低プロファイル順行性禁制洗腸ボタン。
- LV** MiniACE® Button. Zema profila anterogrādās kontinences klizmas Button ierīce.
- LT** „MiniACE®“ Sagutė. Mažo profilio išmatų nelaikymo klizmos sagutė.
- PT** Botão MiniACE®. Botão Enema para Continência Anterógrada de Baixo Perfil.
- PL** Zgłębnik MiniACE® Button. Niskoprofilowy zstępujący zgłębnik Button do lewatywy.
- RO** Butonul MiniACE®. Buton pentru irigare continentă antegradă cu profil scăzut.
- SK** MiniACE® Tlačidlo. Nízko profilové tlačidlo klystíra antegrádnej kontinencie.
- SL** MiniACE® Button. Antegradni kontinenčni gumb z nizkim profilom za klistiranje.
- TR** MiniACE® Butonu. Düşük Profilli Antegrad Kontinans Lavman Butonu.

STERILE EO

MD

Rx ONLY


DEHP

C4628-I 09/2024

Made in the U.S.A.

Figure 1: MiniACE® Button  
 Figura 1: MiniACE® Butón  
 Abbildung 1: MiniACE® Button  
 Figura 1: MiniACE® Bottone  
 Figure 1: Bouton MiniACE®  
 Figure 1: MiniACE®-knapp  
 Figure 1: MiniACE®-knapp  
 Figura 1: Botão MiniACE®  
 Afbeelding 1: MiniACE®-knop  
 Figur 1: MiniACE®-knapp  
 Kuva 1: MiniACE®-nappi  
 MiniACE®: 1 الشغل  
 Фигура 1: Бутон MiniACE®

Obrazek 1: MiniACE® knoflik  
 Joonis 1: MiniACE® Button  
 Εικόνα 1: MiniACE® Κομπτί  
 1. ábra: MiniACE® gomb  
 図1: MiniACE®ボタン  
 1. att: MiniACE® Button  
 1 pav. „MiniACE®“ sagutė  
 Figura 1: Botão MiniACE®  
 Rysunek 1: MiniACE® Button  
 Figura 1: Butonul MiniACE®  
 Obrazok 1: MiniACE® Tlačidlo  
 Slika 1: MiniACE® Button  
 Şekil 1: MiniACE® Butonu

1. Balloon Inflation Port  
 Puerto de inflado del balón  
 Ballon-Befüllungs port  
 Porta per gonfiare il palloncino  
 Ballonnet Orifice de gonflage  
 Ballong-oppblåsnings-port  
 Ballongupplåsningsport  
 Porta de enchimento do balão  
 Ballon Inflatiepoort  
 balloninflaterings-port  
 Pallon täyttöportti  
 مدخل نفخ البالون  
 Балон Порт за запълване на балона

Balónek Nafukovací otvor  
 Balloon Täiteport  
 Μπαλόνι Θύρα διόγκωσης  
 Ballon Felfújó nyílás  
 バルーン 膨張ポート  
 Balons Piepūšanas ports  
 Balionėlis Pripūtimo jungtis  
 Porta de Insuflação do Balão  
 Balon Port do napeiniania  
 Balon Port umflare  
 Balónik naplnenie balónika  
 Balonček Inflačijski port  
 Balon Şişirme Girişi


2. Tubing  
 Tubo  
 Ballon-Befüllungs port  
 Sonda  
 Tubulure  
 Slange  
 Rör  
 Tubo  
 Slangen  
 Sonde  
 Letku  
 الأنبوب  
 Тръби  
 Trubice

Torud  
 Σωλήνωση  
 Csővezés  
 チューブ  
 Caurule  
 Vamzdeliai  
 Sonda  
 Przewód  
 Tub  
 Hadičky  
 Cevka  
 Түр

3. Interlock (feed/med port)  
 Enclavamiento (puerto alimentación/medicación)  
 Kupplung (Port zur Verabreichung von Nahrung/  
 Medikamenten)  
 Interblocco (porta nutr./farm.)  
 Verrouillage (orifice d'alimentation/médication)  
 Forrigling (mat/med. port)  
 Spärrias (närlings/med port)  
 Engate (porta de alim./medic.)  
 Interlock (voedings-/med.port)  
 Interlock (ern/med-port)  
 Lukkoliitin (ruokinta/lääkitysportti)  
 وشاحة (مدخل التغذية / الدواء)  
 Заключване (порт за хранене/лекарства)


Západka (otvor vyživování/léky)  
 Lukustus (toidu/ravimi port)  
 Μαντάλωμα (θύρα σίτισης/ενδίαμεση)  
 Reteszelés (tápláló/gyógyszernyílás)  
 連結 (栄養/薬剤ポート)  
 Starpsavienojums (barošanas/med. ports)  
 Perjungėjas (maistimo / vaistų jungtis)  
 Interlock (porta de alimentação/méd.)  
 Blokada (port do żywienia/podawania leków)  
 Blocare (port hrănire/med)  
 Zámok (port na podávanie výživy/liekov)  
 Zaklep (port za dovajanje hrane/zdravil)  
 Güvenlik kilidi (beslenme/ilac girişi)

|   |  |   |   |   |  |
|---|--|---|---|---|--|
| <p>4. Strap Plug<br/>Tapón de correa<br/>Stöpsel<br/>Tappo copri porta<br/>Bouchon de sangle<br/>Stropp-plugg<br/>Fästpropp<br/>Bujão da alça<br/>Aftakningsplug<br/>Fjederpropp<br/>Tulppa<br/>سدادة بطوق<br/>Лентов щепсел<br/>Zátka s páskem<br/>Rihma pistik<br/>Βύσμα ιμάντα</p> | <p>Heveder dugó<br/>ストラップ・プラグ<br/>Lences spraudnis<br/>Dirželio kaištis<br/>Plugue de alça<br/>Korek paska<br/>Dop cu bandă<br/>Uzávěr s popruhom<br/>Vtič za priključitev<br/>na trak<br/>Aski Tampon</p> | <p>5. Anti-Reflex Valve<br/>Válvula antirreflujo<br/>Anti-Rückflussventil<br/>Valvola antireflusso<br/>Valve anti-reflux<br/>Anti-tilbakeslagsventil<br/>Antireflux ventil<br/>Válvula antirrefluxo<br/>Antirefluxklep<br/>Antitilbageløbsventil<br/>Refluxinesteventiili<br/>صمام مضاد للارتداد<br/>Клапан против<br/>рефлюкс<br/>Antirefluxní ventil<br/>Ühesuunaline ventiil</p> | <p>Αντίρροφη βαλβίδα<br/>Refluxgátó szelep<br/>逆止め弁<br/>Anti-refluxus vārsts<br/>Antirefluxinis vožtuvas<br/>Válvula antirrefluxo<br/>Zavór antyrefluxowsy<br/>Supapán anti-reflux<br/>Ventil proti refluxu<br/>Protirefluxusni ventil<br/>Anti Reflüj Kapağı</p> | <p>6. Balloon<br/>Balón<br/>Ballon<br/>Palloncino<br/>Ballonnet<br/>Ballong<br/>Ballong<br/>Ballon<br/>Ballon<br/>Pallo<br/>بالون<br/>БАЛОН<br/>Balónek<br/>Balloun<br/>Μπαλόνι</p> | <p>Ballon<br/>バルーン<br/>Balons<br/>Balionėlis<br/>Balão<br/>Balon<br/>Balon<br/>Balão<br/>Balónik<br/>Balonček<br/>Balon<br/>Ballon<br/>Μπαλόνι</p> |
|---|--|---|---|---|--|

**Figure 2:** 12" RA Irrigation Set  
**Figura 2:** Juego de Irrigación de Ángulo ecto de 12 pulgadas  
**Abbildung 2:** 12"-Irrigationsset rechtwinklig  
**Figure 2:** Set di irrigazione ad angolo retto da 12 pollici  
**Figure 2:** 12" Kit d'irrigation à angle droit de 12 pouces  
**Figure 2:** 12-tommers høyrevinklet irrigasjonsett  
**Figure 2:** 12" Right Angle spolenhet  
**Figure 2:** Conjunto de Irrigação de Ângulo Reto de 12 Pol.  
**Abbeildung 2:** Irrigationsset met rechte hoek van 12 inch  
**Figure 2:** 12" høyrevinklet udsykningsssæt  
Kuva 2: 30 cm huuhetuletku kulmaliittimellä  
الشكل 2: مجموعة ارواء 12 بوصة قائمة الزاوية  
**Фигура 2:** 12 инчов комплект за иригация с десен ъгъл  
**Obrazek 2:** 12 palcová pravouhlá  
vyplachovací sada

**Joonis 2:** 12-tolline täisnurk-irrigatsioonimoodul  
**Εικόνα 2:** 12 Ιντσών Δεξιάς γωνιάς σετ έκπλυσης  
**2 ábra:** 12" derékszögű irrigáló készlet  
**图2:** 12英寸直角洗淨セット  
**2 att:** 12 collu taisna lenġa iriġċijas komplekts  
**2 pav.** 12 col. dešiniojo kampo iriġacijos rinkinys  
**Figure 2:** Conjunto de Rego Ângulo Direito de 12 Polegadas  
**Rysunek 2:** 12-calowy prawy zestaw do irygacji  
**Figure 2:** Set de irigare cu unghi drept de 12 inci  
**Obrazok 2:** 12 palcová pravouhlá zavlažovacia súprava  
**Slika 2:** 12-palčni pravokotni komplet za izpiranje  
**Şekil 2:** 12 İnc Dik Açılı İrriğasyon Seti

**Figure 3:** Balloon Inflation Port  
**Figura 3:** Puerto de inflado del balón  
**Abbildung 3:** Ballon-Befüllungs-Port  
**Figure 3:** Porta per gonfiare il palloncino  
**Figure 3:** Ballonnet Orifice de gonflage  
**Figure 3:** Ballong-opplåsningsport  
**fig 3:** Ballongupplåsningsport  
**Figure 3:** Porta de enchimento do balão  
**Abbeildung 3:** Ballon Inflatieport  
**Figure 3:** balloninflaterings-port  
**Kuva 3:** Pallon täyttöportti  
الشكل 3: مدخل نفخ البالون  
**Фигура 3:** Балон Порт за запълване на балона

**Obrazek 3:** Balónek Nafukovací otvor  
**Joonis 3:** Balloon Täiteport  
**Εικόνα 3:** Μπαλόνι Ούρα διόγκωσης  
**3. ábra:** Ballon Felfújó nyílás  
**图3:** 气球膨胀气口  
**3. att.** Balons Pierūšanas ports  
**3 pav.** Balionėlis Pripūtimo jungtis  
**Figure 3:** Porta de Insuflação do Balão  
**Rysunek 3:** Balon Port do napełniania  
**Figure 3:** Balon Port umflare  
**Obrazok č. 3:** Port na naplnění balónika  
**Slika 3:** Balonček Inflačijski port  
**Şekil 3:** Balon Şişirme Girişi

# MiniACE<sup>®</sup> Button

## Low Profile Antegrade Continence Enema Button

### INSTRUCTIONS FOR USE

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. **Notice:** Peel off ID sticker from tray and store for future use in patient chart or other convenient location. The device is provided STERILE for single use. Please inspect all contents of the kit for damage. If damaged, do not use.

### INTENDED USE

The MiniACE<sup>®</sup> Button is intended to be used as an effective device for channeling irrigation fluid through a secured (initial placement) or formed (replacement) stoma. The MiniACE<sup>®</sup> Button is intended to be used by clinicians and trained patients/caregivers.

### INDICATIONS FOR USE

The Bowel Management Device is intended to instill fluids through a stoma into the colon to promote evacuation of the contents of the lower bowel through the anus and is intended to be an aid in the management of fecal incontinence. The catheter is placed and maintained in a percutaneously prepared opening, such as a cecostomy or appendicostomy. The Bowel Management Device is intended to be used in children and adults.

### CONTRAINDICATIONS

Contraindications for placement of the MiniACE<sup>®</sup> Button include but are not limited to:

**Initial Placement:** Colonic interposition • Ascites • Portal hypertension • Peritonitis • Uncorrected coagulopathy • Infection around stoma site • Uncertainty as to cecostomy tract direction and length (colon wall thickness)

**Replacement:** Lack of adherence of the cecum/appendix/colon to the abdominal wall • Lack of a well established cecostomy site • Evidence of infection • Uncertainty as to cecostomy tract direction and length (colon wall thickness) • Presence of multiple stoma fistulous tracts

**NOTICE:**

- Please contact a professional health care giver or physician for explanation of the warnings, care, and use of the device.
- The MiniACE<sup>®</sup> Button has been designed to provide irrigation access into the colon. Other applications are not advised.

### COMPLICATIONS

Potential complications when using the MiniACE<sup>®</sup> Button include but are not limited to:

Peristomal pain • Abscess, wound infection and skin breakdown • Pressure necrosis • Hypergranulation tissue • Intraperitoneal leakage • Buried bumper syndrome • Peristomal leakage • Balloon failure or dislodgement • Tube clog • Colonic bleeding and/or ulcerations • Bowel volvulus

**WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, RE-STERILIZE OR REPROCESS THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.**

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

### CLINICAL BENEFITS

Clinical Benefits to be expected when using the MiniACE<sup>®</sup> Button include but are not limited to:

- Facilitates antegrade continence enemas for patients who are in need of flushing the bowels
- Cleared for use in both Appendicostomy (MACE) and Cecostomy procedures
- Device may be placed either during an initial placement procedure or as a replacement device
- Device can be easily exchanged in the exam room without anesthesia
- Balloon design aids in easy, at home replacement when necessary
- Patient-specific sizing minimizes stoma dilation and leads to less movement within the stoma, which may minimize pain and discomfort at the stoma site
- Ease of use and replacement may improve quality of life

## PERFORMANCE CHARACTERISTICS

Performance Characteristics of the MiniACE<sup>®</sup> Button include but are not limited to:

- Unique “apple” shape balloon designed to reduce leakage and reduce the risk of accidental pull-outs from the stoma site
- Anti-leakage valve prevents backflow of irrigation fluids
- Protected tip is gentle on intestinal walls
- Balloon fill volume can be easily adjusted for a custom-fit
- Low-profile device sits flush against the skin
- Soft external bolster
- Interlock allows safe and secure attachment of the irrigation set

## DEVICE MATERIALS

The MiniACE<sup>®</sup> Button is made of the following materials: Medical-grade silicone (81%) • Medical-grade thermoplastic (16%) • Stainless steel spring (2%) • Medical-grade silicone pad print ink (1%)

## KIT CONTENTS

- (1) MiniACE<sup>®</sup> Button (Fig. 1)
- (2) 4 x 4 inch Gauze
- (1) Introducer/Stiffener
- (1) Water Soluble Lubricant
- (1) Luer-Slip Syringe (5ml)
- (1) 12" Right Angle Irrigation Set (Fig. 2)

### **Optional Accessories for Device Placement (Not Included):**

Dilators, Stoma Measuring Device, Introducer Needle, Scalpel, Guide-wire.

### **Additional Components for Use (Not Included):**

Replacement Irrigation Sets, Irrigation Bag

## TYPE OF USAGE – INITIAL PLACEMENT VS REPLACEMENT

The MiniACE<sup>®</sup> Button may be placed either during an initial placement procedure or as a replacement device.

### **PERFORMING AN INITIAL PLACEMENT**

If the patient does not currently have a stoma tract for placement of the MiniACE<sup>®</sup> Button device, a new stoma tract will need to be created. This process can only be completed by a healthcare professional per proper cecostomy/appendicostomy surgical procedures.

**WARNING: INITIAL PLACEMENT OF THE MINIACE<sup>®</sup> BUTTON REQUIRES THAT A PROCEDURE IS PERFORMED TO AFFIX THE INTESTINAL WALL TO THE ANTERIOR ABDOMINAL WALL. DO NOT USE THE RETAINING BALLOON OF THE DEVICE AS AN ATTACHMENT/ANCHORING DEVICE. AN EARLY BALLOON FAILURE MAY PREVENT THE INTESTINAL WALL FROM ATTACHING TO THE ANTERIOR ABDOMINAL WALL.**

**CAUTION: It is recommended to perform a three point anchoring/securement in a triangle configuration to ensure attachment of the intestinal wall to the anterior abdominal wall.**

### **REPLACING A DEVICE IN AN ESTABLISHED STOMA SITE**

If it is determined that a device placed in an existing stoma needs replaced (due to non-optimal functionality or as a pro-active replacement), the current device can be removed from the stoma and the new device can be inserted in the same stoma site. If performing a device replacement, skip straight to the **DEVICE PLACEMENT PROCEDURE** section for the proper method of placing the new device.

**NOTE:** Device replacement can be performed by the healthcare professional or at home by the patient/caregiver. Do not attempt to replace the device until first discussing the procedure with your healthcare professional.

**CAUTION: Selection of the correct size device is critical for safety and comfort of the patient. An inappropriately sized device can cause necrosis, buried bumper syndrome and/or hypergranulation tissue. If patient size/weight has changed since the device was placed, fit is too tight or too loose, or it has been over six months since the device was last measured, it is recommended that you speak to your healthcare professional to ensure that device size does not need changed.**

## MEASURE THE STOMA LENGTH

**CAUTION:** Selection of the correct size device is critical for the safety and comfort of the patient. Measure the length of the patient's stoma with a stoma measuring device. The shaft length of the selected device should be the same as the length of the stoma. An inappropriately sized device can cause necrosis, buried bumper syndrome and/or hypergranulation tissue.

1. Please refer to Manufacturer's Directions for Use for the stoma measuring device being used.
2. Be sure to select the appropriate size MiniACE<sup>®</sup> Button for the abdominal wall thickness measured. If the measurement appears to be between two sizes, always select the next larger size MiniACE<sup>®</sup> Button. Once placed, the external flange should rotate easily.

**WARNING: UNDER-SIZING THE DEVICE MAY CAUSE EMBEDDING WITH EROSION INTO THE COLON WALL, TISSUE NECROSIS, INFECTION, SEPSIS AND ASSOCIATED SEQUELAE.**

## DEVICE PLACEMENT PROCEDURE

**CAUTION:** Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

**CAUTION:** Do not attempt to replace the device or check for placement verification until first discussing the procedure with your healthcare professional.

1. Select the correct MiniACE<sup>®</sup> Button size for placement.

**NOTE:** When replacing a device, stoma length should be periodically measured to ensure the correct MiniACE<sup>®</sup> Button size is being used. If patient size/weight has changed since device was placed, fit is too tight or loose, or it has been over six months since the device was last measured, it is recommended that you speak to your healthcare professional to ensure device length does not need changed.

2. Prior to placing the MiniACE<sup>®</sup> Button, inflate balloon through the balloon inflation port (see Fig 3) using a Luer slip syringe with distilled water or sterile water to the recommended fill volume. The recommended fill-volume can be found on Table 1 or printed above the balloon inflation port of the device. Remove the syringe and verify balloon integrity by gently squeezing the balloon to check for leaks. Visually inspect the balloon to verify symmetry. Verify that sizing information is appropriate for measured length. Reinsert the syringe and deflate all of the water from the balloon after inspection.
3. Lubricate the tip of the tube with water soluble lubricant. Do not use mineral oil or petroleum jelly. Insert optional introducer into the irrigation port if increased stiffness is desired during placement.
4. Gently guide the tube through the stoma and into the colon until the external flange is flush to the skin.
5. Remove Introducer (if used in step 3).
6. Inflate the balloon with sterile or distilled water according to the fill volume in the chart in Table 1.
7. Gently lift the tabs and check for signs of intestinal leakage.
8. To reduce risk of misconnection, place the provided tag on the button strap after device placement.

**NOTE:** If leakage is observed, increase balloon volume in increments of 0.5-1 ml. Do not exceed maximum fill volume.

**Table 1: Balloon Inflation Volumes**

| Fr Size | Minimum Fill Volume | Recommended Fill Volume | Maximum Fill Volume |
|---------|---------------------|-------------------------|---------------------|
| 10 F    | 1.5 ml              | 2 ml                    | 2.5 ml              |
| 12 F    | 2 ml                | 2.5 ml                  | 3 ml                |
| 14 F    | 3 ml                | 4 ml                    | 5 ml                |

## PLACEMENT VERIFICATION

1. If there is any trouble with placement, or any pain, blood, or site trauma occur during device removal or placement, consult with a physician before using the device to confirm proper device placement.
2. If no issues occur with placement, attach the Irrigation Set to the MiniACE<sup>®</sup> Button by lining up the dark line on the irrigation set connector with the dark line on the interlock of the MiniACE<sup>®</sup> Button. Press irrigation set connector to insert into the MiniACE<sup>®</sup> Button. Turn 3/4 to the right (clock-wise) to lock the irrigation set into place.
3. Attach a catheter tip syringe with 5ml of water into the Irrigation set and flush it through the device.
4. If no pain or discomfort occurs during the flush the placement is confirmed.

**WARNING: NEVER INJECT AIR INTO THE MINIACE<sup>®</sup> BUTTON.**

**WARNING: NEVER CONNECT THE IRRIGATION SET TO THE BALLOON INFLATION PORT.**

5. When placement is confirmed, colon irrigation can begin.

**NOTE:** In the case of an initial placement, consult with your physician regarding proper flushing and irrigation procedures immediately following placement.

## PLACEMENT CARE INSTRUCTIONS

1. The MiniACE® Button should be replaced periodically for optimal performance. If in place for three months or longer, frequent performance checks are recommended. Clogging and/or reduced flow are indications of diminished performance. Refer to the **TROUBLESHOOTING** section for additional signs of reduced performance or failure.
2. Be sure the balloon is within the colon and the MiniACE® Button rotates freely **before the irrigations begin**. Inject water or saline into the irrigation port to determine proper function of the valve.
3. The stoma area should be cleansed daily with a mild soap and water. The stoma site should be clean and dry at all times. The MiniACE® Button should be rotated daily for site hygiene.

## IRRIGATION SET INSTRUCTIONS FOR USE

**WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. ONLY USE THIS DEVICE TO CONNECT TO COMPATIBLE IRRIGATION DEVICES.**

1. Inspect contents for damage. If damaged, do not use. Obtain another package.
2. The irrigation set may be used for gravity irrigation/enema.
3. Ensure the clamp is closed and attach the irrigation set to the MiniACE® Button by lining up the dark line on the irrigation set connector with the dark line of the interlock on the MiniACE® Button. Fully press irrigation set connector to insert into the MiniACE® Button. Turn 3/4 to the right (clock-wise) to lock the irrigation set into place.
4. Attach the opposite end of the irrigation set to the connector being used. If using a bolus or luer connector, firmly insert the connector while using a slight rotation into the irrigation set to lock in place. If using a rotating connector, firmly rotate the connector clockwise into the irrigation set while avoiding excessive force to lock in place. Once connected, open the clamp to allow flow.

**WARNING: ONLY TIGHTEN BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.**

**WARNING: ENSURE DEVICE IS CONNECTED TO AN IRRIGATION PORT ONLY AND NOT TO AN IV SET.**

**WARNING: IF THE IRRIGATION SET IS NOT PROPERLY PLACED AND LOCKED, LEAKAGE MAY OCCUR. WHEN TWISTING THE CONNECTOR, HOLD THE MINIACE® BUTTON IN PLACE TO AVOID ROTATION WHILE PLACING THE IRRIGATION SET. NEVER INSERT IRRIGATION SET CONNECTOR INTO THE BALLOON INFLATION PORT. THIS MAY CAUSE THE INTERNAL BALLOON TO DEFLATE OR LEAD TO BALLOON FAILURE IF IMPROPER CONTENTS ARE INSERTED INTO THE BALLOON.**

5. When irrigation is complete, flush with 5-10ml water (if not using water for irrigation). Close the clamp to prevent leakage while removing the irrigation set. Remove the irrigation set by holding the MiniACE® Button in place, grasp the irrigation connector, twist 3/4 turn to the left (counter clock-wise) to unlock, and gently remove the irrigation set from the MiniACE® Button. Snap the MiniACE® Button plug in place to keep the lumen clean.
6. The irrigation set should be cleansed with soap and water and thoroughly rinsed. Store in a clean, dry place. Never use a dishwasher for cleaning.

**WARNING: IN CASE OF FEVER, ABDOMINAL DISTENTION, INFECTION, BLOCKAGE, OR TISSUE NECROSIS, PATIENTS SHOULD SEE THEIR PHYSICIAN IMMEDIATELY.**

**WARNING: DO NOT PLACE ANY FOREIGN OBJECTS INTO THE IRRIGATION PORT.**

**WARNING: IF USING A TRANSITION CONNECTOR, ENSURE THAT IT IS NOT CONNECTED TO AN IV SET PRIOR TO USE.**

**WARNING: WHEN USING A BOLUS STYLE IRRIGATION SET, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO BREATHING SYSTEM, LIMB CUFF, AND NEURAXIAL CONNECTORS.**

**NOTE:** The ENFi® connection is intended to reduce misconnections with other non-enteral healthcare applications.

## UNCLOGGING A DEVICE

First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a catheter tip syringe to an irrigation set and attach into the interlock connector. Fill the syringe with warm water and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced.

**CAUTION: Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.**

## REMOVAL OF THE MINIACE® BUTTON

1. Deflate the balloon by attaching the Luer slip tip syringe into the balloon inflation port (Fig. 3) to withdraw the sterile or distilled water.
2. Gently remove the device from the stoma site.
3. Replace as necessary following the MiniACE® Button directions for use.

**NOTE:** Spontaneous closure of the stoma may occur as early as 24 hours after removal. Insert new device if irrigation by this route is still intended.

**WARNING: DO NOT CUT OFF THE INTERNAL BOLSTER OR TUBING, OR ALLOW TO PASS THROUGH THE INTESTINAL TRACT.**

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

## DEVICE LONGEVITY

Low-profile balloon cecostomy devices are meant to be periodically replaced for optimal performance and functionality.

Exact device longevity cannot be predicted. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-6 months. Some factors that can lead to reduced longevity include: intestinal pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, incorrect stoma length measurement, and overall tube care.

For optimal performance, it is recommended that MiniACE® Button devices be changed every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure. If devices are failing or performance is degrading earlier than the typical range for device longevity, it is recommended that you speak with your healthcare professional regarding eliminating common factors that can lead to early device degradation. Also reference the **TROUBLESHOOTING** section for additional information regarding device performance issues.

**NOTE:** To help prevent un-needed hospital visits, it is recommended that a spare device is kept on hand at all times for replacement in case device failure occurs prior to scheduled replacement.

## MRI SAFETY INFORMATION

Non-clinical testing demonstrated that the MiniACE® Button is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 1,000-gauss/cm (10-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the MiniACE® Button is expected to produce a maximum temperature rise of 1.6°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the MiniACE® Button extends approximately 20-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

## TROUBLESHOOTING

Long-term device performance and functionality is dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your irrigation device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

**A tear has formed:** Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.

**Balloon inflation valve leakage:** Leakage from this valve typically occurs due to residue in the fill-valve. Always use a clean syringe when inflating the balloon and only inflate water through the valve. Leakage can also occur due to valve becoming stuck if syringe is inserted too hard into valve. Insert the syringe into the balloon inflation port if leakage is occurring to attempt to reset valve. Several attempts may be necessary before valve resets.

## TROUBLESHOOTING

**Leaking of balloon volume:** If balloon is deflating, remove balloon from stoma and inflate balloon with recommended fill volume. Check balloon for leakage by gently massaging tubing and balloon. If no leakage is noted, re-insert balloon back into stoma and re-inflate the balloon to desired fill-volume. Do not exceed maximum fill volume. Only access the balloon inflation port for balloon inflation and deflation. Other uses can result in balloon leakage or device failure. **NOTE:** The balloon is made of a semi-permeable material and can lose a small amount of fill-volume over time depending on environment and usage conditions.

**Anti-reflux valve leakage or blockage:** Leakage/blockage of the anti-reflux valve typically occurs due to residue (irrigation solution, intestinal contents, etc) becoming stuck in the valve area, preventing the valve from fully closing. Make sure device is flushed after each use. The valve can also become inverted in rare cases. Insert irrigation set into port to reset valve if this occurs.

**Tubing has reduced flow or has become clogged:** Tubing can become blocked due to not properly flushing after each use, use of thick irrigation solutions, and/or fungus growth. If clogged, refer to the **UNCLOGGING A DEVICE** section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced.

**Device fit is too tight or too loose:** Balloon fit can be adjusted by adjusting balloon inflation within the balloon inflation range in **Table 1**. If balloon fit is too loose, increase balloon fill-volume above recommended but not above max fill-volume. If balloon fit is too tight, decrease fill-volume below recommended but not below min fill-volume. If device does not properly fit with the fill-volume range, a new device length will be needed.

**Balloon is misshaped:** Be sure to inflate and inspect balloon prior to placement to check balloon symmetry. Balloons becoming misshaped typically occur due to excessive force or tension on the device (pulling device out of stoma while balloon is inflated). Balloons can be gently massaged with fingers back into symmetry if not excessively misshaped. Devices should be replaced if balloon is found excessively misshaped.

**Device has become discolored:** The device can become discolored over days to months of use. This is normal depending on the types of solutions being used with the device.

**Balloon will not inflate or deflate:** Inflation and deflation problems typically occur due to residue blocking the fill-lumen. Always use a clean syringe when inflating the balloon. In some rare occurrences, fungus can grow inside of the tubing and block the fill-lumen. Fungus growth can occur depending on patient environment and irrigation fluids being administered through the device. If balloon will not deflate, contact your healthcare professional for removal. If deflation problem is the result of fungal growth, eliminating source of fungal growth or anti-fungal medication may be required.

**Interlock failure or cracked:** The interlock has been designed to withstand extreme forces without detaching or cracking. However, the strength of the bond and material can reduce over pro-longed use depending on solutions used through the device. The device should be replaced if interlock is found cracked, leaking, or separating from the device.

**Foul smell coming from the device:** Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed and stoma site should be gently cleaned with soap and warm water. If foul smell does not go away, it is recommended that you contact your healthcare professional.

**Balloon failure:** Early balloon failure can occur due to a number of patient or environment factors, including but not limited to: Intestinal pH, diet, certain medications, improper balloon fill volume, placement of device, improperly administering irrigation solution into balloon port, trauma, contact with a sharp or abrasive material, incorrect stoma length measurement, and overall care of the device.

**Plug will not stay closed:** Ensure that plug is being firmly and fully pressed into the interlock connector. If plug is not staying closed, check the plug and irrigation-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.

## THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

Australian users please visit us at the link below for the MiniACE® Balloon Button Patient Information Leaflet:  
[www.appliedmedical.net/resources/miniace\\_aus\\_pil.pdf](http://www.appliedmedical.net/resources/miniace_aus_pil.pdf)

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Applied Medical Technology, Inc.