

AMT PRODUCT CATALOG

A Complete List of Products

For a complete list of products visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



AppliedMedical.net



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
MR Conditional | MR SAFE
Made In USA



© 2021 Applied Medical Technology, Inc. **C-4338-E**

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

Tyvek® is a registered trademark of E.I. du Pont de Nemours and Company or its affiliates | Salem Sump® is a registered trademark of Cardinal Health Ireland Unlimited Company. | Chait Trapdoor™ is a trademark of Cook® Medical | ENFit® is a registered trademark of GEDSA. | Mic-key® is a registered trademark of Avent, Inc.

TABLE OF CONTENTS:

	PAGE
NutraGlide™ Nasal Feeding Tube	02 - 03
AMT Bridle™ Nasal Tube Retaining System Family	04 - 05
Bridle Pro®	
AMT Bridle™	
MiniONE® Gastrostomy Tube Family	06 - 11
MiniONE® Balloon Button	
MiniONE® Capsule Non-Balloon Button	
MiniONE® Non-Balloon Button	
Traditional Length Gastrostomy Tube Family	12 - 17
Balloon G-Tube	
Capsule Dome G-Tube	
Capsule Monarch® G-Tube	
G-JET® Gastric-Jejunal Feeding Tube Family	18 - 27
micro G-JET® Button	
G-JET® Button	
Traditional G-JET®	
MiniACE® Enema Button	28 - 29
Accessories	30 - 37
Extension Sets: Legacy + ENFit®	
CINCH® Tube Securement Device	
AMT Clamp™ Tube Securement Device	
AMT Syringes	
Guidewires	
Transition Adapters: ENFit® + Legacy	
Replacement Y-Port Adapters	
Surgical Line	38 - 43
Initial Placement Gastrostomy Kits	
TLC® Retractors + Stays	
Specialized Surgical Kits	
Remora® Smoke + Fluid Evacuation Device	
Dedication to Education	44 - 45
Contact Us and References	46 - 48

MEDICAL DEVICE SOLUTIONS:

Applied Medical Technology, Inc. (AMT) is a global leader in the development, design, and manufacture of enteral feeding devices, accessories, and surgical products. We are honored to serve the global community as we continue to introduce innovative products and materials, designed to increase patient comfort while improving patient outcomes. In our state-of-the-art facility near Cleveland, Ohio, we believe in manufacturing the highest quality medical devices to meet the needs of our customers.

AMT has been a pioneering force within the medical device market, revolutionizing the standard of care for both pediatric and adult patients. Our Founder/CEO, Dr. Picha, co-invented and manufactured the first low profile gastrostomy button, changing the enteral nutrition community forever.

Over time, AMT's product portfolio has evolved to include nasal tube retention systems, innovative solutions for pediatric jejunal feeding, safe tube securement options, bowel management, and a complete line of ENFit® products. AMT's patient focus has even led to the development of an enteral nutrition mobile app: AMT ONE Source®.

At AMT, we stand by our motto: Innovating. Educating. Changing Lives.™ We are committed to improving the quality of life for patients, caregivers and clinicians around the world, by providing educational resources, support, and continuous device innovation.

ALL AMT PRODUCTS ARE PROUDLY MADE IN THE USA!



Dr. George Picha

M.D., Ph.D., F.A.C.S.
 Founder & CEO of Applied Medical Technology, Inc.



Innovating. Educating. Changing Lives.™



NutraGlide™ + Bridle Pro®
Essential Nutrition, Uninterrupted

Improve Nutrition

- Bridled patients receive 16% more of their caloric goals.¹
- Consistent nutrition delivery can reduce length of stay.²

Safe + Secure

The NutraGlide™ Nasal Feeding Tube is compatible with the AMT Bridle™ Family of Nasal Tube Retaining Systems.

AMT Offers the Largest Range of Pediatric + Adult Bridle Sizes³


Patient Comfort

The dissolvable tip of the NutraGlide™ and the small probes of the AMT Bridle™ Family allow for a more comfortable experience overall.

Decrease Costs

Using the NutraGlide™ in combination with the AMT Bridle Family can DECREASE costs by reducing the risk of clog formation and tube pullouts.

NutraGlide™ Nasal Feeding Tube




Reduced Clogging

The NutraGlide™ features a unique expanding internal diameter which is designed to reduce the risk of clogging.

Dissolvable Tip (8F, 10F)

The NutraGlide™ features a dissolvable tip created with OptiSol™ technology. The hydrophilic nature of the tip aids in placement as it self-lubricates once moisture is introduced. Our rounded, dissolvable tip offers comfort during placement and the superior flow of an open-ended tube during feeding.



Changing the Standard of Care for Nasal Feeding

NutraGlide™ | Order Information

FRENCH SIZE	LENGTH (cm)	LEGACY ORDER #		ENFIT® ORDER #	
		STYLET	NO STYLET	STYLET	NO STYLET
5F	38 cm	NFT-05038	NFT-05038-N	NFT-05038-I	NFT-05038-N-I
	55 cm	NFT-05055	NFT-05055-N	NFT-05055-I	NFT-05055-N-I
8F	110 cm	NFT-08110	NFT-08110-N	NFT-08110-I	NFT-08110-N-I
	140 cm	NFT-08140	NFT-08140-N	NFT-08140-I	NFT-08140-N-I
10F	110 cm	NFT-10110	NFT-10110-N	NFT-10110-I	NFT-10110-N-I
	140 cm	NFT-10140	NFT-10140-N	NFT-10140-I	NFT-10140-N-I

Add the prefix "S" to order Sterile Nasal Feeding Tubes, e.g. SNFT-05038. NutraGlide™ is sold 10/box.

Bridle Pro®

AMT Bridle™ Nasal Tube Retaining System Family



An AMT Exclusive. The New Standard of Care.



1 Flexible Tubing

Easy to clean and comfortable for patients.

2 Alignment Markings

Assures device alignment for proper placement.

3 Pre-Attached Clip

Improved patient safety and improved handling.

4 Ease of Placement

Magnets at the distal end help facilitate a quick and easy connection.

Bridle Family Clips

Clips accommodate several nasal tube sizes.



AMT Bridle™ + Bridle Pro® | Order Information

AMT BRIDLE PRO® PEDIATRIC SIZE NASAL TUBE RETAINING SYSTEM		
DESCRIPTION	ORDER #	BOX/QTY
5F-6F Pediatric Bridle Pro® Range (Purple)	4-420506P	5
8F-10F Pediatric Bridle Pro® Range (Teal)	4-420810P	5
8F Pediatric Bridle Pro® (White)	4-4208P	5
Kit contents include: flexible bridle catheter, stylet guide, retrieval probe, pre-attached clip, clip opening tool and water soluble lubricant.		
AMT BRIDLE PRO® NASAL TUBE RETAINING SYSTEM		
8F-10F Bridle Pro® Range (Teal)	4-420810	5
12F-14F Bridle Pro® Range (Dark Blue)	4-421214	5
16F-18F Bridle Pro® Range (Yellow)	4-421618	5
8F Bridle Pro® (White)	4-4208	5
10F Bridle Pro® (Teal)	4-4210	5
12F Bridle Pro® (Dark Blue)	4-4212	5
14F Bridle Pro® (Yellow)	4-4214	5
16F Bridle Pro® (Yellow)	4-4216	5
18F Bridle Pro® (Light Blue)	4-4218	5
Kit contents include: flexible bridle catheter, stylet guide, retrieval probe, pre-attached clip, clip opening tool and water soluble lubricant.		
AMT BRIDLE PRO® SPARE CLIPS		
5F-6F Bridle Pro® Range (Purple)	4-460506	5
8F-10F Bridle Pro® Range (Teal)	4-460810	5
12F-14F Bridle Pro® Range (Dark Blue)	4-461214	5
16F-18F Bridle Pro® Range (Yellow)	4-461618	5
8F Bridle Pro® (White)	4-4608	5
10F Bridle Pro® (Teal)	4-4610	5
12F Bridle Pro® (Dark Blue)	4-4612	5
14F Bridle Pro® (Yellow)	4-4614	5
16F Bridle Pro® (Yellow)	4-4616	5
18F Bridle Pro® (Light Blue)	4-4618	5
AMT BRIDLE™ NASAL TUBE RETAINING SYSTEM		
8F AMT Bridle™ (White)	4-4108	5
10F AMT Bridle™ (Teal)	4-4110	5
12F AMT Bridle™ (Dark Blue)	4-4112	5
14F AMT Bridle™ (Yellow)	4-4114	5
16F AMT Bridle™ (Yellow)	4-4116	5
18F AMT Bridle™ (Light Blue)	4-4118	5
AMT BRIDLE™ SPARE CLIPS		
8F AMT Bridle™ (White)	4-4508	5
10F AMT Bridle™ (Teal)	4-4510	5
12F AMT Bridle™ (Dark Blue)	4-4512	5
14F AMT Bridle™ (Yellow)	4-4514	5
16F AMT Bridle™ (Yellow)	4-4516	5
18F AMT Bridle™ (Light Blue)	4-4518	5

MiniONE® Gastrostomy Tube Family



INNOVATIVE

MiniONE® Balloon Button

Low Profile Gastrostomy Tube

Legacy + ENFit® | 12F-24F

Page 08-09



AppliedMedical.net




EXCLUSIVE

MiniONE® Capsule Non-Balloon

Low Profile Gastrostomy Tube

Legacy + ENFit® | 14F, 18F, 20F, 24F

Page 10-11



AppliedMedical.net



DEPENDABLE

MiniONE® Non-Balloon

Low Profile Gastrostomy Tube

Legacy + ENFit® | 14F, 18F, 20F, 24F

Page 10-11



AppliedMedical.net

Customer Success Stories

AMT's Commitment: Real Customers. Real Experiences.

They make it possible for Sam to do more than ever before!



Sam's Mom
AMT Customer



MiniONE® Family


Connect with Us Online!



MiniONE® Balloon Button

Low Profile Balloon Button Gastrostomy Tube

STUDY CONFIRMS: Over 70% of families prefer the AMT MiniONE® over the Avanos Mic-Key®.⁴



1 Soft Silicone
Flexible external bolster.

2 Glow Green™ Technology
Easier nighttime feeds.

3 Balloon Vol. Clearly Labeled
The Balloon fill volume printed directly on balloon fill port.

4 Exclusive AMT Balloon
Apple-shaped, expands horizontally to create a superior gastric seal.⁷

MiniONE® Balloon Button | Order Information

LEGACY ORDER #						
LENGTH (cm)	12F	14F	16F	18F	20F	24F
0.8	M1-5-1208	M1-5-1408				
1.0	M1-5-1210	M1-5-1410	M1-5-1610	M1-5-1810	M1-5-2010	
1.2	M1-5-1212	M1-5-1412	M1-5-1612	M1-5-1812	M1-5-2012	
1.5	M1-5-1215	M1-5-1415	M1-5-1615	M1-5-1815	M1-5-2015	M1-5-2415
1.7	M1-5-1217	M1-5-1417	M1-5-1617	M1-5-1817	M1-5-2017	M1-5-2417
2.0	M1-5-1220	M1-5-1420	M1-5-1620	M1-5-1820	M1-5-2020	M1-5-2420
2.3	M1-5-1223	M1-5-1423	M1-5-1623	M1-5-1823	M1-5-2023	M1-5-2423
2.5	M1-5-1225	M1-5-1425	M1-5-1625	M1-5-1825	M1-5-2025	M1-5-2425
2.7	M1-5-1227	M1-5-1427	M1-5-1627	M1-5-1827	M1-5-2027	M1-5-2427
3.0	M1-5-1230	M1-5-1430	M1-5-1630	M1-5-1830	M1-5-2030	M1-5-2430
3.5	M1-5-1235	M1-5-1435	M1-5-1635	M1-5-1835	M1-5-2035	M1-5-2435
4.0	M1-5-1240	M1-5-1440	M1-5-1640	M1-5-1840	M1-5-2040	M1-5-2440
4.4		M1-5-1444	M1-5-1644	M1-5-1844	M1-5-2044	M1-5-2444
5.0		M1-5-1450	M1-5-1650	M1-5-1850	M1-5-2050	M1-5-2450
5.5		M1-5-1455	M1-5-1655	M1-5-1855	M1-5-2055	M1-5-2455
6.0		M1-5-1460	M1-5-1660	M1-5-1860	M1-5-2060	M1-5-2460
6.5		M1-5-1465				
7.0		M1-5-1470				
8.0		M1-5-1480				
9.0		M1-5-1490				
10.0		M1-5-14100				

ENFit® ORDER #						
LENGTH (cm)	12F	14F	16F	18F	20F	24F
0.8	M1-5-1208-I	M1-5-1408-I				
1.0	M1-5-1210-I	M1-5-1410-I	M1-5-1610-I	M1-5-1810-I	M1-5-2010-I	
1.2	M1-5-1212-I	M1-5-1412-I	M1-5-1612-I	M1-5-1812-I	M1-5-2012-I	
1.5	M1-5-1215-I	M1-5-1415-I	M1-5-1615-I	M1-5-1815-I	M1-5-2015-I	M1-5-2415-I
1.7	M1-5-1217-I	M1-5-1417-I	M1-5-1617-I	M1-5-1817-I	M1-5-2017-I	M1-5-2417-I
2.0	M1-5-1220-I	M1-5-1420-I	M1-5-1620-I	M1-5-1820-I	M1-5-2020-I	M1-5-2420-I
2.3	M1-5-1223-I	M1-5-1423-I	M1-5-1623-I	M1-5-1823-I	M1-5-2023-I	M1-5-2423-I
2.5	M1-5-1225-I	M1-5-1425-I	M1-5-1625-I	M1-5-1825-I	M1-5-2025-I	M1-5-2425-I
2.7	M1-5-1227-I	M1-5-1427-I	M1-5-1627-I	M1-5-1827-I	M1-5-2027-I	M1-5-2427-I
3.0	M1-5-1230-I	M1-5-1430-I	M1-5-1630-I	M1-5-1830-I	M1-5-2030-I	M1-5-2430-I
3.5	M1-5-1235-I	M1-5-1435-I	M1-5-1635-I	M1-5-1835-I	M1-5-2035-I	M1-5-2435-I
4.0	M1-5-1240-I	M1-5-1440-I	M1-5-1640-I	M1-5-1840-I	M1-5-2040-I	M1-5-2440-I
4.4		M1-5-1444-I	M1-5-1644-I	M1-5-1844-I	M1-5-2044-I	M1-5-2444-I
5.0		M1-5-1450-I	M1-5-1650-I	M1-5-1850-I	M1-5-2050-I	M1-5-2450-I
5.5		M1-5-1455-I	M1-5-1655-I	M1-5-1855-I	M1-5-2055-I	M1-5-2455-I
6.0		M1-5-1460-I	M1-5-1660-I	M1-5-1860-I	M1-5-2060-I	M1-5-2460-I
6.5		M1-5-1465-I				
7.0		M1-5-1470-I				
8.0		M1-5-1480-I				
9.0		M1-5-1490-I				
10.0		M1-5-14100-I				

MiniONE® Family

Feed Sets Included in Kit | Full Extension Set Portfolio on Page 30-31



MiniONE® Button comes packaged as 1/box (sterile). Kit contents include: device (1), bolus extension set (1), y-port extension set (1), 5cc/10cc balloon syringe (1), 35cc bolus syringe (1), optional introducer/stiffener (12F-18F only), gauze, and water soluble lubricant (1).

MiniONE® Non-Balloon + Capsule Non-Balloon Button

Low Profile Non-Balloon Gastrostomy Tube



Exclusive Capsule Technology

Material facilitates easy and comfortable insertion.

Glow Green™ Technology

Easier nighttime feeds.

Improve Comfort

Combines the advantage of a low profile device with the longer life of a non-balloon device.

Exclusive AMT Capsule Technology



More Durable

Internal bolster is smaller than most devices, durable and less likely to block the pylorus.¹²


MiniONE® Non-Balloon Button | Order Information

LENGTH (cm)	14F		18F		20F		24F	
	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #
1.0	M1-1-1410	M1-1-1410-I						
1.2	M1-1-1412	M1-1-1412-I	M1-1-1812	M1-1-1812-I				
1.5	M1-1-1415	M1-1-1415-I	M1-1-1815	M1-1-1815-I				
1.7	M1-1-1417	M1-1-1417-I	M1-1-1817	M1-1-1817-I	M1-1-2017	M1-1-2017-I	M1-1-2417	M1-1-2417-I
2.0	M1-1-1420	M1-1-1420-I	M1-1-1820	M1-1-1820-I				
2.4			M1-1-1824	M1-1-1824-I	M1-1-2024	M1-1-2024-I	M1-1-2424	M1-1-2424-I
2.5	M1-1-1425	M1-1-1425-I						
3.0	M1-1-1430	M1-1-1430-I	M1-1-1830	M1-1-1830-I				
3.4			M1-1-1834	M1-1-1834-I	M1-1-2034	M1-1-2034-I	M1-1-2434	M1-1-2434-I
4.4			M1-1-1844	M1-1-1844-I	M1-1-2044	M1-1-2044-I	M1-1-2444	M1-1-2444-I

MiniONE® Capsule Non-Balloon Button | Order Information

LENGTH (cm)	14F		18F		20F		24F	
	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #
1.0	M1-2-1410	M1-2-1410-I						
1.2	M1-2-1412	M1-2-1412-I	M1-2-1812	M1-2-1812-I				
1.5	M1-2-1415	M1-2-1415-I	M1-2-1815	M1-2-1815-I				
1.7	M1-2-1417	M1-2-1417-I	M1-2-1817	M1-2-1817-I	M1-2-2017	M1-2-2017-I	M1-2-2417	M1-2-2417-I
2.0	M1-2-1420	M1-2-1420-I	M1-2-1820	M1-2-1820-I				
2.4			M1-2-1824	M1-2-1824-I	M1-2-2024	M1-2-2024-I	M1-2-2424	M1-2-2424-I
2.5	M1-2-1425	M1-2-1425-I						
3.0	M1-2-1430	M1-2-1430-I	M1-2-1830	M1-2-1830-I				
3.4			M1-2-1834	M1-2-1834-I	M1-2-2034	M1-2-2034-I	M1-2-2434	M1-2-2434-I
4.4			M1-2-1844	M1-2-1844-I	M1-2-2044	M1-2-2044-I	M1-2-2444	M1-2-2444-I

Feed Sets Included in Kit | Full Extension Set Portfolio on Page 30-31



MiniONE® Button comes packaged as 1/box (sterile). Kit contents include: device (1), bolus extension set (1), y-port extension set (1), 10cc syringe (1), 35cc bolus syringe (1), Placement T-Handle/ Removal T-Handle*, snap arm assembly, and gauze.

Spare Removal Tools | For Obturated Non-Balloon and Capsule Non-Balloon Buttons

FRENCH SIZE	DESCRIPTION	ORDER #	BOX/QTY
14F	1.0, 1.2 cm	R-S1012	1
14F	1.5, 1.7 cm	R-S1517	1
14F	2.0 cm	R-S2023	1
14F	2.5 cm	R-S2427	1
14F	3.0 cm	R-S3033	1
18F	1.2 cm	R-L1012	1
18F, 20F, 24F	1.5, 1.7 cm	R-L1517	1
18F	2.0 cm	R-L2023	1
18F	2.4 cm	R-L2427	1
18F	3.0 cm	R-L3033	1
18F, 20F, 24F	3.4 cm	R-L3437	1
18F, 20F, 24F	4.4 cm	R-L4044	1

AMT Obturator for Distending button during insertion and removal.



*Non-Balloon Kit: Placement and Removal T-Handle
Capsule Non-Balloon Kit: Removal T-Handle

Traditional Length Gastrostomy Tube Family

Offering Long Term Performance



COMFORT + SAFETY

Balloon G-Tube

Legacy and ENFit® | 12F – 24F

Page 14-15



AppliedMedical.net



EASE OF PLACEMENT

Capsule Dome G-Tube

Legacy and ENFit® | 20F – 24F

Page 16-17



AppliedMedical.net



EXCLUSIVE

Capsule Monarch® G-Tube

Legacy and ENFit® | 12F, 14F, 18F, 20F

Page 16-17



AppliedMedical.net

Customer Success Stories

AMT's Commitment: Real Customers. Real Experiences.

I am fully convinced that he is here today solely because of his feeding tube. It saved his life.



Eli's Mom
AMT Customer



Traditional Family

Connect with Us Online!




Traditional G-Tube, Balloon

Adjustable Length Gastrostomy Tube



Adjustable External Bolster

Can accommodate various stoma lengths and has stabilizer feet that allow for air flow.



Traditional G-Tube, Balloon

Adjustable Length Gastrostomy Tube

Reduce Leaks

- Exclusive AMT Balloon expands horizontally to provide superior stoma coverage.⁷

Improve Comfort + Safety

- Balloon protects the distal tip of the tube, lowering the risk for gastric ulceration.
- Enhanced External Retention Ring minimizes contact and allows skin to breathe.

Ease of Placement

- Often placed in the office or on-site.
- Tapered tip design for easier placement.
- Guidewire compatible.
- Radiopaque for x-ray confirmation.

G-Tube, Balloon | Order Information

DESCRIPTION	LEGACY ORDER #	ENFit® ORDER #
12F Balloon G-Tube	7-1210	7-1210-I
14F Balloon G-Tube	7-1410	7-1410-I
16F Balloon G-Tube	7-1610	7-1610-I
18F Balloon G-Tube	7-1810	7-1810-I
20F Balloon G-Tube	7-2010	7-2010-I
24F Balloon G-Tube	7-2410	7-2410-I

AMT Balloon G-Tube is packaged 1/Box.
All ENFit® Balloon G-Tube kits are pouched with an ENFit® Single-Port to Y-Port Adapter and Pinch Clamp.



REPLACEMENT Y-PORT ADAPTERS ORDER INFORMATION - SEE PAGE 37		
DESCRIPTION	ORDER #	BOX/QTY
ENFit® Single-Port to Y-Port Adapter	4-0005-ISOSAF	10

Traditional G-Tube, Non-Balloon

Adjustable Length Gastrostomy Tubes


Capsule Dome G-Tube



Capsule Monarch® G-Tube



Stages of Capsule Bolster Deployment



- 1: Encapsulated Bolster Inserted into Established Stoma
- 2: Tether Pulled to Split Capsule Material, Stiffener Removed
- 3: Capsule Material Falls Away and Dissolves Rapidly

Traditional G-Tube, Non-Balloon

Adjustable Length Gastrostomy Tube

Exclusive AMT Capsule Technology
Ease of Placement. Improve Comfort + Durability. Reduce Leaks.

G-Tube, Capsule Dome Non-Balloon | Order Information

DESCRIPTION	LEGACY ORDER #	ENFit® ORDER #
20F	9-C2020	9-C2020-I
24F	9-C2420	9-C2420-I
20F w/Guidewire (0.038" x 90 cm)	9-CG2020	9-CG2020-I
24F w/Guidewire (0.038" x 90 cm)	9-CG2420	9-CG2420-I


AMT Capsule Dome G-Tube is Packaged 1/box.
All ENFit® Capsule Dome kits are pouched with an ENFit® Single-Port to Y-Port Adapter and Pinch Clamp.
For Guidewire Kits: Guidewire is included in a separate pouch, sterile, 90 cm length.

G-Tube, Capsule Monarch® Non-Balloon | Order Information

DESCRIPTION	LEGACY ORDER #	ENFit® ORDER #
12F	9-C1210	9-C1210-I
14F	9-C1410	9-C1410-I
18F	9-C1810	9-C1810-I
20F	9-C2010	9-C2010-I
12F w/Guidewire (0.018" x 70 cm)	9-CG1210	9-CG1210-I
14F w/Guidewire (0.018" x 70 cm)	9-CG1410	9-CG1410-I
18F w/Guidewire (0.038" x 90 cm)	9-CG1810	9-CG1810-I
20F w/Guidewire (0.038" x 90 cm)	9-CG2010	9-CG2010-I

AMT Capsule Monarch® G-Tube is Packaged 1/box.
All ENFit® Capsule Monarch® kits are pouched with an ENFit® Single-Port to Y-Port Adapter and Pinch Clamp.
For Guidewire Kits: Guidewire is included in a separate pouch, sterile, length varies by French size.

Replacement Y-Port Adapters | Order Information - See Page 37



LEGACY
Purple Y-Port Adapter

LEGACY
Clear Y-Port XL Adapter

ENFit®
Single-Port Adapter

ENFit®
Single-Port to Y-Port Adapter

G-JET® Gastric-Jejunal Feeding Tube Family



PEDIATRIC

micro G-JET®

Low Profile Gastric-Jejunal Button
Legacy + ENFit® | 14F Gastric – 8F Jejunal
Page 20-21



AppliedMedical.net



LOW PROFILE

G-JET® Button

Low Profile Gastric-Jejunal Button
Legacy + ENFit® | 14F, 16F, 18F
Page 22-24



AppliedMedical.net



TRADITIONAL

Traditional G-JET®

Adjustable Length Gastric-Jejunal Enteral Feeding Tube
Legacy + ENFit® | 14F, 16F, 18F, 22F
Page 26-27



AppliedMedical.net

Customer Success Stories

AMT's Commitment: Real Customers. Real Experiences.

In fact, we've had virtually zero issues since switching to the micro G-JET®. We are incredibly thankful to have found this feeding tube when we did. It has been the perfect solution!



Mimi's Dad,
micro G-JET® User



G-JET® Family

Connect with Us Online!



micro G-JET® Button

Pediatric G-J Enteral Feeding Tube

One of a Kind!

Tubing transitions from a 14F Gastric segment to an 8F Jejunal segment.



1 Low Profile Bolster

Clearly labeled G + J Ports, Featuring Glow Green™ Technology.

2 Exclusive AMT Balloon

The balloon you trust, adapted for smaller stomachs.

3 8F + Anti-Kink Protection

Anti-Kink coil spans the entire 8F jejunal segment.

4 Soft Tip


Gentle, atraumatic.

micro G-JET® Button

Pediatric G-J Enteral Feeding Tube

Designed for Pediatric Patients⁵

The Only Low Profile Gastric-Jejunal Feeding Device with an 8F Jejunal Segment




micro G-JET® Button | Order Information

GASTRIC 14F	JEJUNAL 8F	LEGACY ORDER #	ENFit® ORDER #
0.8 cm	10 cm	MGJ-1408-10	MGJ-1408-10-I
	15 cm	MGJ-1408-15	MGJ-1408-15-I
	22 cm	MGJ-1408-22	MGJ-1408-22-I
1.0 cm	10 cm	MGJ-1410-10	MGJ-1410-10-I
	15 cm	MGJ-1410-15	MGJ-1410-15-I
	22 cm	MGJ-1410-22	MGJ-1410-22-I
1.2 cm	10 cm	MGJ-1412-10	MGJ-1412-10-I
	15 cm	MGJ-1412-15	MGJ-1412-15-I
	22 cm	MGJ-1412-22	MGJ-1412-22-I
	30 cm	MGJ-1412-30	MGJ-1412-30-I
1.5 cm	10 cm	MGJ-1415-10	MGJ-1415-10-I
	15 cm	MGJ-1415-15	MGJ-1415-15-I
	22 cm	MGJ-1415-22	MGJ-1415-22-I
	30 cm	MGJ-1415-30	MGJ-1415-30-I
1.7 cm	10 cm	MGJ-1417-10	MGJ-1417-10-I
	15 cm	MGJ-1417-15	MGJ-1417-15-I
	22 cm	MGJ-1417-22	MGJ-1417-22-I
	30 cm	MGJ-1417-30	MGJ-1417-30-I
2.0 cm	22 cm	MGJ-1420-22	MGJ-1420-22-I
	30 cm	MGJ-1420-30	MGJ-1420-30-I
2.3 cm	22 cm	MGJ-1423-22	MGJ-1423-22-I
	30 cm	MGJ-1423-30	MGJ-1423-30-I
2.5 cm	22 cm	MGJ-1425-22	MGJ-1425-22-I
	30 cm	MGJ-1425-30	MGJ-1425-30-I

micro G-JET® Button comes packaged as 1/box (sterile). | Kit contents include: device (1), gastric extension set (1), jejunal extension set (1), guidewire adapter (1), 5cc/10cc balloon syringe (1), 10cc contrast syringe (1), 35cc bolus syringe (1), winged luer adapter (1), and gauze.

Extension + Feed Sets Included in Kit | Full Extension Set Portfolio on Pages 30-33

*AMT Jejunal feed sets are compatible with Mic-Key® devices.



G-JET® Family

G-JET® Button

Low Profile Gastric-Jejunal Enteral Feeding Tube

An Advantage in Enteral Nutrition
Anti-Kink. Exclusive AMT Balloon. Glow Green™ Technology.



1 Guidewire Adapter

For ease of placement.

2 Glow Green™ Technology

Easier nighttime feeds.

3 Exclusive AMT Balloon

Apple-shaped, expands horizontally creating a superior gastric seal.⁷

4 Anti-Kink Protection

16F and 18F sizes offer spring reinforcement to maintain flow.


5 Single Lumen

Internal geometry designed to reduce clogging.

G-JET® Button

Low Profile Gastric-Jejunal Enteral Feeding Tube

Same French Size, Larger Lumens
Designed to Reduce Clogging



GASTRIC Lumen Comparison⁶

AMT 14F is **58% Larger** than Competitor 14F

AMT 16F is **41% Larger** than Competitor 16F

AMT 18F is **34% Larger** than Competitor 18F

JEJUNAL Lumen Comparison⁶


AMT 14F is **56% Larger** than Competitor 14F
56% Larger than Competitor 16F

AMT 16F is **131% Larger** than Competitor 16F
81% Larger than Competitor 18F

AMT 18F is **172% Larger** than Competitor 18F
57% Larger than Competitor 22F

Mutually Exclusive Ports Leading the Way in Patient Safety

Accidentally feeding the stomach instead of the jejunum can be dangerous! Mutually exclusive **GASTRIC** and **JEJUNAL** ports are designed to reduce the risk of misconnections. Different feed/extension sets are required for each port, so it is imperative to maintain adequate stock of both sets. Always remember to **Feed the Green!**



Complete **GASTRIC** set offering found on pages 32-33 and complete **JEJUNAL** set offering found on pages 30-31.

G-JET® Button

Low Profile Gastric-Jejunal Enteral Feeding Tube


G-JET® Button | Order Information

GASTRIC	JEJUNAL	14F		16F		18F	
		LEGACY ORDER #	ENFit® ORDER #	LEGACY ORDER #	ENFit® ORDER #	LEGACY ORDER #	ENFit® ORDER #
1.0	15	GJ-1410-15	GJ-1410-15-I	GJ-1610-15	GJ-1610-15-I		
	22	GJ-1410-22	GJ-1410-22-I	GJ-1610-22	GJ-1610-22-I		
1.2	15	GJ-1412-15	GJ-1412-15-I	GJ-1612-15	GJ-1612-15-I		
	22	GJ-1412-22	GJ-1412-22-I	GJ-1612-22	GJ-1612-22-I	GJ-1812-22	GJ-1812-22-I
	30	GJ-1412-30	GJ-1412-30-I	GJ-1612-30	GJ-1612-30-I	GJ-1812-30	GJ-1812-30-I
1.5	15	GJ-1415-15	GJ-1415-15-I	GJ-1615-15	GJ-1615-15-I		
	22	GJ-1415-22	GJ-1415-22-I	GJ-1615-22	GJ-1615-22-I	GJ-1815-22	GJ-1815-22-I
	30	GJ-1415-30	GJ-1415-30-I	GJ-1615-30	GJ-1615-30-I	GJ-1815-30	GJ-1815-30-I
	45	GJ-1415-45	GJ-1415-45-I	GJ-1615-45	GJ-1615-45-I	GJ-1815-45	GJ-1815-45-I
1.7	15	GJ-1417-15	GJ-1417-15-I	GJ-1617-15	GJ-1617-15-I		
	22	GJ-1417-22	GJ-1417-22-I	GJ-1617-22	GJ-1617-22-I	GJ-1817-22	GJ-1817-22-I
	30	GJ-1417-30	GJ-1417-30-I	GJ-1617-30	GJ-1617-30-I	GJ-1817-30	GJ-1817-30-I
	45	GJ-1417-45	GJ-1417-45-I	GJ-1617-45	GJ-1617-45-I	GJ-1817-45	GJ-1817-45-I
2.0	22	GJ-1420-22	GJ-1420-22-I	GJ-1620-22	GJ-1620-22-I	GJ-1820-22	GJ-1820-22-I
	30	GJ-1420-30	GJ-1420-30-I	GJ-1620-30	GJ-1620-30-I	GJ-1820-30	GJ-1820-30-I
	45	GJ-1420-45	GJ-1420-45-I	GJ-1620-45	GJ-1620-45-I	GJ-1820-45	GJ-1820-45-I
2.3	22	GJ-1423-22	GJ-1423-22-I				
	30	GJ-1423-30	GJ-1423-30-I	GJ-1623-30	GJ-1623-30-I	GJ-1823-30	GJ-1823-30-I
	45	GJ-1423-45	GJ-1423-45-I	GJ-1623-45	GJ-1623-45-I	GJ-1823-45	GJ-1823-45-I
2.5	22	GJ-1425-22	GJ-1425-22-I				
	30	GJ-1425-30	GJ-1425-30-I	GJ-1625-30	GJ-1625-30-I	GJ-1825-30	GJ-1825-30-I
	45	GJ-1425-45	GJ-1425-45-I	GJ-1625-45	GJ-1625-45-I	GJ-1825-45	GJ-1825-45-I
2.7	30	GJ-1427-30	GJ-1427-30-I	GJ-1627-30	GJ-1627-30-I	GJ-1827-30	GJ-1827-30-I
	45	GJ-1427-45	GJ-1427-45-I	GJ-1627-45	GJ-1627-45-I	GJ-1827-45	GJ-1827-45-I
3.0	45	GJ-1430-45	GJ-1430-45-I	GJ-1630-45	GJ-1630-45-I	GJ-1830-45	GJ-1830-45-I
3.5	45	GJ-1435-45	GJ-1435-45-I	GJ-1635-45	GJ-1635-45-I	GJ-1835-45	GJ-1835-45-I
4.0	45	GJ-1440-45	GJ-1440-45-I	GJ-1640-45	GJ-1640-45-I	GJ-1840-45	GJ-1840-45-I
4.4	45	GJ-1444-45	GJ-1444-45-I	GJ-1644-45	GJ-1644-45-I	GJ-1844-45	GJ-1844-45-I
5.0	45	GJ-1450-45	GJ-1450-45-I	GJ-1650-45	GJ-1650-45-I	GJ-1850-45	GJ-1850-45-I
5.5	45	GJ-1455-45	GJ-1455-45-I	GJ-1655-45	GJ-1655-45-I	GJ-1855-45	GJ-1855-45-I
6.0	45	GJ-1460-45	GJ-1460-45-I	GJ-1660-45	GJ-1660-45-I	GJ-1860-45	GJ-1860-45-I
6.5	45	GJ-1465-45	GJ-1465-45-I	GJ-1665-45	GJ-1665-45-I	GJ-1865-45	GJ-1865-45-I
7.0	45	GJ-1470-45	GJ-1470-45-I	GJ-1670-45	GJ-1670-45-I	GJ-1870-45	GJ-1870-45-I

G-JET® Button comes packaged as 1/box (sterile). | Kit contents include: device (1), gastric extension set (1), jejunal extension set (1), guidewire adapter (1), 5cc/10cc balloon syringe (1), 10cc contrast syringe (1), 35cc bolus syringe (1), winged luer adapter (1), and gauze.

Extension + Feed Sets Included in Kit | Full Extension Set Portfolio on Page 30-33

*AMT Jejunal feed sets are compatible with Mic-Key® devices.



Customer Success Stories

AMT's Commitment: Real Customers. Real Experiences.

The AMT G-JET® would be the best fit for me and my lifestyle. It was flexible and so comfortable I could forget I even had a feeding tube! The anti-kink technology was important to me as well because I didn't want my tube to coil back up in my stomach. Bonus points for glowing in the dark too!

Thanks to my G-JET®, I am a 3rd degree black belt, a 9-time Martial Arts Taekwondo World Champion, trophy fisherman, blogger, and skier.



Meghan
AMT Customer




G-JET® Family

Connect with Us Online!



Traditional Length G-JET®

Adjustable Length Gastric-Jejunal Enteral Feeding Tube



Comprehensive French Sizes

Designed for both pediatrics and adults ranging from 14F - 22F.

1 Anti-Kink Technology


Jejunal section offers spring reinforcement to maintain flow (offered in 16F and 18F).

2 Exclusive AMT Balloon

Apple-shaped, expands horizontally creating a superior gastric seal.⁷

3 Adjustable Bolster

Accommodates changes in stoma length.



Traditional Length G-JET®


Adjustable Length Gastric-Jejunal Enteral Feeding Tube

The G-JET® You Trust, Now Offered in Traditional Length Designed to Reduce Clogging

GASTRIC Profile



JEJUNAL Profile



GASTRIC Lumen Comparison⁶

AMT 14F is **Exclusive**, Not Offered by Competitor

AMT 16F is **41% Larger** than Competitor 16F

AMT 18F is **34% Larger** than Competitor 18F

JEJUNAL Lumen Comparison⁶

AMT 14F is **56% Larger** than Competitor 16F
22% Larger than Competitor 18F

AMT 16F is **131% Larger** than Competitor 16F
81% Larger than Competitor 18F

AMT 18F is **172% Larger** than Competitor 18F
57% Larger than Competitor 22F

Traditional Length G-JET® | Order Information

JEJUNAL	GASTRIC 14F (Without Spring)		GASTRIC 16F (With Spring)		GASTRIC 18F (With Spring)		GASTRIC 22F (Without Spring)	
	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #
15	TGJ-14-15	TGJ-14-15-I	TGJ-16-15	TGJ-16-15-I				
22	TGJ-14-22	TGJ-14-22-I	TGJ-16-22	TGJ-16-22-I	TGJ-18-22	TGJ-18-22-I		
30	TGJ-14-30	TGJ-14-30-I	TGJ-16-30	TGJ-16-30-I	TGJ-18-30	TGJ-18-30-I	TGJ-22-30	TGJ-22-30-I
45	TGJ-14-45	TGJ-14-45-I	TGJ-16-45	TGJ-16-45-I	TGJ-18-45	TGJ-18-45-I	TGJ-22-45	TGJ-22-45-I

Traditional G-JET® comes packaged as 1/box (sterile). | Kit contents include: device (1), 5cc/10cc balloon syringe (1), 35cc bolus syringe (1), and gauze. NOTE: G-JET® devices should never be cut to custom lengths.


G-JET® Family

MiniACE® Enema Button

Low Profile Antegrade Continenence Enema Button



A New Standard of Care in Bowel Management
The Only FDA Cleared Balloon Button Device for ACE Procedures



1 10F Tubing

Minimizes stoma dilation associated with conversion from other devices to MiniACE®.

2 Exclusive AMT Balloon


Apple-shaped, expands horizontally to create a superior seal which is designed to reduce leakage.¹³

3 Protected Tip

Balloon bolster internally shields the patient with its soft silicone design, protecting the cecum.

The Difference is Clear:


MiniACE® Enema Button: Low Profile Antegrade Continenence Enema Button



MiniACE® Enema Button | Order Information

LENGTH (cm)	10F		12F		14F	
	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #
0.8	ACE-1008	ACE-1008-I	ACE-1208	ACE-1208-I	ACE-1408	ACE-1408-I
1.0	ACE-1010	ACE-1010-I	ACE-1210	ACE-1210-I	ACE-1410	ACE-1410-I
1.2	ACE-1012	ACE-1012-I	ACE-1212	ACE-1212-I	ACE-1412	ACE-1412-I
1.5	ACE-1015	ACE-1015-I	ACE-1215	ACE-1215-I	ACE-1415	ACE-1415-I
1.7	ACE-1017	ACE-1017-I	ACE-1217	ACE-1217-I	ACE-1417	ACE-1417-I
2.0	ACE-1020	ACE-1020-I	ACE-1220	ACE-1220-I	ACE-1420	ACE-1420-I
2.3	ACE-1023	ACE-1023-I	ACE-1223	ACE-1223-I	ACE-1423	ACE-1423-I
2.5	ACE-1025	ACE-1025-I	ACE-1225	ACE-1225-I	ACE-1425	ACE-1425-I
2.7	ACE-1027	ACE-1027-I	ACE-1227	ACE-1227-I	ACE-1427	ACE-1427-I
3.0	ACE-1030	ACE-1030-I	ACE-1230	ACE-1230-I	ACE-1430	ACE-1430-I
3.5	ACE-1035	ACE-1035-I	ACE-1235	ACE-1235-I	ACE-1435	ACE-1435-I
4.0	ACE-1040	ACE-1040-I	ACE-1240	ACE-1240-I	ACE-1440	ACE-1440-I
4.4	ACE-1044	ACE-1044-I	ACE-1244	ACE-1244-I	ACE-1444	ACE-1444-I
5.0	ACE-1050	ACE-1050-I	ACE-1250	ACE-1250-I	ACE-1450	ACE-1450-I
5.5	ACE-1055	ACE-1055-I	ACE-1255	ACE-1255-I	ACE-1455	ACE-1455-I
6.0	ACE-1060	ACE-1060-I	ACE-1260	ACE-1260-I	ACE-1460	ACE-1460-I
6.5	ACE-1065	ACE-1065-I	ACE-1265	ACE-1265-I	ACE-1465	ACE-1465-I
7.0	ACE-1070	ACE-1070-I	ACE-1270	ACE-1270-I	ACE-1470	ACE-1470-I
8.0	ACE-1080	ACE-1080-I	ACE-1280	ACE-1280-I	ACE-1480	ACE-1480-I
9.0	ACE-1090	ACE-1090-I	ACE-1290	ACE-1290-I	ACE-1490	ACE-1490-I
10.0	ACE-10100	ACE-10100-I	ACE-12100	ACE-12100-I	ACE-14100	ACE-14100-I

Irrigation Set Included in Kit



MiniACE® Button comes packaged as 1/box (sterile). | Kit contents include: device (1), irrigation set (1), 5cc syringe (1), introducer/stiffener (1), gauze, and water soluble lubricant (1).

www.AppliedMedical.net

MiniACE® Button


Legacy Extension Sets

Feeding – Decompression – Medication






8-XXXX extension sets are compatible with:

MiniONE® Gastric Port (G-Port) + G-JET® Jejunal Port (J-Port)



Legacy Extension Sets | Order Information

LEGACY			
DESCRIPTION	ORDER #	BOX/QTY	
 2" Straight Purple Bolus Medication Set	8-0211	10	
 12" Straight Purple Bolus Feed Set	8-1211	10	
 24" Straight Purple Bolus Feed Set	8-2411	10	
 12" Right Angle Purple Bolus Feed Set	8-1221	10	
 24" Right Angle Purple Bolus Feed Set	8-2421	10	
LEGACY			
DESCRIPTION	ORDER #	BOX/QTY	
 2" Straight Clear Y-Port XL Medication Set	8-0212	10	
 12" Straight Clear Y-Port XL Feed Set	8-1212	10	
 2" Right Angle Clear Y-Port XL Medication Set	8-0222	10	
 12" Right Angle Clear Y-Port XL Feed Set	8-1222	10	
 24" Right Angle Clear Y-Port XL Feed Set	8-2422	10	
LEGACY			
DESCRIPTION	ORDER #	BOX/QTY	
 12" Right Angle Purple Y-Port Feed Set	8-1255	10	
 12" Right Angle Purple Y-Port (w/ Enhanced Med-Port)	8-1257G	10	
 24" Right Angle Purple Y-Port Feed Set	8-2455	10	


ENFit® Extension Sets

Feeding – Decompression – Medication







8-XXXX-ISOSAF extension sets are compatible with:

MiniONE® Gastric Port (G-Port) + G-JET® Jejunal Port (J-Port)



ENFit® Extension Sets | Order Information


ENFIT®			
DESCRIPTION	ORDER #	BOX/QTY	
 2" Straight Single ENFit® Medication Set w/ Press Cap	8-0211-ISOSAF	10	
 12" Straight Single ENFit® Feed Set w/ Press Cap	8-1211-ISOSAF	10	
 12" Right Angle Single ENFit® Feed Set w/ Press Cap	8-1221-ISOSAF	10	
 24" Straight Single ENFit® Feed Set w/ Press Cap	8-2411-ISOSAF	10	
 24" Right Angle Single ENFit® Feed Set w/ Press Cap	8-2421-ISOSAF	10	
 12" Right Angle Single ENFit® Feed Set w/ Twist Cap	8-1226-ISOSAF	10	
 12" Right Angle Single ENFit® Feed Set	8-1223-ISOSAF	10	
 24" Right Angle Single ENFit® Feed Set	8-2423-ISOSAF	10	
ENFIT®			
DESCRIPTION	ORDER #	BOX/QTY	
 2" Straight Dual ENFit® Y-Port Medication Set	8-0212-ISOSAF	10	
 12" Straight Dual ENFit® Y-Port Feed Set	8-1212-ISOSAF	10	
 2" Right Angle Dual ENFit® Y-Port Medication Set	8-0222-ISOSAF	10	
 12" Right Angle Dual ENFit® Y-Port Feed Set	8-1222-ISOSAF	10	
 24" Right Angle Dual ENFit® Y-Port Feed Set	8-2422-ISOSAF	10	
 12" Right Angle Dual ENFit® Y-Port Feed Set	8-1255-ISOSAF	10	
 24" Right Angle Dual ENFit® Y-Port Feed Set	8-2455-ISOSAF	10	
NON-PVC ENFIT®			
DESCRIPTION	ORDER #	BOX/QTY	
 12" Straight Single ENFit® Non-PVC Feed Setw/ Press Cap	8-1211-ISOSAF-NP	10	
 12" Right Angle Dual ENFit® Y-Port Non-PVC Feed Set	8-1255-ISOSAF-NP	10	
 24" Right Angle Dual ENFit® Y-Port Non-PVC Feed Set	8-2455-ISOSAF-NP	10	

Legacy Extension Sets

Decompression – Medication



6-XXXX extension sets are compatible with:
G-JET® Gastric Port (G-Port)



Legacy Extension Sets | Order Information


LEGACY			
	DESCRIPTION	ORDER #	BOX/QTY
	2" Straight Purple Bolus Medication Set	6-0211	10
	12" Straight Purple Bolus Extension Set	6-1211	10
	24" Straight Purple Bolus Extension Set	6-2411	10
	12" Right Angle Purple Bolus Extension Set	6-1221	10
	24" Right Angle Purple Bolus Extension Set	6-2421	10
LEGACY			
	DESCRIPTION	ORDER #	BOX/QTY
	2" Straight Clear Y-Port Medication Set	6-0212	10
	12" Straight Clear Y-Port Extension Set	6-1212	10
	12" Right Angle Clear Y-Port Extension Set	6-1222	10
	24" Right Angle Clear Y-Port Extension Set	6-2422	10

ENFit® Extension Sets

Decompression – Medication



6-XXXX-ISOSAF extension sets are compatible with:
G-JET® Gastric Port (G-Port)





ENFit® Extension Sets | Order Information

ENFIT®			
	DESCRIPTION	ORDER #	BOX/QTY
	2" Straight Single ENFit® Medication Set w/ Press Cap	6-0211-ISOSAF	10
	12" Straight Single ENFit® Extension Set w/ Press Cap	6-1211-ISOSAF	10
	24" Straight Single ENFit® Extension Set w/ Press Cap	6-2411-ISOSAF	10
	12" Right Angle Single ENFit® Extension Set w/ Press Cap	6-1221-ISOSAF	10
	24" Right Angle Single ENFit® Extension Set w/ Press Cap	6-2421-ISOSAF	10
	12" Right Angle Single ENFit® Extension Set w/ Twist Cap	6-1226-ISOSAF	10
	12" Right Angle Single ENFit® Extension Set	6-1223-ISOSAF	10
	24" Right Angle Single ENFit® Extension Set	6-2423-ISOSAF	10
ENFIT®			
	DESCRIPTION	ORDER #	BOX/QTY
	2" Straight Clear Dual ENFit® Y-Port Medication	6-0212-ISOSAF	10
	12" Straight Clear Dual ENFit® Y-Port Extension Set	6-1212-ISOSAF	10
	12" Right Angle Clear Dual ENFit® Y-Port Extension Set	6-1222-ISOSAF	10
	24" Right Angle Clear Dual ENFit® Y-Port Extension Set	6-2422-ISOSAF	10
NON-PVC ENFIT®			
	DESCRIPTION	ORDER #	BOX/QTY
	12" Right Angle Dual ENFit® Y-Port Non-PVC Feed Set	6-1222-ISOSAF-NP	10


CINCH®
Tube Securement Device




From **4F to 30F**, The CINCH® Has You Covered



The CINCH® comes fully assembled with an easy-peel split backing.



AMT CINCH® | Order Information



MEDIUM CINCH® TUBE SECUREMENT DEVICE		
DESCRIPTION	ORDER #	BOX/QTY
Medium CINCH® 4F – 17F	CINCH417M-10	10
Medium CINCH® 4F – 17F	CINCH417M-50	50
LARGE CINCH® TUBE SECUREMENT DEVICE		
DESCRIPTION	ORDER #	BOX/QTY
Large CINCH® 9F – 30F	CINCH930L-10	10
Large CINCH® 9F – 30F	CINCH930L-50	50

The AMT CINCH® was designed with the most delicate skin in mind. Our hypoallergenic silicone adhesive provides a gentle yet secure attachment. Secures 4F to 30F medical tubes.

1 Versatile

Medium CINCH® secures tubes 4F to 17F,
Large CINCH® secures tubes 9F to 30F.

2 Comfortable

Soft, low profile, flexible.

3 Breathable Skin


In two convenient sizes, Medium (3 cm x 4 cm)
and Large (4 cm x 6 cm).

4 Reusable

Elastomer strap can be repeatedly opened
and closed.

Hypoallergenic

Gentle adhesive for sensitive skin.



Customer Success Stories

AMT's Commitment: Real Customers. Real Experiences.

I fought my enteral company for 11 months to get the correct feed sets for my button. I contacted AMT and had the correct feed attachments in 1.5 weeks! They are GREAT at working with customers.



April
AMT Customer



Connect with Us Online!



AMT Accessories

Order Information



AMT CLAMP™ SECUREMENT DEVICE		
DESCRIPTION	ORDER #	BOX/QTY
Tube Securement Device	4-3000	10


The AMT Clamp™ device can reduce disconnects between a “Christmas Tree”/Stepped adapter and a feed set.



AMT SYRINGES		
DESCRIPTION	ORDER #	BOX/QTY
5 ml Slip Tip Syringe	B1081	24
10 ml Slip Tip Syringe	B1380	24
35 ml Slip Catheter Syringe	B1752	24
60 ml Slip Catheter Syringe	B1147	24

WINGED LUER ADAPTER		
DESCRIPTION	ORDER #	BOX/QTY
Winged Luer Adapter on 10 cc Syringe	4-7000	10

The Winged Luer Adapter allows for direct G and J access in order to inject contrast dye to confirm proper placement of the GJ device.



GUIDEWIRE		
DESCRIPTION	ORDER #	BOX/QTY
0.018” x 45 cm	4-GW-1	5
0.038” x 180 cm	4-GW-5	5


AMT Accessories

Order Information



TRN101 TRN102 TRN201 TRN202 TRN203 TRN204

ENFit® TRANSITION ADAPTERS		
DESCRIPTION	ORDER #	BOX/QTY
Male ENFit® Transition Adapter	TRN101	10
ENFit® Transition Stepped Adapter	TRN102	10
Female ENFit® Transition Adapter	TRN201	10
ENFit® Bolus Transition Adapter	TRN202	10
ENFit® Y-Port Transition Adapter	TRN203	10
ENFit® Y-Port Adapter w/ Enhanced Med-Port	TRN204	10



LEGACY Purple Y-Port Adapter LEGACY Clear Y-Port XL Adapter ENFit® Single-Port Adapter ENFit® Single-Port to Y-Port Adapter


TRADITIONAL LENGTH G-TUBE REPLACEMENT ADAPTERS		
DESCRIPTION	ORDER #	BOX/QTY
Capsule Dome, Capsule Monarch®		
12F Purple Y-Port Adapter	4-2012	10
14F Purple Y-Port Adapter	4-2014	10
15F Clear Y-Port XL Adapter	4-2015	10
20F Clear Y-Port XL Adapter	4-2020	10
24F Clear Y-Port XL Adapter	4-2024	10
12F ENFit® Single-Port Adapter	4-2012-ISOSAF	10
14F ENFit® Single-Port Adapter	4-2014-ISOSAF	10
20F ENFit® Single-Port Adapter	4-2020-ISOSAF	10
24F ENFit® Single-Port Adapter	4-2024-ISOSAF	10
Balloon, Capsule Dome, Capsule Monarch®		
ENFit® Single-Port to Y-Port Adapter	4-0005-ISOSAF	10

AMT Surgical: Initial Placement Kits

Initial Placement Gastrostomy Sets



AMT Quality, Right from the Start.




1 Multiple Size Dilators

For use with 10F –18F devices.




2 Ease of Placement

An AMT exclusive, size-specific, guidewire compatible Introducer/Stiffener.



3 Ensure a Proper Fit

Use AMT's Balloon Stoma Measuring Device with each placement.



AMT Surgical: Initial Placement Kits

Initial Placement Gastrostomy Sets

Designed to Minimize Trauma Associated with Stoma Tract Development

- Tapered Dilators Help Ensure Proper Dilatation
- Requires Only One Endoscopic Pass
- Compatible with Laparoscopic, Endoscopic, and Fluoroscopic Placement
- For Use with 10F - 18F Devices

Initial Placement Kits | Order Information

COMPREHENSIVE I.P. KITS*	
ORDER #	KIT CONTENTS
IP-10-S	(3) Dilators (7F, 12F, 14F-S)
	(1) Guidewire, 70 cm (0.018")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(1) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
	(4) 4" x 4" Gauze
ORDER #	KIT CONTENTS
IP-1218	(5) Dilators (7F, 12F, 16F, 18F, 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
(1) 5 cc Syringe	
(4) 4" x 4" Gauze	
ORDER #	KIT CONTENTS
IP-1218-S	(5) Dilators (7F, 12F, 16F, 18F-S , 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
(1) 5 cc Syringe	
(4) 4" x 4" Gauze	

I.P. Kits are packaged 1/box (sterile). | *NOTE: Kits do not contain T-Fasteners or U-Stitches. "-S" I.P. Kits include one sheathed dilator, noted in bold above.

BASIC I.P. KITS*	
ORDER #	KIT CONTENTS
IP-DIL-10	(3) Dilators (7F, 12F, 14F)
	(1) Guidewire, 45 cm (0.018")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible Introducer/Stiffener
ORDER #	KIT CONTENTS
IP-DIL	(4) Dilators (7F, 12F, 16F, 20F)
	(1) Guidewire, 45 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible Introducer/Stiffener


I.P. Kits are packaged 1/box (sterile). | *NOTE: Kits do not contain T-Fasteners or U-Stitches.

I.P. KIT ACCESSORIES – AVAILABLE SEPARATELY	
ORDER #	DESCRIPTION
1804	Balloon Stoma Measuring Device Guidewire Compatible (up to 0.038")
1704	Balloon Stoma Measuring Device Non-Guidewire Compatible
4-6010	Introducer/Stiffener Guidewire Compatible (up to 0.018") Device Compatibility: 10F and larger Low Profile Devices
4-6000	Introducer/Stiffener - 4.6" Length Guidewire Compatible (up to 0.038") Device Compatibility: 12F and larger Low Profile Devices
4-6200*	Introducer/Stiffener - 9.6" Length Guidewire Compatible (up to 0.038") Device Compatibility: 12F and larger Low Profile Devices
4-5000	Introducer/Stiffener Non-Guidewire Compatible Device Compatibility: 12F and larger MINI™ Button Devices

I.P. Kit accessories are individually pouched and sold 10/box (sterile). *Note: 4-6200 is not CE marked.

AMT Surgical: TLC® Retractors and Stays

Self-Retaining Retractor System



Faster Deployment

Flexible hinges adjust to the patient's anatomy in an instant. There are no thumbscrews to adjust!




Improve Field of Vision

Low profile design allows for greater visualization of the surgical site.



Eliminate Slippage

TLC® ribbed stays provide a secure fit, virtually eliminating slippage and the risk of popping out.



Increase Control

Uniform tension allows better control of the tissue as the retractor is precisely positioned around the surgical site.

AMT Surgical: TLC® Retractors and Stays

Self-Retaining Retractor System

TLC® Retractors and Stays | Order Information

AMT SURGICAL LINE – TLC® RETRACTORS									
43 cm	LARGE	8	6	8	MEDIUM	8	6	SMALL	PEDIATRIC
36 cm									
22 cm									
15 cm									
0 cm									
ORDER #	TLC900	TLC300	TLC200	TLC600	TLC400	TLC800	TLC500	TLC700	TLC920

Disposable retractor rings are pre-sterilized and packaged individually in Tyvek® pouches for easy delivery to the sterile surgical field.

TLC® RETRACTOR STAYS		
ORDER #	DESCRIPTION	
TLC090	3 mm Sharp	
TLC110	5 mm Sharp	
TLC100	5 mm Blunt	
TLC150	12 mm Blunt	
TLC203	20 mm Rake	

Stays are available in 4, 8, and 50 packs.
The 20 mm 3-Finger Rake is available in a 4 pack.
4 and 8 packs: Stays are assembled in a tray and pouched (sterile).
50 pack: Stays are individually pouched (sterile).

TLC® UROGYNECOLOGY RETRACTOR SYSTEM		
ORDER #	DESCRIPTION	BOX/QTY
TLC5015	Retractor Frame, TLC200	1
	5 mm Sharp Stay	8

TLC® COLORECTAL RETRACTOR SYSTEM		
ORDER #	DESCRIPTION	BOX/QTY
TLC5017	Retractor Frame, TLC700	1
	5 mm Sharp Stay	8

WILSON™ ENHANCED RETRACTOR SYSTEMS		
ORDER #	DESCRIPTION	BOX/QTY
TLC5011	Retractor Frame, TLC900	1
	12 mm Blunt Stay	6
	5 mm Sharp Stay	1
	Small Deavers (17 mm Wide)	2
	20 mm 3-Finger Rake	2
	Penile Elevation Strap	1
ORDER #	DESCRIPTION	BOX/QTY
TLC5013	12 mm Blunt Stay	6
	5 mm Sharp Stay	1
	Small Deavers (17 mm Wide)	2
	20 mm 3-Finger Rake	2
	Penile Elevation Strap	1

WILSON™ RETRACTOR SYSTEMS		
ORDER #	DESCRIPTION	BOX/QTY
TLC5042	Retractor Frame, TLC200	1
	12 mm Blunt Stay	6
	5 mm Sharp Stay	1
	Penile Elevation Strap	1
ORDER #	DESCRIPTION	BOX/QTY
TLC5040	12 mm Blunt Stay	6
	5 mm Sharp Stay	1
	Penile Elevation Strap	1


AMT Surgical: AMT Remora®

Smoke and Fluid Evacuation System



Safety Is Our Highest Priority!

Each Year **500,000+** Healthcare Workers Are Exposed to Surgical Smoke⁸



- Water Vapor
- Contaminants/Carcinogens
- Electrosurgical Smoke Environment

ONE Day of **O.R.** Surgical Smoke Exposure is Equivalent to Smoking **27-30** Unfiltered Cigarettes a Day.^{10, 11}

What's in Surgical Smoke?

- Nearly 150 Chemical Classes⁹
- Alkanes
- Alkenes
- Aldehydes
- Ketones
- Alkyl aromatics
- Pyrroles
- Furans
- Pyridines
- Thiols
- Thiocyanates
- Pyrimidines
- Nitriles
- 600+ Unidentified Compounds

Breathe Easy with Remora®

Smoke + Fluid Evacuation

- Removes smoke and fluid at the source.
- Low profile design allows for clear visualization of the surgical site.
- No impact to the surgeon's dexterity.
- Complies with A.O.R.N. and OSHA guidelines.

Universal Compatibility

- The Remora® simply slips over almost any electrosurgical pencil.
- Universal suction adapter which plugs into most O.R. suction sources.
- No capital equipment investment required to incorporate the Remora®.

AMT Surgical: AMT Remora®

Smoke and Fluid Evacuation System




ST9000 Remora®
Smoke and Fluid Evacuator with Pinch Clamp (Shown below in Universal Compatibility).




ST9100 Remora® Micro*
Smoke and Fluid Evacuator for Tight, Detailed Procedures.

Universal Compatibility



Remora® | Order Information



ST9000 REMORA®	
ORDER #	BOX/QTY
ST9000	10
ST9100 REMORA® MICRO*	
ORDER #	BOX/QTY
ST9100*	10

Electrosurgical unit is not included.
 Remora® includes: Evacuator (1), needle tip adapter (1), and pinch clamp (1).
 Remora® Micro includes: Evacuator (1), and needle tip adapter (1).
 Remora® can accommodate many different tips (needle, blade, or extended) and electrosurgical pencil handles.
 *ST9100 is not licensed for sale in Canada.

Dedication to Education

Educational Resources and Support



Online Resources + Support:

- AMT ONE Source® Mobile App
- Patient Education Guides
- How-To Videos
- Downloadable Content
- FAQs and MORE!

Educational Training Tools:

- Acrylic Comparison Stand
- Mini BUDDY® Bears
- AMT Training Dolls
- Bridle Pro® Skull Kit
- AMT Education Kit (Order# P-EDU KIT)*



Download the **AMT ONE Source®** App Today!
Simply Scan QR Code with Smartphone

App Store
DOWNLOAD
Scan QR Code For iPhone



Google Play
DOWNLOAD
Scan QR Code For Android



"We are here for you."

- Dr. Picha




Connect with Us Online!



Contact AMT

Need Help with a Product or Have a Question?



GET IN TOUCH WITH AMT, WE'RE HERE TO HELP!

At AMT, our team is here for YOU! We are committed to empowering patients, caregivers, and medical professionals through educational resources and support. AMT offers a broad range of online educational resources for both healthcare providers and patients, including a vast video library.

Our trained and certified Enteral Product Specialists are available to offer in-person or virtual training and technical support. Additionally, our Customer Service team and Nurse Practitioners are ready to answer your questions.

Reach Out Today!



CALL

TOLL FREE: 800-869-7382
LOCAL: 440-717-4000



FAX

TOLL FREE: 440-717-4200



EMAIL

USA: CS@AppliedMedical.net
INTERNATIONAL:
ICS@AppliedMedical.net

Contact AMT

Setting Up an Account – Placing an Order

Healthcare professionals can place an order via fax or email. To place an order, follow the simple three-step process outlined below.



SET UP AN ACCOUNT:

New customers will need to submit a New Account Application prior to placing their first order. Interested parties can contact AMT's Customer Service team to obtain the necessary documents.



CHOOSE PAYMENT METHOD:

AMT accepts payment via NET 30 terms or Prepayment/Credit Card. To establish terms, customers will need to complete a Credit Application and include three favorable references.



PLACE YOUR ORDER:

Items that are in stock will ship the next business day. Overnight and Second Day Air shipping are available for orders received before 2:00 p.m. Eastern Time (only applicable to in-stock items). All orders are shipped via UPS FOB Origin. No third-party shipping.

AMT's knowledgeable Customer Service team is available from 8:00 a.m. to 5:00 p.m. Eastern Time to assist with the purchasing process.

INSURANCE BILLING AND REIMBURSEMENT

AMT does not process medical insurance billing paperwork. For more information on coding options for AMT devices, please visit: <https://www.appliedmedical.net/order>

Contact AMT

Product Education – Product Questions – Who is My AMT Representative?



NEED PRODUCT EDUCATION OR TRAINING?

We are here to help! Our live Customer Service team can assist with placing an order, answering product questions, and getting you in contact with your local AMT Rep.



HAVE A PRODUCT-SPECIFIC QUESTION OR NEED ASSISTANCE?

At AMT, our team is here for YOU! Our team can help with device troubleshooting questions and provide additional educational resources.



WANT TO MEET YOUR LOCAL AMT REPRESENTATIVE?

Our fully-credentialed Enteral Product Specialists are available to offer in-person or virtual training and technical support. Reach out to AMT and we'll get you connected!

WHAT'S YOUR STORY? WE'D LOVE TO HEAR IT!

Take a moment to share your personal experience and testimonial with us!

www.AppliedMedical.net/Stories/Share

References:

1. Seder, C.W., et al. Crit Care Med 2010;38:797-801
2. Badosa, E.L. Nutr Hosp. 2017;34(4):907-913
3. Claim based on internal testing of 5F – 18F devices.
4. Abdelhadi, Ruba A., Dekonenko, C., Dorman, Robert M., et al. A Prospective, Randomized Cross-over Trial of Two Low-profile Gastrostomy Buttons to Determine Family Preference. JPGN 2020;70:386 – 388.
5. The micro G-JET® is cleared for patients > 6 kg.
6. Claim based on internal testing of 14-22F devices.
7. Claim based on internal testing of 16F devices.
8. Hill D.S., O'Neill, J.K., Powell, R.J., et al. Surgical smoke—a health hazard in the operating theatre: a study to quantify exposure and a survey of the use of smoke extractor systems in UK plastic surgery units. Journal of Plastic Reconstructive & Aesthetic Surgery. 2012; 65(7):911-916.
9. Pierce J.S., Lacey S.E., Lippert, J.F., et al. Laser-Generated Air Contaminants from Medical Laser Applications: A State-of-the-Science Review of Exposure Characterization, Health Effects, and Control. JOEH. 2011; 8:447-466
10. Tomita, Y., Mihashi, S., Nagata, K., et al. Mutagenicity of smoke condensates induced by CO2-laser irradiation and electrocauterization. Mutation research 89 2 (1981):145-149.
11. Hill D.S., O'Neill, J.K., Powell, R.J., et al. Surgical smoke—a health hazard in the operating theatre: a study to quantify exposure and a survey of the use of smoke extractor systems in UK plastic surgery units. Journal of Plastic Reconstructive & Aesthetic Surgery. 2012; 65(7):911-916.
12. Based on internal testing of 14F, 18F, 20F, and 24F devices.
13. Claim based on internal testing of 10F devices.

Customer Success Stories

AMT's Commitment: Real Customers. Real Experiences.

Huge shout-out to the amazing team at AMT! Can't say enough wonderful things about this company! Not only are their products top notch, but their customer service is unparalleled.



Amanda,
AMT G-JET® User



Connect with Us Online!



NOTES

Lined area for notes on page 50.



Lined area for notes on page 51.



NOTES

Lined notes area on page 52, featuring horizontal lines and a light blue shaded background that tapers from left to right.

Lined notes area on page 53, featuring horizontal lines and a light blue shaded background that tapers from right to left.