



# MiniACE<sup>®</sup> Button Cecostomy + Appendicostomy

# MiniACE<sup>®</sup> Button Low Profile Antegrade Continence Enema Button



## MiniACE<sup>®</sup> Button Sizes

New! 10F Stoma Size and Longer-Length 12F Devices!

AMT	10F	12F	14F
LENGTH	LEGACY ORDER NUMBERS*		
0.8 cm	ACE-1008	ACE-1208	ACE-1408
1.0 cm	ACE-1010	ACE-1210	ACE-1410
1.2 cm	ACE-1012	ACE-1212	ACE-1412
1.5 cm	ACE-1015	ACE-1215	ACE-1415
1.7 cm	ACE-1017	ACE-1217	ACE-1417
2.0 cm	ACE-1020	ACE-1220	ACE-1420
2.3 cm	ACE-1023	ACE-1223	ACE-1423
2.5 cm	ACE-1025	ACE-1225	ACE-1425
2.7 cm	ACE-1027	ACE-1227	ACE-1427
3.0 cm	ACE-1030	ACE-1230	ACE-1430
3.5 cm	ACE-1035	ACE-1235	ACE-1435
4.0 cm	ACE-1040	ACE-1240	ACE-1440
4.4 cm	ACE-1044	ACE-1244	ACE-1444
5.0 cm	ACE-1050	ACE-1250	ACE-1450
5.5 cm	ACE-1055	ACE-1255	ACE-1455
6.0 cm	ACE-1060	ACE-1260	ACE-1460
6.5 cm	ACE-1065	ACE-1265	ACE-1465
7.0 cm	ACE-1070	ACE-1270	ACE-1470
8.0 cm	ACE-1080	ACE-1280	ACE-1480
9.0 cm	ACE-1090	ACE-1290	ACE-1490
10.0 cm	ACE-10100	ACE-12100	ACE-14100

\*Add "-I" for ENFit<sup>®</sup> MiniACE<sup>®</sup> Button kits. MiniACE<sup>®</sup> Kits are packaged with a 5 cc syringe, an introducer, gauze, and a lube packet. A sterile pouched 2-1221 bolus irrigation set is also included in the kit.

### Irrigation Set Included in Kit



### Complete Portfolio – Irrigation Sets

Adapter Options	Right Angle Connector Options**	
	12"	24"
Bolus	2-1221	2-2421

\*\*Add "-ISOSAF" for ENFit<sup>®</sup> Irrigation Sets.

**Note:** Irrigation sets are packaged as 10/box (non-sterile)



## Guts of a Warrior SUCCESS STORIES

River-Kai proudly rocks two AMT devices, the MiniONE<sup>®</sup> and the MiniACE<sup>®</sup>!!


[Our surgeons have commented] that they've never seen a stoma look so good. I owe that to AMT for such a great product! I recommend AMT products to any parent who complains about granulation. – River-Kai's Mom

**River-Kai's MiniACE<sup>®</sup>**  
AppliedMedical.net/River-Kais-MiniONE-and-MiniACE

### Balloon Fill Volumes

French Size	Min	Rec.	Max	
10F	1 ml	1.5 ml	2 ml	
12F	2 ml	2.5 ml	3 ml	
14F	3 ml	4 ml	5 ml	

Recommended water fill volume is clearly printed on the balloon fill valve of each MiniACE<sup>®</sup> Button. Consult your healthcare provider before adjusting fill volumes.



**Cecostomy + Appendicostomy**

**Cleared** for Initial Placement

**Superior** Leakage Control

**Improved** Stoma Retention

For a complete list of products visit us at [www.AppliedMedical.net](http://www.AppliedMedical.net)  
800 869 7382 | [CS@AppliedMedical.net](mailto:CS@AppliedMedical.net) | [ICS@AppliedMedical.net](mailto:ICS@AppliedMedical.net)

Initial Placement Cleared

MR Conditional



Applied Medical Technology, Inc.  
8006 Katherine Boulevard  
Brecksville, OH 44141  
[www.AppliedMedical.net](http://www.AppliedMedical.net)  
An ISO 13485 Registered Company





Innovating. Educating. Changing Lives.™

© 2020 Applied Medical Technology, Inc. C-4642-D  
Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>  
Chait Trapdoor™ is a trademark of Cook® Medical | ENFit® is a registered trademark of GEDSA. | 1.Claim based on internal testing of 10F-14F devices.



### MiniACE® Button

- Cleared for Initial Placement
- **NOW IN 10F!** Minimizes stoma dilation associated with conversion from Chait Trapdoor™ to MiniACE®
- Cecostomy + Appendicostomy




### Patient Comfort First

- Medical grade silicone, gentle to the skin
- Flexible + Soft, easy to clean around the button
- Exclusive AMT Balloon, designed to reduce pullouts and leakage
- Balloon bolster shields the internal tip of the MiniACE®, protecting the cecum

### Reduce Leaks, Superior Retention

- AMT's MiniACE® Button is 450% more secure within the stoma when compared to the Chait Trapdoor™<sup>1</sup>
- Exclusive AMT Balloon, along with tapered tubing at the external bolster, reduce the risk of leaks.
- Anti-leakage valve helps prevent backflow.



### Safe and Secure Irrigation Set Attachment

- The connection between the MiniACE® Button and its irrigation set is 10x stronger than the competition.<sup>1</sup>

